## NPPTL COVID-19 Response: Beyond Shelf Life/Stockpiled Respirator Assessment

Manufacturer: Inovel, LLC

Model Tested: MFRN95-ML

Date Tested: May 27, 2020

Report Prepared: May 28, 2020

These findings pertain to the Inovel, LLC, model MFRN95-ML (TC-84A-4100) submitted for testing, and may not be applicable to other stockpile facilities and/or under different environmental storage conditions. The maximum and minimum filter efficiency was 99.61% and 96.99%, respectively. All twenty respirators measured more than 95% efficiency. No shelf life was designated for this model by the previous NIOSH approval holder.

NIOSH regulation sets the minimum quality and performance requirements for the approval of respirators (42 CFR 84). **NIOSH does not have requirements for shelf life or storage conditions for particulate-only APRs.** The approval holder (i.e. the entity that is granted the approval from NIOSH) is responsible for understanding how their products' design or performance may be affected by various use or storage conditions and must provide instruction for establishing the proper use, storage, and maintenance procedures for their approved products, which may include designating a shelf life. FFR or particulate filter packaging (such as the box) often includes NIOSH-approved user instructions, label information, and recommendations on shelf life. Additionally, some approval holders also disseminate recommendations related to storage and shelf life through resources such as user and web notices. The respirators tested in this study were generally not designed for long-term storage.

Based on research conducted by NIOSH and this limited testing, NIOSH does not have enough information to definitively know the level of protection that may be provided by respirators that 1) are stored for prolonged periods of times; 2) are stored under various storage conditions; or 3) have exceeded the approval holder's designated shelf life. Users of respirators that have exceeded the designated shelf life should be forewarned to avoid a false sense of confidence; these devices may not provide the same level of protection as those that have not exceeded the designated shelf life. We recommend contacting the approval holder(s) of the respirators in the stockpile with specific questions regarding the use of product beyond the manufacturer- designated shelf life.

The results provided in this letter are specific to the subset of NIOSH-approved N95s, past their designated shelf life, that were provided to NPPTL for evaluation.

These results will be added to the CDC guidance for [Stockpiled N95 Filtering Facepiece Respirators Beyond the Manufacturer-Designated Shelf Life](#).

## Evaluation of Stockpiled and Beyond Manufacturer-Designated Shelf Life N95s

**Test:** TEB-APR-STP-0059

**Manufacturer:** Inovel, LLC

**Item Tested:** MFRN95-ML

**Expiration Date:** None Provided

**Manufacture Date:** None Provided

Filter	Flow Rate (LPM)	Initial Filter Resistance (mmH <sub>2</sub> O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency (%)
1	85	8.6	0.97	0.97	99.03
2	85	8.2	0.95	0.95	99.05
3	85	9.4	0.90	0.90	99.10
4	85	8.3	1.42	1.42	98.58
5	85	8.4	1.03	1.03	98.97
6	85	8.3	1.97	1.97	98.03
7	85	10.0	0.52	0.52	99.48
8	85	8.7	0.87	0.87	99.13
9	85	8.3	1.32	1.32	98.68
10	85	10.2	0.39	0.39	99.61
11	85	8.2	1.23	1.23	98.77
12	85	7.8	3.01	3.01	96.99
13	85	8.2	1.00	1.00	99.00
14	85	8.8	0.98	0.98	99.02
15	85	8.3	0.94	0.94	99.06
16	85	9.5	0.64	0.64	99.36
17	85	8.0	1.34	1.34	98.66
18	85	10.3	0.67	0.67	99.33
19	85	8.3	2.01	2.01	97.99
20	85	8.2	1.00	1.00	99.00
<b>Minimum Filter Efficiency: 96.99</b>			<b>Maximum Filter Efficiency: 99.61</b>		

REF MFRN95-ML *MTT-2020-198.1* Qty. 20

# Secure-Gard® N95 Healthcare Particulate Respirator and Surgical Mask

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- ENG**
- **Extra Comfort** – Contoured Shape, Soft Foam Nose Cushion
  - **Easy to fit** – Molded Nose Bridge
  - **Color coded sizes**
  - **Fluid Resistant** – 160 mm Hg'
  - **BFE >99.9%**
  - **Dura-Mesh®** – Prevents glove/hand contact on filter media
  - **Meets CDC guidelines for TB exposure control**

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## WARNING

Failure to follow instructions and warnings on the proper use of this respirator may reduce respirator effectiveness, wearer protection and may result in sickness or death. For additional information consult your supervisor, read all information enclosed and phone 866-546-6835.

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Latex Free  
Sans Latex  
Sin Látex

Distributed by:  
**Cardinal Health**  
McGaw Park, IL 60085-6787 USA  
DMFRN95-ML-60

Manufactured by:

# NPPTL COVID-19 Response: Beyond Shelf Life/Stockpiled Respirator Assessment

## DESCRIPTION:

The Inovel/Cardinal Health MFRN95 Series healthcare N95 Particulate Respirator and Surgical Mask is designed to help provide respiratory protection for the wearer. This product has been tested<sup>1</sup> and certified by NIOSH as an N95 respirator and as having a filter efficiency level of 95% or greater against particulate aerosols free of oil. It is fluid resistant,<sup>2</sup> disposable and may be worn in surgery or throughout the hospital.

This product contains no components made from natural rubber latex.

## INTENDED USE:

The various models of Inovel/Cardinal Health Type N95 Healthcare Particulate Respirators and Surgical Masks meet CDC Guidelines for TB Exposure Control within healthcare facilities. These devices are also intended to be worn by healthcare personnel during surgical procedures to protect both the patient and healthcare personnel from the transfer of microorganisms, body fluids and particulate material.

## CONTRAINDICATIONS:

1. Not for use with beards, or other facial hair that prevents direct contact between the face and sealing surface of the respirator.
2. Eyewear must not prevent direct contact between the face and sealing surface of the respirator.
3. Not to be used on children.

- 1 Tested in accordance with NIOSH 42 CFR 84
- 2 Passed ASTM F 1862 @ 160mm Hg

## LIMITED WARRANTY – IMPORTANT NOTICE TO PURCHASER

This limited warranty is made in lieu of the warranties of merchantability, fitness for particular purposes and all other warranties, express or implied. There are no other warranties which extend beyond the description on the face hereof. The physical standards and specifications of Inovel will be met by products sold. Exclusive Remedies: damages for the breach of this limited warranty are limited to the replacement of such quantity of Inovel products proved to be defectively manufactured. Except as provided above, Inovel shall not be liable or responsible for any loss, damage, or liability, direct, indirect, incidental, special, or consequential, arising out of sale, use, or misuse, or the inability to use products by the user.

## INOVEL LLC

10111 W. Jefferson Blvd., Culver City, CA 90232 Phone 866-546-6835 310-204-5237, ext. 523 Fax 310-204-4337. E-mail: sales@inovelmedical.com www.inovelmedical.com

U.S. Patents #0285,374, #4,850,347 and #4,384,577 Trademarks and Patents of Moldex-Metric, Inc., Culver City, CA 90232 USA MFRN95ML-802 REV A 1/07 Made in U.S.A.

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INOVEL LLC  
CULVER CITY, CALIFORNIA, U.S.A.  
866-546-6835

These respirators are approved only in the following configurations:

TC	PROTECTION <sup>1</sup>	RESPIRATOR				CAUTIONS AND LIMITATIONS <sup>2</sup>
		MFRN95-XS	MFRN95-S	MFRN95-A	MFRN95-ML	
84A-4T00	N95	X	X	X	X	A B C J M N O P

- 1. Protection**  
N95 - Particulate Filter (95% filter efficiency level) efficient against particulate aerosols free of oil; time use restrictions may apply
- 2. Cautions and Limitations**
  - A. Not for use in atmospheres containing less than 19.5 percent oxygen.
  - B. Not for use in atmospheres immediately dangerous to life or health.
  - C. Do not exceed maximum use concentrations established by regulatory standards.
  - J. Failure to properly use and maintain this product could result in injury or death.
  - M. All approved respirators shall be selected, fitted, used and maintained in accordance with MSHA, OSHA, and other applicable regulations.
  - N. Never substitute, modify, add or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
  - O. Refer to User's Instructions and/or maintenance manuals for information on use and maintenance of these respirators.
  - P. NIOSH does not evaluate respirators as surgical masks.

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## WARNING TO USER:

Non-latex straps are used in the headbands of these respirators. However, individuals highly sensitive to natural rubber latex may have an allergic reaction.

This product does not eliminate the risk of contracting any disease or infection. Change immediately if the respirator comes in contact with blood or fluids, including body fluids. OSHA has not set a permissible exposure level for airborne biohazards.

This product has not been sold with warning or use instructions for personnel involved in industrial or related situations.

If you are considering such uses, first contact the Inovel Technical Services Department at 866-546-6835 or 310-204-5237, ext. 523.

### RESTRICTIONS:

1. Before use, a written respiratory protection program in accordance with 29 CFR 1910.134(f) must be implemented.
2. This respirator does not supply oxygen and must not be used in atmospheres containing less than 19.5% oxygen.
3. Prior to each use, inspect filter media and strap attachment area for tears and damage. Staple or strap attachment perforations do not affect NIOSH approval.
4. If respirator is damaged, distorted, a proper fit cannot be obtained or breathing becomes difficult, leave contaminated area as soon as possible and replace respirator.
5. If the respirator comes in contact with blood or fluids, including body fluids, leave contaminated area as soon as possible and discard and replace the respirator.
6. Do not alter, modify or abuse this respirator.
7. Store unused respirators in box in a clean, dry, non-contaminated area.
8. Dispose of respirator according to your employer's policy and local regulations.
9. Do not reuse or store for reuse or hang around neck.
10. When used for surgical procedures, discard after every use.
11. Use respirator before the "use by" expiration date printed on box.

For technical assistance call Inovel Technical Service Department, 866-546-6835 or 310-204-5237, ext. 523

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
B0149ISOJ2

Secure-Gard N95 Particulate Respirator  
Surgical Mask Blue 20/Box - New

ie,  
partement  
ovel,  
4-5237,

**FITTING**


- If you cannot obtain a proper fit do not enter the contaminated area. OSHA regulation 29 CFR 1910.134(f) requires that the user must be fit tested before being assigned a respirator. Use Inovel Qualitative Fit Test Kit or an OSHA approved Quantitative fit test such as TSI's PortaCount® with N95 Companion™.
- Users must follow instructions each time respirator is worn. It's the wearer's responsibility to obtain a proper fit.
- To choose the correct MFRN95 Series size/shape you must fit test and user seal check.
- When changing from one Inovel/Cardinal Health model or size to another or when changing from another manufacturer's product to Inovel/Cardinal Health, you must fit test with the new respirator.
- Not for use with beards, or other facial hair that prevents direct contact between the face and sealing surface of the respirator.
- Eyewear must not prevent direct contact between the face and sealing surface of the respirator.

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**MFRN95 SERIES DONNING INSTRUCTIONS**

1. Hold respirator in hand with molded nose contour (narrow end) at finger tips, allowing headstraps to fall below hand.
2. Place respirator under chin with molded nose contour (narrow end) up. Cushion must not be creased inside respirator. Raise top strap to top back of head. Pull shorter bottom strap over head, below ears, to around neck. Do not wear with only one strap because it may affect fit.
3. Adjust respirator for comfortable fit.
4. Each time user enters contaminated area seal should be fit checked. Cover front of respirator by cupping both hands. **INHALE SHARPLY.** A negative pressure should be felt inside respirator. If any leakage is detected at respirator edges, adjust straps by pulling back along the sides and/or reposition respirator. Repeat until sealed properly, otherwise see your supervisor. Entry into a contaminated area with an improper fit may result in sickness or death.

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


NPPTL COVID-19 Response: Beyond Shelf Life/Stockpiled Respirator Assessment

NPPTL COVID-19 Response: Beyond Shelf Life/Stockpiled Respirator Assessment

NPPTL COVID-19 Response: Beyond Shelf Life/Stockpiled Respirator Assessment