NPPTL COVID-19 Response: Beyond Shelf Life/Stockpiled Respirator Assessment

Manufacturer: 3M

Model Tested: 1860

Date Tested: July 29, 2020

Report Prepared: July 29, 2020

These findings pertain to the 3M, model 1860 submitted for testing; and may not be applicable to other stockpile facilities and/or under different environmental storage conditions. The maximum and minimum filter efficiency was 99.35% and 98.57%, respectively. All nineteen respirators measured more than 95% efficiency. There was no designated expiration date.

NIOSH regulation sets the minimum quality and performance requirements for the approval of respirators (42 CFR 84). **NIOSH does not have requirements for shelf life or storage conditions for particulate-only APRs.** The approval holder (i.e. the entity that is granted the approval from NIOSH) is responsible for understanding how their products' design or performance may be affected by various use or storage conditions and must provide instruction for establishing the proper use, storage, and maintenance procedures for their approved products, which may include designating a shelf life. FFR or particulate filter packaging (such as the box) often includes NIOSH-approved user instructions, label information, and recommendations on shelf life. Additionally, some approval holders also disseminate recommendations related to storage and shelf life through resources such as user and web notices. The respirators tested in this study were generally not designed for long-term storage.

Based on research conducted by NIOSH and this limited testing, NIOSH does not have enough information to definitively know the level of protection that may be provided by respirators that 1) are stored for prolonged periods of times; 2) are stored under various storage conditions; or 3) have exceeded the approval holder's designated shelf life. Users of respirators that have exceeded the designated shelf life should be forewarned to avoid a false sense of confidence; these devices may not provide the same level of protection as those that have not exceeded the designated shelf life. We recommend contacting the approval holder(s) of the respirators in the stockpile with specific questions regarding the use of product beyond the manufacturer- designated shelf life.

The results provided in this letter are specific to the subset of NIOSH-approved N95s, past their designated shelf life, that were provided to NPPTL for evaluation.

These results will be added to the CDC guidance for [Stockpiled N95 Filtering Facepiece Respirators Beyond the Manufacturer-Designated Shelf Life](#).

Evaluation of Stockpiled and Beyond Manufacturer-Designated Shelf Life N95s

Test: TEB-APR-STP-0059

Manufacturer: 3M

Item Tested: 1860

Expiration Date: None Provided

Manufacture Date: None Provided

Filter	Flow Rate (LPM)	Initial Filter Resistance (mmH ₂ O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency (%)
1	85	8.7	0.54	1.09	98.91
2	85	10.2	0.84	1.22	98.78
3	85	7.7	0.27	0.83	99.17
4	85	8.8	0.51	1.15	98.85
5	85	9.4	0.20	0.65	99.35
6	85	7.3	0.51	1.23	98.77
7	85	7.5	0.41	1.09	98.91
8	85	8.3	0.47	1.43	98.57
9	85	8.9	0.34	1.06	98.94
10	85	7.7	0.27	1.03	98.97
11	85	8.8	0.33	1.12	98.88
12	85	8.8	0.23	0.72	99.28
13	85	7.3	0.34	0.98	99.02
14	85	9.2	0.29	1.03	98.97
15	85	8.0	0.35	0.82	99.18
16	85	9.4	0.26	0.96	99.04
17	85	9.1	0.44	0.93	99.07
18	85	9.4	0.48	0.91	99.09
19	85	6.9	0.42	1.20	98.80
Minimum Filter Efficiency: 98.57			Maximum Filter Efficiency: 99.35		

3M FLUID RESISTANT/RÉSISTE AUX FLUIDES
RESISTENTE A LOS FLUIDOS/RESISTENTE A FLUIDOS

**Health Care Particulate Respirator
and Surgical Mask** 

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Particulate Respirator N95

NIOSH Approved: N95

At least 95% filtration efficiency against solid and liquid aerosols that do not contain oil.

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⚠ WARNING

This respirator helps protect against certain particulate contaminants but does not eliminate exposure to or the risk of contracting any disease or infection. **Misuse may result in sickness or death.** For proper use, see supervisor, or *User Instructions*, or call 3M in U.S.A., 1-800-247-3941. In Canada, call Technical Service at 1-800-267-4414.

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Exported by/Exporté par/Exportado por:
3M Occupational Health and
Environmental Safety Division
3M Center, Building 0235-02-W-70
St. Paul, MN 55144-1000
Made in the U.S.A. from U.S. and Imported Materials

Imported by:
3M Occupational Health and
Environmental Safety Division
3M Canada Company
P.O. Box 5757
London, Ontario N6A 4T1
Made in the U.S.A. from U.S. and Imported Materials

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1860 (regular, courant, mediano)

FLUID RESISTANT/RÉSISTE AUX FLUIDES

3M™ Health Care Particulate Respirator and Surgical Mask

IMPORTANT: Before use, the wearer must read and understand the *User Instructions* inside.

DESCRIPTION: The 3M™ 1860 Health Care Particulate Respirator and Surgical Mask is designed to help provide respiratory protection for the wearer. This product has a filter efficiency level of 95% or greater against particulate aerosols free of oil¹. It is fluid resistant, disposable and may be worn in surgery. It can fit a wide range of face sizes.

This product contains no components made from natural rubber latex.

INTENDED USE: This product meets CDC guidelines for *Mycobacterium tuberculosis* exposure control. As a respirator, it is intended to help reduce wearer exposure to certain airborne particles including those generated by electrocautery, laser surgery, and other powered medical instruments. As a surgical mask, it is designed to be fluid resistant to splash and spatter of blood and other infectious materials; when worn properly and in combination with protective eyewear, it complies with the OSHA Bloodborne Pathogens Standard. It also provides >99% BFE² against wearer generated micro-organisms.

CONTRAINDICATIONS: Not for industrial use. Not for use with beards or other facial hair that prevents direct contact between the face and the sealing surface of the respirator. This respirator can help reduce inhalation exposures to certain airborne biological particles (e.g. mold, avian influenza viruses, *Bacillus anthracis*, *Mycobacterium tuberculosis*, etc.) but cannot eliminate the risk of contracting infection, illness or disease. OSHA and other government agencies have not established safe exposure limits for these contaminants. This respirator is designed for occupational/professional use by adults who are properly trained in their use and limitations. This respirator was not designed to be used by children.

USE INSTRUCTIONS:

1. Before use for respiratory protection, a written respiratory protection program must be implemented meeting all requirements of OSHA 29 1910.134 such as medical evaluation, training and fit testing. In Canada, CSA standard Z94.4 requirements must be met. The 3M FT-10 (sweet solution) or FT-30 (bitter solution) Qualitative Fit Test Apparatus is recommended for this respirator. When used only as a surgical mask, fit testing is not required.
2. Respirator may be used until damaged, breathing becomes difficult, or contaminated with blood or body fluids. Otherwise, it may be stored and reused according to the facility's infection control policy. Discard after every use when used for surgical procedures.
3. Inspect respirator before each use to ensure that it is in good operating condition. Examine all the respirator parts for signs of damage including the two headbands, staples, noseclip, and nosefoam. The respirator should be disposed of immediately upon observation of damaged or missing parts.
4. Filtering facepieces are to be inspected prior to each use to assure there are no holes in the breathing zone other than the punctures around staples and no damage has occurred. Enlarged holes resulting from ripped or torn filter material around staple punctures are considered damage. Immediately replace respirator if damaged. Staple perforations do not affect NIOSH approval.

FOR MORE INFORMATION and assistance on 3M occupational health and environmental safety products, contact your local 3M representative or call 3M OH&ESD Technical Service toll free in U.S.A., 1-800-243-4630. In Canada, call Technical Service at 1-800-267-4414.

¹ Tested against a 0.3 micron particle (mass median aerodynamic diameter) per 42 CFR 84.

² Bacterial Filtration Efficiency (BFE) determined by the modified Greene and Vesley test method.

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©3M 2006