



**BACKUP DATA REPORT
NIOSH Method No. 6011**

Title: Bromine and Chlorine

Analytes: Bromine and Chlorine

Author/developer: Mary Ellen Cassinelli (Issues 1 & 2), Method revised by
Dawn Farwick (Issue 3)

Date: February 18, 2022

Disclaimer: *Mention of any company or product does not constitute endorsement by the National Institute for Occupational Safety and Health, Centers for Disease Control and Prevention. In addition, citations to websites external to NIOSH do not constitute NIOSH endorsement of the sponsoring organizations or their programs or products. Furthermore, NIOSH is not responsible for the content of these websites. All web addresses referenced in this document were accessible as of the publication date.*

Table 8. Method validation for bromine at low humidity (~20%)

Concentration Level (x PEL)	n	Concentration taken ^a (mg/m ³)	Concentration found (mg/m ³)	s _r	Recovery (%)
0.1	6	0.076	0.077	0.039	101
0.5	6	0.351	0.358	0.033	102
1	6	0.696	0.673	0.026	96.7
2	6	1.42	1.50	0.028	106

Mean recovery = 102%; Precision (pooled s_r) = 0.032

^aConcentration determined from generator source and dilution flows where sampling rate=0.3 L/min and sampling time=4 hr.

All data were within the accuracy criterion of ±25% at the 95% confidence level and coefficient of variation (precision) recommended for a valid NIOSH method [10]. Bartlett's test for homogeneity was applied to the data before pooling. The mean recovery for the overall method (Cl₂ and Br₂) was 98.7% (-1.3% bias) with a sampling and analytical precision of 5.2%, and a total coefficient of variation estimate (CV_T) of 7.2%.

Instrument LOD/LOQ and linear range

The limit of detection (LOD), the smallest amount of analyte that can be distinguished from background, was defined in this study as three times the standard error of the calibration graph divided by its slope. The limit of quantitation (LOQ) was defined as ten times the standard error of the calibration graph divided by its slope [9]. The LOD of Cl⁻ in the thiosulfate matrix was 0.06 µg/mL with a limit of quantitation (LOQ) of 0.19 µg/mL. For Br⁻, the LOD and LOQ were 0.16 and 0.54 µg/mL, respectively. The linear range for Cl⁻ at the sensitivity setting used in this work (10 µS) was 0.05 to 5 µg/mL, and for Br⁻, 0.20 to 15 µg/mL. The LOD and LOQ may be reduced and linear range proportionally extended at the lower end by using a more sensitive instrumental setting, e.g., 3 µS or 1 µS full scale. Similarly, the linear range can be extended at the upper end by using a less sensitive instrumental setting or by diluting the sample.

The lower limit for the working range of air concentrations for a given sampling volume can be calculated based on the limit of quantitation (LOQ) of the analytical method and the liquid sample volume. The lower limit for the working range of air concentrations (mg/m³) for a given sampling volume = LOQ (µg/mL) x liquid sample volume (mL) ÷ air sample volume (L):

For Cl₂: Lower limit of the working range for 45 L sampling volume = 0.19 µg/mL x 10 mL ÷ 45 L = 0.04 mg/m³.

For Br₂: Lower limit of the working range for 360 L sampling volume = 0.54 µg/mL x 10 mL ÷ 360 L = 0.015 mg/m³.

Cl₂ field study

Field samples for Cl₂ were collected at a plant that produces trichloroisocyanuric acid. The compound is prepared by treating cyanuric acid with Cl₂. Paired silver filter samples and sulfamic acid bubblers (OSHA Method ID-101) [12] were collected at rates of 0.5 and 1.0 L/min, respectively. The results are presented in Table 9.

Table 9. Chlorine field study results using paired silver filters and bubbler^a samples

Pair number	Sampler numbers	Air volume (L)	Mass Cl ₂ (µg)	Concentration (mg/m ³)	Difference ^b (mg/m ³)
1	Filter – 1	189	134	0.71	0.10
1	Bubbler – 1	377	230	0.61	
2	Filter – 3	211	122	0.578	0.12
2	Bubbler – 4	305	140	0.459	
3	Filter – 13	167	368	2.20	-0.02
3	Bubbler – 11	239	530	2.22	
4	Filter – 11	157	257	1.64	-0.10
4	Bubbler – 13	270	470	1.74	
5	Filter – 10	152	41.4	0.272	0.14
5	Bubbler – 12	313	43.0	0.137	
6	Filter – 4	218	460	2.11 ^c	
6	Bubbler – 2	361	1700	4.71	
7 ^d	Filter – 2	185	11.5	0.062	
7 ^d	Bubbler – 3	342	4.8	0.014	
8 ^d	Filter – 12	158	ND ^e	ND ^e	
8 ^d	Bubbler - 14	293	1.6	0.006	

^aOSHA Method ID-101^bDifferences not statistically significant at 95% confidence level^cProbable breakthrough of silver filter^dBackground levels in plant^e None detected, detection limit = 0.6 µg/sample

Pairs 1 through 6 were area samples collected at the production site. Pairs 7 and 8 were collected in the plant but away from the chlorination operation to ascertain background levels of Cl₂. In pair 6, possible breakthrough of the silver filter occurred, as evidenced by the large mass of Cl₂ found in the bubbler. The mean difference for pairs 1 to 5 was 0.048 mg/m³. When a paired t-test for determining differences in paired observations [19] was applied to these data, the results showed that the difference between the two monitoring methods was not statistically significant at the 95% confidence level. However, since the number of samples taken at this field site was limited, more field studies must be done to verify the soundness of the significance statement.

References

1. Cassinelli ME [1991]. Development of a solid sorbent monitoring method for chlorine and bromine in air with determination by ion chromatography. *Appl Occ Environ Hyg* 6:215-226.
2. Moeller T [1952]. *Inorganic Chemistry*. John Wiley & Sons, New York.
3. Cotton FA, Wilkinson G [1980]. *Advanced Inorganic Chemistry*, 4th ed. New York, NY: John Wiley & Sons.
4. Glasstone ES [1946]. *Textbook of Physical Chemistry*, 2nd ed. New York, NY: D. Van Nostrand Company, Inc.
5. Hi-Q Environmental Products Co [1979]. 27 Catalog K. La Jolla, CA. La Jolla Scientific Co, Inc.
6. Mercer TT [1979]. Adsorption of mercury vapor by gold and silver. *Anal Chem* 51:1026.

7. Long SJ, Scott DR, Thompson RJ [1973]. Atomic absorption determination of elemental mercury collected from ambient air on silver wool. *Anal Chem* 45(13):2227.
8. NIOSH [1984]. Crystalline Silica by XRD: Method 7500. In: Eller PM, Cassinelli ME, eds. NIOSH manual of analytical methods, 3rd ed. Cincinnati, OH: Department of Health and Human Services, Centers for Disease Control, National Institute for Occupational Safety and Health, DHHS (NIOSH) Publication No. 84-100.
9. NIOSH [1984]. Development and Evaluation of Methods Chapter. In: Eller PM, Cassinelli ME, eds. NIOSH manual of analytical methods, 3rd ed. Cincinnati, OH: Department of Health and Human Services, Centers for Disease Control, National Institute for Occupational Safety and Health, DHHS (NIOSH) Publication No. 84-100.
10. NIOSH [1981]. By Busch KA, Taylor DG. Statistical Protocol for the NIOSH Validation Tests. In: Chemical Hazards in the Workplace -- Measurement and Control. Choudhary, G., Ed., ACS Symposium Series 149, pp 503-517, American Chemical Society, Washington, DC.
11. OSHA [1982] Bromine in Workplace Atmospheres, OSHA ID-108. Salt Lake City, UT: Department of Labor, Occupational Safety and Health Administration (OSHA), <https://www.osha.gov/dts/sltc/methods/inorganic/id108/id108.html>.
12. OSHA [1982]. Chlorine in Workplace Atmospheres, OSHA ID-101. Salt Lake City, UT: Department of Labor, Occupational Safety and Health Administration (OSHA), <https://www.osha.gov/dts/sltc/methods/inorganic/id101/id101.html>.
13. Lodge JP, Pate JB, Huitt HA [1963] The use of impregnated filters to collect traces of gases in the atmosphere, I. Suitability of membrane filters. *Am Ind Hyg Assoc J.* 24:380.
14. Chemical Rubber Co. [1970]. Handbook of Chemistry and Physics, 51st edition. Weast RC, ed. Cleveland, OH: Chemical Rubber Company.
15. Dillon HK, Fowler WK [1983]. Methods development for sampling and analysis of Chlorine, Chlorine Dioxide, Bromine, and Iodine. Southern Research Institute. Prepared for NIOSH under Contract No. 210 -80-0067. NTIS No. PB83-245282.
16. NIOSH [2014] Method 7907: Volatile acids by Ion Chromatography. In: Andrews R, O'Connor P, eds. NIOSH manual of analytical methods, 5th ed. Cincinnati, OH: Department of Health and Human Services, Centers for Disease Control and Prevention, National Institute for Occupational Safety and Health, DHHS (NIOSH) Publication No. 2014-151, www.cdc.gov/niosh/nmam.
17. NIOSH [1984]. Inorganic Acids: Method 7903. In: Eller PM, Cassinelli ME, eds. NIOSH manual of analytical methods, 3rd ed. Cincinnati, OH: Department of Health and Human Services, Centers for Disease Control, National Institute for Occupational Safety and Health, DHHS (NIOSH) Publication No. 84-100.
18. Kotrly S, Sucha L [1985]. Handbook of chemical equilibria in analytical chemistry. Ellis Horwood, Chichester, UK.
19. Bethea RM, Duran SS, Boullion TL [1985]. Statistical methods for engineers and scientists, 2nd edition. New York, NY: Marcel Dekker.