

Warning:

1. This product does not eliminate the risk of contracting any disease or infection.
2. Failure to follow all instructions and limitations on the use of this product could reduce the effectiveness of the respirator and result in sickness or death.
3. A properly selected respirator is essential to protect your health. Before using this respirator, consult an Industrial Hygienist, Occupational Safety Professional or Infection Control Practitioner to determine suitability for your intended use.
4. This product does not supply oxygen. Use only in adequately ventilated areas containing sufficient oxygen to support life. Do not use this respirator when oxygen concentration is less than 19.5%.
5. Do not use when concentrations of contaminants are immediately dangerous to health or life.
6. Leave work area and return to fresh air if (a) breathing becomes difficult or (b) dizziness or other distress occurs.
7. Facial hairs or beards, as well as certain facial characteristics may reduce the effectiveness of this respirator.
8. Never alter or modify the respirator in any way.
9. Only for single use. No maintenance is necessary. Discard the used respirator after single use.
10. Keep respirators in the display box and away from direct sunlight until use.

Fitting Instruction:

1. Hold the respirator in your hand with its nosepiece touching your fingertips, allowing the head straps to hang loosely.
2. Cover your mouth and nose with the respirator, while retaining straps around back of hand. Pull lower portion of the strap over head and fit it around your neck.
3. Stretch and pull the upper strap over your head, to the back of your head, above your ears.
4. Using both hands, form the nosepiece around your nose. Adjust the tension in the upper and lower portions of the continuous loop head strap to obtain a positive face seal. For a tighter facial fit, increase tension by pulling the straps at either side of the buckle. For a looser fit, decrease tension by pushing out the back of the buckle.
5. To test fit: Cup both hands over the respirator and exhale sharply. If air flows around your nose, tighten the nosepiece. If air flows around the edges, repeat steps 3 and 4.

Change respirator immediately if breathing becomes difficult, if respirator becomes damaged or distorted, or if a proper face fit can not be maintained. Careful observance of these instructions is an important step in safe respirator use.