## INSTRUCTIONS - IMPORTANT

Fitting instructions to be followed each time the respirator is worn:

1. Hold the respirator in your hand, allowing the headbands to hang freely below your hands.
2. Position the respirator comfortably under your chin and against the bridge of your nose.
3. Pull the top headband over your head, placing it at the crown of your head. Pull the bottom headband over your head and place it below your ears.
4. USING BOTH HANDS, mold the metal nosepiece to the shape of your nose and cheek area. Using one hand to pinch the nosepiece may cause a bad fit, resulting in lower effective respirator performance.
5. FIT CHECK: Place both hands completely over the respirator. Inhale sharply. Be careful not to disturb the position of the respirator. A negative pressure (mask wanting to collapse) should be felt. If not, or if air leaks around the edges of the respirator, reposition the respirator and headbands for a better fit. If unable to adequately seal filter with hands, obtain the optional respirator fit check cup to perform this test.

If you CANNOT achieve a proper fit, DO NOT enter the contaminated area. See your supervisor or call 1-888-MSA-HOW2 (1-888-672-4692) for proper fitting respiratory protection.

### Applications

Particulate aerosols free of oils meeting the requirements described herein.

### Preparations for Use:

1. Appropriately selected respirator is essential to protect your health. Before using this product, consult your supervisor and review the instructions. If uncertain about the suitability of this respirator for your intended use, call MSA Safety Works at 1-888-MSA-HOW2 (1-888-672-4692).
2. Before use, the wearer must be trained by the employer in proper respirator use.
3. OSHA (29 CFR 1910.134) requires that each wearer of this respirator have a qualitative or quantitative fit test performed to determine the amount of protection it will provide.

### Respirator Use Limitations Do not exceed the following:

#### 1. MAXIMUM USE CONCENTRATION:

- a) 10 times the OSHA permissible exposure limit (PEL), or other governmental limitation, for the contaminant(s) present.
  - b) Immediately dangerous to life or health (IDLH) concentration for any contaminant present.
2. 8 hours of use (continuous or intermittent). [Service time can be extended by performing an evaluation in the specific workplace setting that demonstrates (a) that the extended use will not degrade the filter efficiency below 95%, or (b) that the total mass loading of the filter is less than 200 mg.]
  3. Any applicable regulatory or industry standard limitation or any limitation contained in the above warning.

### Replace the respirator when:

- Unable to pass a fit check.
- Breathing resistance is too high.
- Time use limitation is reached.

## ▲ WARNING

This respirator provides LIMITED protection. It may help reduce exposure to airborne biological agents, including avian (bird) flu virus, other types of influenza, SARS, or other bacterial or viral biological agents, but will not eliminate the risk of exposure, infection, illness, or death. • This respirator can filter these contaminants at the designated efficiency; however, the government has NOT established a safe level of exposure to these substances. Therefore, breathing a particle of any of these substances may be unsafe. • The Center for Disease Control and Prevention (CDC) recommends health-care workers exposed to persons under evaluation for avian influenza use a respirator at least as protective as a NIOSH approved N-95 respirator. They specify that respirators MUST be used in accordance with a complete respiratory protection program as required by OSHA, which includes proper selection, training, fit-testing, and fit-checking. Detailed information on a respiratory protection program is available by contacting OSHA or visiting [www.osha.gov](http://www.osha.gov). • Do NOT remove respirator in contaminated areas. The outer surface of the respirator MUST be treated as if it is contaminated at all times. A respirator will NOT prevent one from catching the flu in other ways such as by touching the mouth, nose, or eyes with contaminated hands or objects. Biological agents, such as the flu virus, can be transmitted when infected individuals cough or sneeze and spread virus particles through the air to exposed surfaces which are touched. Tight-fitting safety goggles, or a full-facepiece respirator, MUST be used to help prevent transmission of the disease. • The CDC recommends frequent hand washing and wearing gloves to help prevent transmission of the disease due to exposure to surfaces where contaminants may be present, and also immediately following removal of the respirator. • Do NOT reuse or share maintenance-free respirators. ALWAYS clean cartridge-style respirators before reuse in accordance with the instructions provided. • This respirator was NOT evaluated for use by children or by people with a medical condition that might prevent the use of a respirator. These individuals MUST consult a physician before using. • Failure to follow all warnings and instructions can result in serious personal injury or death.

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