

JFY 4155S, 4155ML (N95)

## Donning Procedures

1. Hold the respirator in one hand with the nosepiece at your fingertips, allowing the straps to hang below your hand.
2. Press the respirator against your face with the nosepiece on the bridge of your nose.
3. Place the top strap high on the back of your head. Move the bottom strap over your head and position it below your ears.
4. Adjust the strap tension by pulling the straps as shown in the Figure.
5. Test the fit. Cup both hands over the respirator and exhale vigorously. If air flows around your nose, tighten the nosepiece. If air leaks around the edges, reposition the straps for better fit.

If you CANNOT achieve a proper fit, DO NOT enter the contaminated area.  
See your supervisor for proper fitting respiratory protection.

