

着用方法



① 図のように鼻あてが前にくるようにゴムバンドが下にたれるように持ちます。

② マスクをアゴで包むように当てます。

③ 片手でマスクを押さえながら下のゴムを首のまわりにつけ、次に上のゴムを頭上後頭部につけて下さい。

④ 両手の指で鼻あてが鼻に密着するように軽く押します。

⑤ 両手でマスク全体をおおい、息を吐き空気のもれをチェックして密着のいい場所にマスクを合わせます。

注意：

マスクは着用の都度にしっかりとフィットしていることを確認しなければなりません。両手をマスクの上にと置き、息吐き出します。鼻または周辺から空気が漏れる場合は、マスクがしっかりとフィットするまで鼻当てとヘッドバンドを調整します。

Respirator Fitting Instructions



① Hold respirator in your hand as show, with the nosepiece at your fingertips and the headbands hanging freely below your hand.

② Place respirator under your chin with the nosepiece up. Pull the top headbands to top back of head.

③ Pull the bottom headband overhead and place around neck below ears.

④ Adjust headbands for comfortable fit. Push the nosepiece inward with fingertips to conform to the shape of your nose.

⑤ Check the facepiece-to-face seal by cupping both hands over the respirator and exhale vigorously. If air escapes around your nose, reform the metal strip. If air escapes around the edges, reposition the straps. Repeat as necessary until no air escapes, and you have obtained a satisfactory seal.

Caution :

The respirator should be fit checked for each wearing. Place both hands completely over the respirator and exhale. If air leaks around nose or the edges, adjust the nosepiece and / or headband until a good fit is achieved.