



Routine Microscopy Procedures

Student Laboratory Exercises


Introduction

After you have completed the Basic Microscopy eLearning course, it is strongly recommended that you complete the following laboratory exercises to transfer the didactic content of the course to experiential knowledge gained through hands-on laboratory exercises with your equipment in your laboratory.

Your supervisor/mentor should work with you to develop these laboratory skills as well as confirm that these exercises have been completed. The number and types of exercises you will complete will be at the discretion of your supervisor/mentor based on procedures followed within your laboratory.

Included in the laboratory exercises portion of this course are the objectives of the exercises as well as the prepared exercises.

After the laboratory exercises are completed and discussed with your supervisor/mentor, your supervisor/mentor should then follow-up the exercises with instruction related to your laboratory's specific procedures or guidelines.



This job aid is a component of the free, on-demand CDC training course "Routine Microscopy Procedures." Find the course at <https://www.cdc.gov/labtraining>.

Student Laboratory Exercise Objectives

Laboratory Exercise Objectives:

After completing the routine microscopy procedures laboratory exercises, you will be able to:

- Outline the steps of preparing a smear.
- Express the purpose of the Gram stain procedure.
- Identify the types of reagents used in the Gram stain procedure.
- Sequence the steps in the Gram stain procedure.
- Interpret the results seen in the bacterial cells, with the effects of the various reagents during the Gram stain procedure.
- Outline the potassium hydroxide (KOH) procedure and its uses.
- Identify how to prepare and interpret a wet mount.
- Identify the steps and results obtained in the India Ink procedure.
- Identify and resolve commonly encountered problems during routine microscopy procedures.

Note: Be sure to review the proper use of personal protective equipment (PPE) and laboratory equipment according to your laboratory's procedures and safety manual.

Preparing for Student Laboratory Exercises

Supply List

1. Personal protective equipment
2. Brightfield (Compound) microscope with 10X, ,20X, 40X, and 100X objectives
3. Immersion oil
4. Lens paper
5. Lens cleaning solution
6. Microscope slides, frosted-edge
7. Cover slips
8. Loops (sterile plastic or metal)
9. Sterile pipettes
10. Slide rack
11. Slide warmer or Bunsen burner (optional for heat fixing the smear)
12. Absorbent paper, such as bibulous paper
13. Agar plate containing isolated colonies
14. Specimen or sample containing Trichomonads
15. Specimen containing Yeast
16. Specimen containing Clue cells
17. Pencil or wax pencil
18. Biohazard waste container: used for personal protective equipment, alcohol swabs and lens paper.
19. Sharps container: For microscope slides if they will be discarded after the examination is completed.

Reagent List

1. Sterile saline or water
2. Methanol (optional for fixing the smear)
3. Crystal violet
4. Gram's iodine
5. Decolorizer
6. Safranin (or carbol fuschin)
7. Potassium Hydroxide (KOH)

Student Laboratory Exercise II: Gram Stain

Laboratory Exercise II Objective

After completing this laboratory exercise, the participant will be able to:

- Utilize the Gram stain procedure to correctly perform a Gram stain.

Laboratory Exercise II: Performing a Gram Stain

1. Use the fixed smear from Exercise I to perform the Gram stain.
2. Place the prepared fixed smear on a slide rack then flood the slide with crystal violet.
3. Wait _____ seconds then rinse the slide with water.
4. Flood the slide with Gram’s iodine.
5. After _____ seconds rinse the slide with water.
6. Apply the decolorizer to the slide.
7. Rinse the slide immediately with water.
8. Flood the slide with counterstain.
9. Wait _____ seconds then rinse the slide with water.
10. Blot the slide with absorbent paper. Be careful not to wipe the cells off the slide.
11. Allow the newly stained slide to air dry completely.
12. View the slide under oil using the oil immersion objective for a total magnification of 1000X.
13. Record your results.

Notes

Mentor/Supervisor /Date

Student Laboratory Exercise IV: KOH Procedure

Laboratory Exercise IV Objective

After completing this laboratory exercise, the participant will be able to:

- Perform the potassium hydroxide (KOH) procedure.

Laboratory Exercise IV: Performing a KOH Procedure

1. Place a drop or 10 μ m of a specimen on a labeled microscope slide.
2. Add a drop of KOH to the specimen on the slide.
3. Add a coverslip.
4. Observe using the 40X objective. Look for:
 - a. Yeast
 - i. Round or oval
 - ii. Non-motile
 - iii. Vary in size (2-6 μ m)
 - iv. Some show buds
 - b. Pseudohyphae
 - i. Filaments with rounded ends
 - ii. Vary in length (20-100 μ m)

Remember: You are able to use your KOH job aid for this exercise.

Notes

Mentor/Supervisor /Date