Understanding Your Elastomeric Respirator

Elastomeric Quarter Facepiece Respirator

- Covers your nose and mouth, and seals between your chin and lower lip
- Assigned Protection Factor (APF) = 5
- Can have an exhalation valve or no valve
- Protects against particles, gases, or vapors when equipped with the proper filter or cartridge*







EHMR with an exhalation valve



EHMR with an exhalation valve filter



EHMR with no exhalation valve

Elastomeric Half Mask Respirator (EHMR)

- Covers your nose and mouth, and seals under your chin
- APF = 10
- Available in multiple sizes to fit a diverse user population
- EHMR with an exhalation valve
 - Your inhaled breath enters through the filter or cartridge material and your *unfiltered* exhaled breath exits through the exhalation valve
 - Exhalation valve increases comfort by reducing breathing resistance, removing moisture, and reducing CO₂ concentrations inside the respirator
 - Protects against particles, gases, or vapors when equipped with the proper filter or cartridge

• EHMR with no exhalation valve


- Your inhaled breath enters through the filter material and your exhaled breath exits through the same filter
- Protects against particles only
- Cannot be used with gas or vapor cartridges

· EHMR with an exhalation valve filter


- Your inhaled breath enters through the filter material and your exhaled breath is filtered after it passes through the exhalation valve
- Protects against particles, gases, or vapors when equipped with the proper filter or cartridge*

Elastomeric Full Facepiece Respirator

- Covers your whole face, from roughly the hairline to below the chin, and provides eye protection
- APF = 50[†]
- Available in multiple sizes to fit a diverse user population
- · Has an exhalation valve
- Your inhaled breath enters through the filter or cartridge material and your unfiltered exhaled breath exits through the exhalation valve
- Protects against particles, gases, or vapors when equipped with the proper filter, cartridge, or canister



All illustrations by NIOSH



NIOSH Approved is a certification mark of the U.S. Department of Health and Human Services (HHS) registered in the United States and several international jurisdictions.

For more information, visit cdc.gov/niosh/npptl/.

