

Assessment of Non-NIOSH Approved Innovative FFRs

Organization: JInJJ LLC

Date Tested: 3.26.2021

Respirator Model(s): Bal-Blk-A (Balaclava-style Mask), Bnd-Blk-A (Bandana-style Mask), Loop-Blk (Ear Loop-style mask), Filters

Tests: Filtration with NaCl (modified version of STP-0059)

This test result is for a filtering facepiece mask not approved by NIOSH or by any other international standards (e.g., China-KN95, Korea KFN94, UK FFP2 & FFP3, Australia P2, P3, Japan DS/DL2 & DS/DL3). The test conducted for this report is for filtration efficiency and is not equivalent to the test for NIOSH respirator certification per NIOSH test procedure STP-0059.

This assessment provides the filtration performance of an innovative respirator using a modified version of NIOSH test procedure STP-0059. This assessment does not include an evaluation of the fit of these respirators. Use of these innovative respirators in occupational settings would be contingent upon approval by the appropriate regulatory agency.

Filtration Efficiency Results: The minimum and maximum filter efficiencies were as follows: Bal-Blk-A (38.6% and 49.7%), Bnd-Blk-A (35% and 39.3%), Loop-Blk (33.8% and 36.3%), and Filter (20.3% and 21.5%). None of the masks or filters measured more than 95% efficiency. See Table 1.



Figure 1: Masks evaluated. A) Bal-Blk-A, Balaclava-style Mask. B) Bnd-Blk-A, Bandana-style Mask. C) Loop-Blk, Ear Loop-style Mask. D) Filter only as provided for testing.



Figure 2: Laboratory Test Photo (TSI, Inc. 8130 filter tester). A) Bal-Blk-A, Balaclava-style Mask. B) Bnd-Blk-A, Bandana-style Mask. C) Loop-Blk, Ear Loop-style Mask. D) Filter only as provided for testing.

Table 1. Filter Efficiency Evaluation

Respirator Model	Sample #	Flow Rate (Lpm)	Initial Filter Resistance (mmH ₂ O)	Initial Percent Leakage (%) (1 min)	Final Percent Leakage (%) (5 min)	Percent Leakage Variance (%)	Maximum Percent Leakage (%)	Filter Efficiency (%)
Bal-Blk-A Min: 38.6% Max: 49.7%	1	85	2.5	60.000	60.300	-0.300	61.300	38.700
	2	85	2.5	46.500	49.400	-2.900	50.300	49.700
	3	85	2.9	55.600	59.000	-3.400	59.200	40.800
	4	85	3.0	58.600	61.400	-2.800	61.400	38.600
	5	85	2.5	51.900	53.500	-1.600	53.500	46.500
Bnd-Blk-A Min: 35.0% Max: 39.3%	1	85	2.5	63.000	64.700	-1.700	65.000	35.000
	2	85	2.5	63.100	63.900	-0.800	64.000	36.000
	3	85	2.5	61.700	63.000	-1.300	64.400	35.600
	4	85	2.2	59.400	60.900	-1.500	60.900	39.100
	5	85	2.5	60.300	60.100	0.200	60.700	39.300
Loop-Blk Min: 33.8% Max: 36.3%	1	85	3.0	63.000	65.700	-2.700	66.200	33.800
	2	85	2.8	63.500	63.900	-0.400	64.000	36.000
	3	85	2.9	65.600	65.200	0.400	66.100	33.900
	4	85	2.2	62.300	64.800	-2.500	64.800	35.200
	5	85	2.9	63.200	63.700	-0.500	63.700	36.300
Filters Min: 20.3% Max: 21.5%	1	85	2.0	77.600	79.000	-1.400	79.200	20.800
	2	85	2.4	76.800	78.600	-1.800	78.600	21.400
	3	85	1.6	76.200	78.500	-2.300	78.500	21.500
	4	85	1.7	77.000	78.800	-1.800	79.300	20.700
	5	85	2.2	79.600	79.700	-0.100	79.700	20.300

Notes:

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not necessarily meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- **BOLD** filter efficiencies < 95%.