Manufacturer: Nanjing Sunny Technology Co., Ltd. Model Tested: Disposable Face Mask Date Tested: June 8, 2020

These findings pertain to the Nanjing Sunny Technology Co., Ltd., Disposable Face Mask. The packaging for this product indicates that it meets GB2626-2006 (the Chinese standard for Respiratory Protective Equipment – Non-Powered Air-Purifying Particle Respirator).

Ten respirators were submitted for evaluation. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found <u>here</u>.

No certificate of approval was provided with the samples received; therefore, the authenticity of the claims cannot be validated.

The maximum and minimum filter efficiency was 52.50% and 25.80%, respectively. All ten respirators measured less than 95%.

While the above-listed product classification has similar performance requirements to NIOSH-approved devices, NIOSH does not have knowledge about the sustained manufacturer quality system and product quality control for these products. NIOSH also does not have knowledge about the product's handling and exposures after leaving its manufacturer's control.

In addition, this product is an ear loop design. Currently, there are no NIOSH-approved products with ear loops; NIOSH-approved N95s have head bands. Furthermore, limited assessment of ear loop designs, indicate difficulty achieving a proper fit. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirator's represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for <u>Crisis Capacity Strategies (during known</u> <u>shortages)</u>.

Evaluation of International Respirators

Test: Modified TEB-APR-STP-0059

Date Tested: June 8, 2020

Report Prepared: June 8, 2020

Manufacturer: Nanjing Sunny Technology Co., Ltd.

Item Tested: Disposable Face Mask

Country of Certification: China (GB2626-2006)

Filter	Flow Rate (Lpm)	Initial Filter Resistance (mmH ₂ O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency	
1	85	8.6	58.4	58.4	41.60	
2	85	6.7	68.7	69.0	31.00	
3	85	5.9	72.4	72.4	27.60	
4	85	7.3	67.6	67.6	32.40	
5	85	5.7	67.2	67.3	32.70	
6	85	7.3	74.2	74.2	25.80	
7	85	6.7	68.4	68.7	31.30	
8	85	7.9	63.7	64.3	35.70	
9	85	10.1	47.5	47.5	52.50	
10	85	8.0	70.2	70.2	29.80	
	Minimum Filter Eff	iciency: 25.80	Maximum Filter Efficiency: 52.50			

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.

NPPTL National Personal Protective Technology Laboratory

Pictures have been added to the end of this report.

		15						
I.E	Item 产品名称							
1.55	Color 颜色	White 白色	Size 规格	140*156mm	1911			
	Standard 执行标准	GB 2626-2006						
E	Material 主要材料	Non-woven and Meltbl 无纺布和熔喷布						
	Pro. Date 生产日期 -	2020-04-22	Shelf Life 保质期	2年 2 years	15-0			
	Batch No 生产批号	A042220200423	Qty 数 量	10	12/18			
	QC No 质检员	01	Check result 检验结果	Qualified 合格				
	Storage and Transport 储存和运输	1.Temperature 0-40℃ 贮存环境,温度0-40℃ 2.Sun and rain shoul 运输过程中应防止日晒						
	委托方	Nanjing Sunny Techno 南京尚益科技有限公司						
10-	MFR H T C	温州鑫博医疗器械有限公司 斯江省温州市平阳县万全镇郑宋大道668-688号 电话0577-63782888						
	sk							
-1. 7	<u>非医用通用口罩</u> MADE IN CHINA 中国制造							
			1.	- 12				









