Manufacturer: CSIC Pacli (Nanjng) Tech Co. Ltd. Model Tested: Filter Protective Mask KN95 Date Tested: April 21, 2020

These findings pertain to the CSIC Pacli (Nanjng) Tech Co. Ltd., Filter Protective Mask KN95. The packaging indicates that this product meets GB2626-2006, the Chinese standard for Respiratory Protective Equipment – Non-Powered Air-Purifying Particle Respirator.

Ten respirators were submitted for evaluation. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found <u>here</u>.

No certificate of approval was provided with the samples received; therefore, the authenticity of the claims cannot be validated.

The maximum and minimum filter efficiency was 97.72% and 96.74%, respectively. All ten respirators measured more than 95%.

While the above-listed product classification has similar performance requirements to NIOSH-approved devices, NIOSH does not have knowledge about the sustained manufacturer quality system and product quality control for these products. NIOSH also does not have knowledge about the product's handling and exposures after leaving its manufacturer's control.

In addition, this product is an ear loop design. Currently, there are no NIOSH-approved products with ear loops; NIOSH-approved N95s have head bands. Furthermore, limited assessment of ear loop designs, indicate difficulty achieving a proper fit. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirator's represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for <u>Crisis Capacity Strategies (during known</u> <u>shortages)</u>.

### **Evaluation of International Respirators**

Test: Modified TEB-APR-STP-0059

Date Tested: April 21, 2020

Report Prepared: April 21, 2020

Manufacturer: CSIC Pacli (Nanjng) Tech Co. Ltd.

Item Tested: Filtered Protective Mask KN95

Country of Certification: China (GB2626-2006)

Filter	Flow Rate (Lpm)	Initial Filter Resistance (mmH <sub>2</sub> O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency	
1	85	7.7	2.71	2.79	97.21	
2	85	8.7	2.28	2.30	97.70	
3	85	8.0	2.95	2.95	97.05	
4	85	7.8	2.60	2.62	97.38	
5	85	8.7	2.28	2.28	97.72	
6	85	8.5	2.57	2.57	97.43	
7	85	7.8	2.60	2.60	97.40	
8	85	8.2	2.84	2.84	97.16	
9	85	8.0	2.73	2.73	97.27	
10	85	8.0	3.26	3.26	96.74	
Minimum Filter Efficiency: 96.74 Maximum Filter Efficiency: 97.72						

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.



Pictures have been added to the end of this report.



item 产品名称	Di		IED CERT 格证		7	
Color 颜色	W	小IEWがF日早 hite 自色	Size 规格	140*156mm		
Standard 执行标准	G	B 2626-2006 呼吸防热	- Andrew Barris	防颗粒物呼吸器		
Standard 执行标准		Non-woven and Meltblown Filter 无纺布和熔喷布				
Pro. Da 生产日第	te	2020. 03. 30	Shelf Life 保质期	2年 2 years		
Batch 生产批		20200330	Qty 数量	50		
QC No 质检		02	Check result 检验结果	Qualified 含格		
Storage a Transport 條存和运 WFR 生产厂 ADD 地址		d 1. Temperature 0-40°C, RH≤80% 此下于"""。 注意 2. 10 - 40°C, 相对 湿度 ≤ 80%				
		OSIC Pacil (Naning) Tech Co. Ltd				
		No 32 Changqing St. Shuige Rd. Nanjng. Jiangsu 江苏省南京市江宁经济技术开发区长青街32号				
		NADE IN	CHINA			







