#### NPPTL COVID-19 Response: International Respirator Assessment

Manufacturer: Lanshan Shendun Technology Co., LTD

Model Tested: SD-KN95 Date Tested: August 7, 2020

These findings pertain to the Lanshan Shendun Technology Co., LTD, Model SD-KN95. The packaging for this product indicates that it meets GB2626-2006 (the Chinese standard for Respiratory Protective Equipment – Non-Powered Air-Purifying Particle Respirator).

Ten respirators were submitted for evaluation. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found here.

No certificate of approval was provided with the samples received; therefore, the authenticity of the claims cannot be validated.

The maximum and minimum filter efficiency was 99.87% and 99.37%, respectively. No respirators measured less than 95%.

While the above-listed product classification has similar performance requirements to NIOSH-approved devices, NIOSH does not have knowledge about the sustained manufacturer quality system and product quality control for these products. NIOSH also does not have knowledge about the product's handling and exposures after leaving its manufacturer's control.

In addition, this product is an ear loop design. Currently, there are no NIOSH-approved products with ear loops; NIOSH-approved N95s have head bands. Furthermore, limited assessment of ear loop designs, indicate difficulty achieving a proper fit. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirators represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for <u>Crisis Capacity Strategies (during known shortages)</u>.

### **Evaluation of International Respirators**



Test: Modified TEB-APR-STP-0059

Date Tested: August 7, 2020

Report Prepared: August 7, 2020

Manufacturer: Lanshan Shendun Technology Co., LTD

Item Tested: SD-KN95

Country of Certification: China (GB2626-2006)

Pictures have been added to the end of this report.

Filter	Flow Rate (LPM)	Initial Filter Resistance (mmH₂O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency (%)
1	85	17.1	0.36	0.36	99.64
2	85	17.1	0.57	0.57	99.43
3	85	13.7	0.17	0.17	99.83
4	85	15.8	0.13	0.13	99.87
5	85	22.1	0.26	0.26	99.74
6	85	16.7	0.26	0.26	99.74
7	85	13.9	0.63	0.63	99.37
8	85	12.9	0.47	0.47	99.53
9	85	15.6	0.26	0.26	99.74
10	85	19.4	0.19	0.19	99.81
Minimum Filter Efficiency: 99.37%			Maximum Filter Efficiency: 99.87%		

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of
  respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and
  therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.

NPPTL COVID-19 Response: International Respirator Assessment



NPPTL COVID-19 Response: International Respirator Assessment





## MASK (KN95 NON-MEDICAL)

**KN95 Protection Level** 

#### Instructions

- 1. The user shall receive necessary training before use and wear it in the
- 2. This product can only be used in a ventilated place with sufficient
- 3. If the mask is blocked and has difficulty breathing or is damaged, please discard the mask or replace it.
- If breathing is difficult or unpleasant while wearing a mask, leave the workplace immediately and breathe fresh air.

#### Important Information

- air pollution is serious or when sleeping.
- 2. Use with caution on patients with asthma and other respiratory
- 3. Please use the mask correctly and do not abuse it. 4. The storage temperature is-28°C-38°C, and the relative humidity is lower than 80%.

#### How to wear

- 1. Open the mask. Hold the ear strips while the nose clip is facing the
- 2. Position the bottom of the mask under your chin and the nose clip against bridge of your nose.
- 3. Pull the ear straps behind the ear.
- 4. Using both hands, mould metal nose piece comfortably to shape of nose. Pinch the nose piece using one hand may result in less effective respirator performance.



Name: Mask (Kn95 Non-medical)

Model: SD-KN95

Category: Disposable, no breathing valve

Standard: GB2626-2006 Quantity: 10pcs/bag

Shelf life: 1 years Production date: See the seal Tel:0746-2999789 Manufacturer: LANSHAN SHENDUN

Technology Co., Ltd

Address: No.5 leather industrial park, Lanshan county, Hunan province









# 合格证

# CERTIFICATE OF QUALITY

产品名称	口罩(KN95非医用)		
PRODUCT NAME	MASK (KN95 NON-MEDICAL)		
产品型号 PRODUCT MODEL	SD-KN95		
产品规格 PRODUCT SPECIFICATION	10.5CM*15.5CM		
材质	熔喷布28% 无纺布72%		
TEXTURE OF MATERIAL	28%MELT BLOWN CLOTH 72%NON-WOVEN		
执行标准 STANDARD	GDZ020 2000		
制造商	蓝山县神盾科技有限公司		
MANUFACTURER	LANSHAN SHENDUN TECHNOLOGY CO., LTD		
地址	湖南省蓝山县湘江源皮具产业园5号		
ADDRESS	NO.5 LEATHER INDUSTRIAL PARK, LANSHAN COUNTY, HUNAN PROVINCE		
数量	10个/包		
QUANTITY	10PCS/BAG		
保质期	12个月		
SHELF LIFE	Parmonths		
联系电话 TEL	0746-2999789 PASS		
生产批号 PRODUCTION BATCH NUMBER	见外包装 03 PRINTED ON PACKAGE		
生产日期	见外包装		
THE DATE OF PRODUCTION	PRINTED ON PACKAGE		

## SD-KN95 mask manual (non-sterile type)

## Please read the instructions carefully

[Product Name] SD-KN95 (non-sterile type)

[Model and ifications] Ear-trap mask, takes 10.5cmx15.5cm as specification

[Product composition] The product is composed of a mask body, a nose bridge clip, and a mask band. The outer and inner layers are made of non-stick nonwoven fabric, and the filter layer is made of polypropylene melt-blown fabric. Non-sterile product supply. [Scope of application] Applicable to the health care of the wearer in the general environment without the risk of body fluids and

splashing.

[Performance structure]

- 1. The appearance of the mask should be neat and intact, and there should be no holes or stains on the surface.
- 2. After wearing the mask, it should cover the mouth, nose and jaw of the wearer. The mask should meet the design size. The maximum deviation should be no more than and less than 5%.
- 3. Nose clips must be provided on the mask. The nose clips can be made of bendable plastic material and the length should be no less than 8.5cm.
- 4. The mask band should be easy to wear and take, and the breaking strength at the connection point between each mask band and the mask body should be no less than 10N.
- 5. The particle filtration efficiency of the mask should be no less than 95%.

[Contraindications] None

[Precautions] 1. Please choose aseptic or non-sterile products according to the needs of your own use. 2. The longest use time of this product is 72H, and the filtration efficiency will gradually decrease with time going by. After use, it should be processed according to the requirements of hospitals and environmental protection departments.

3. Please check whether the packaging is intact before use, if the packaging or mask is damaged, do not use. Confirm the external packaging mark, production date, and validity period before using it, and use it within the validity period. 4. Use with caution if you are allergic to non-woven fabrics. 5. If there is a significant increase in respiratory resistance during use (when breathing is difficult), the mask is damaged, or the mask is contaminated (when contaminated with foreign substances such as blood stains and droplets), it should be replaced in time.

[Instructions] Open the package and take out the mask, check whet-her the mask is intact. The convex side faces outwards, the concave side faces inwards, and the side with the nose clip is facing upwards. There is a silicone adjustment buckle and auxiliary strip in the package. First fold the auxiliary strip through the ear strap rope and then insert the silicone adjustment buckle into the ear strap rope. Hang the ear strap to the ear, adjust the silicone adjustment buckle to the appropriate position (the silicone adjustment buckle adjusts the ear strap to within 17.5cm for the best air tightness). Gently press the nose clip to fit the nose bridge, and then press the nose clip to adjust the lower end of the mask to the jaw. Adjust as comfortably as possible, avoid touching the inside of the mask with your hands, place the mask on the face, and adjust it to the appropriate position.

[Storage conditions] Store it in a dry, ventilated, non-corrosive gas environment, away from fire and flammable materials.

[Transportation and storage] It should be protected from heavy pressure, direct sunlight and rain and snow.

[Packaging] 10 pcs/bag. For specific packaging specifications, please refer to product label or outer packaging.

[Production date] refer to the outer packaging

[Product validity period] 1 year

[Manufacturer] LANSHAN SHENDUN Technology Co., Ltd

[ Address] No.5 leather industrial park, Lanshan county, Hunan province

[Contact number] 0746-2999789

[Executive Standard] GB2626-2006

NPPTL COVID-19 Response: International Respirator Assessment





NPPTL COVID-19 Response: International Respirator Assessment







