Manufacturer: Dongguan Sure-on Industrial Ltd., (Makrite Industries, Inc.) Model Tested: 9500-N95; Lot# 041520 Date Tested: October 5, 2020

These findings pertain to the Dongguan Sure-on Industrial Ltd., (Makrite Industries, Inc.), model 9500-N95; Lot# 041520. The packaging and labeling indicate that it is a NIOSH-approved product under approval number TC-84A-5411. The packaging also indicates it meets GB2626-2006 (the Chinese standard for Respiratory Protective Equipment – Non-Powered Air-Purifying Particle Respirator).

Ten respirators were submitted for evaluation. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found <u>here</u>.

The maximum and minimum filter efficiency was 99.42% and 88.73%, respectively. Seven respirators measured more than 95%. Three respirators measured less than 95%.

This product has head bands/straps. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirators represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for <u>Crisis Capacity Strategies (during known</u> <u>shortages)</u>.

Evaluation of International Respirators

Test: Modified TEB-APR-STP-0059

Date Tested: October 5, 2020

Report Prepared: October 6, 2020

Manufacturer: Dongguan Sure-on Industrial Ltd., (Makrite Industries, Inc.)

Item Tested: 9500-N95; Lot# 041520

Country of Certification: China (GB2626), USA (claimed, 42 CFR 84)

Filter	Flow Rate (LPM)	Initial Filter Resistance (mmH ₂ O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency (%)
1	85	12.7	4.80	4.80	95.20
2	85	13.5	4.97	4.97	95.03
3	85	9.2	11.27	11.27	88.73
4	85	14.8	0.58	0.58	99.42
5	85	12.1	5.50	5.50	94.50
6	85	12.3	5.42	5.42	94.58
7	85	12.7	2.85	2.85	97.15
8	85	14.4	4.74	4.74	95.26
9	85	14.5	2.33	2.33	97.67
10	85	13.4	3.46	3.46	96.54
Minimum Filter Efficiency: 88.73% Maximum Filter					99.42%

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.



Pictures have been added to the end of this report.





LIMITED WARRANTY IMPORTANT NOTICE TO PURCHASER:

This limited warranty is made in lieu of the warranties of merchantability, fitness for particular purposes and all other warranties, expressed or implied. There are no other warranties which extend beyond the description on the face hereof. The physical standards and specifications of Makrite will be met by products sold.

Exclusive Remedies: damages for the breach of this limited warranty are limited to the replacement of such quantity of Makrite, products proved to be defectively manufactured. Except as provided above, Makrite shall not be liable or responsible for any loss, damage, or liability, direct, indirect, incidental, special, or consequential, arising out of sale, use, or misuse, or the inability to use products by the user.

WARNINGS:

i

- 1. Wearer must read fitting instructions and use limitations before use.
- 2. This respirator is not guaranteed to protect against the risk of contracting a disease or infection.
- 3. Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.

Made in China

RESPIRATOR FITTING INSTRUCTION



STEP 1

Hold the respirator in hand with the nosepiece at your fingertips, allowing the headbands to hang freely below your hand.



STEP 2

Press the respirator firmly against your face with the nosepiece on the bridge of your nose.



STEP 3

Stretch and position the top band high on the back of the head. Stretch the bottom band over the head and position below your ears.



STEP 4 Using both hands, mold nosepiece to the shape of your nose.



STEP 5


To test fit: a) Cup both hands over the respirator being careful not to distrub position, and b) inhale vigorously. If air leaks around the edges, reposition the straps or adjust strap tension for better fit.

Please carefully follow these fitting instructions during each use to achieve proper fit.

Made in China







Makrite Qualification 合 过滤口罩 Disposable Respirator 品 (非失用 号: 9500-N95 品尺寸: 144mm*177mm*53mm(生化 Product Size 品颜色: 白色 要成份: 56%涤纶短纤/18%不织布/16%像根 56%PE/18% Polyester/16%Koitted Elastic S 050920 Major Ingredient e Foam 生产日期: Production Date 041520 生产批号: Lot number 效期:三年 有 Period Of Validity 3 years 验员: Q003 检 ctor 生产企业: 东莞迅安塑胶纤维制品有限公司 Dongguan facturer 立地址: 广东省东莞市塘厦镇诸佛岭工业区 No.18, Guangyuan St., 168 Industrial Area, Zho Fo Ling, Tang Xia, Dongguan City, China Address 联系电话: 0769-87723170 亍标准: 美国NIOSH N95/GB 2626-2006 ecutive Standard 本工厂通过ISO 9001:2015国际质量管理体系认证 The factory has passed ISO 9001:2015 international quality management system



