Manufacturer: Dongguan Sure-On Plastic & Fiber Products Co., Ltd.., (Makrite Industries, Inc.) Model Tested: 9500-N95; Lot# 071720 Date Tested: October 7, 2020

These findings pertain to the Dongguan Sure-On Plastic & Fiber Products Co., Ltd.., (Makrite Industries, Inc.), model 9500-N95; Lot# 071720. The packaging and labeling indicate that it is a NIOSH-approved product, under approval number TC-84A-5411. The packaging also indicates that it meets GB2626-2006 (the Chinese standard for Respiratory Protective Equipment – Non-Powered Air-Purifying Particle Respirator).

Ten respirators were submitted for evaluation. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found <u>here</u>.

The maximum and minimum filter efficiency was 99.97% and 99.41%, respectively. All ten respirators measured more than 95%.

This product has head bands/straps. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirators represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for <u>Crisis Capacity Strategies (during known</u> <u>shortages)</u>.

## **Evaluation of International Respirators**

Test: Modified TEB-APR-STP-0059

Date Tested: October 7, 2020

Report Prepared: October 8, 2020

Manufacturer: Dongguan Sure-On Plastic & Fiber Products Co., Ltd.., (Makrite Industries, Inc.)

Item Tested: 9500-N95; Lot#071720

Country of Certification: China (GB2626), USA (claimed, 42 CFR 84)


Filter	Flow Rate (LPM)	Initial Filter Resistance (mmH <sub>2</sub> O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency (%)	
1	85	13.1	0.14	0.14	99.86	
2	85	16.4	0.04	0.04	99.96	
3	85	15.5	0.15	0.15	99.85	
4	85	15.2	0.11	0.11	99.89	
5	85	16.3	0.05	0.05	99.95	
6	85	14.4	0.17	0.17	99.83	
7	85	14.8	0.26	0.26	99.74	
8	85	15.2	0.59	0.59	99.41	
9	85	14.2	0.10	0.10	99.90	
10	85	15.5	0.03	0.03	99.97	
Ν	Ainimum Filter Effi	ciency: 99.41%	Maximum Filter Efficiency: 99.97%			

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.



Pictures have been added to the end of this report.





# **RESPIRATOR FITTING INSTRUCTION**



#### STEP 1

Hold the respirator in hand with the nosepiece at your fingertips, allowing the headbands to hang freely below your hand.



### STEP 2

Press the respirator firmly against your face with the nosepiece on the bridge of your nose.



#### **STEP 3**

Stretch and position the top band high on the back of the head. Stretch the bottom band over the head and position below your ears.



**STEP 4** Using both hands, mold nosepiece to the shape of your nose.




### **STEP 5**

To test fit: a) Cup both hands over the respirator being careful not to distrub position, and b) inhale vigorously. If air leaks around the edges, reposition the straps or adjust strap tension for better fit.

Please carefully follow these fiting instructions during each use to achieve proper fit.

Made in China





					1
THINKS SERVICE	11F-5 Xintai New	RITE INDUSTR 5, No.79, Sec. 5 5th Rd., XiZhi Taipei City 221 886-2-2698-24	1, Dist., , Taiwan	OSH	
	-1				
This respirator	is approved only in Protection <sup>1</sup>	Respirator	Cautions and Limitations <sup>2</sup>		
84A-5411	N95	9500-N95 x	ABCJMNOP		
	1. Prote	N95-Particula filter efficience	te Filter (95% y level)effective ulate aerosols		
		free of oil;time may apply	e use restrictions		
C - Do Not exce J - Failure to pro M - All approved with MSHA, N - Never subst configuration O - Refer to use maintenance	eed maximum use operly use and ma d respirators shall OSHA and other itute, modify, add, as specified by th	concentrations e aintain this produc be selected, fitte applicable regula or omit parts. Us he manufacturer. nd/or maintenanc s.	e only exact replacement par e manuals for information on	th. cordance ts in the	
SERVICE	TIME LI	AITATIO			
	e same respirato laces that could		N 8 hours, continuous or inte er loading.	ermittent,	
in dirty workpl 2. Leave contan proper fit can	laces that could ninated area and not be obtained (	result in high fil remove respirator or the respirator		ifficult or torted.	
in dirty workpl 2. Leave contan proper fit can Do not use th 3. The maximun	laces that could ninated area and not be obtained ne same respirat n service time of	result in high fill d remove respirator or the respirator tor again once it f this respirator i	ter loading. ator if breathing becomes d r is damaged, soiled, or dist	ifficult or torted. nated area.	
in dirty workpl 2. Leave contan proper fit can Do not use th 3. The maximun	laces that could ninated area and not be obtained ne same respirat n service time of r/s of non-use), b	result in high fill d remove respirato or the respirator tor again once it f this respirator i beginning from t	ter loading. ator if breathing becomes d r is damaged, soiled, or dist has been used in contamin s five (5) consecutive calen he first day of use.	ifficult or torted. nated area.	
in dirty workpl 2. Leave contan proper fit can Do not use th 3. The maximun (including day	laces that could ninated area and not be obtained ne same respirat n service time of r/s of non-use), b	result in high fill d remove respirato or the respirator tor again once it f this respirator i beginning from t	ter loading. ator if breathing becomes d is damaged, soiled, or dis has been used in contamir s five (5) consecutive calen he first day of use.	ifficult or torted. nated area.	

Makrite 合格证\_Qualification (非医用) (Non-medical) (一次性) Product Type (Single Use 产品型号: 9500-N95 Product Model 产品尺寸: 144mm\*117mm\*53mm(±10%) Product Size 产品颜色: 白色 Product Colour White 主要成份: 87.6% 涤纶短纤/12.4% 熔 Major Ingredient 87.6% PE/12.4% melt-blown fabric 07-生产日期: Production Date 生产批号: 07-Lot number 有效期:三年 Period Of Validity 3 years 检验员: 0003 Inspector 生产企业:东莞迅安塑胶纤维制品有限 公百 Manufacturer Dongguan Sure\_On Plastic & Fiber Products Co. LTD 生产地址:广东省东莞市塘厦镇诸佛岭工业区 Address ZhuFo Ling Industrial Area, Town, DongGuan City, GuangDong, China 联系电话: 86-0769-87723170 执行标准:美国NIOSH N95/GB 2626-2006 Executive Standard 本工厂通过ISO 9001:2015国际质量管理体系认证 The factory has passed ISO 9001 :2015 international quality management system







