Manufacturer: 3M Model Tested: 9542V

Date Tested: October 7, 2020

These findings pertain to the 3M, model 9542V. The packaging and labeling for this product indicate that it meets GB2626-2006 (the Chinese standard for Respiratory Protective Equipment – Non-Powered Air-Purifying Particle Respirator).

Ten respirators were submitted for evaluation. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found here.

No certificate of approval was provided with the samples received; therefore, the authenticity of the claims cannot be validated.

The maximum and minimum filter efficiency was 99.19% and 95.48%, respectively. All ten respirators measured more than 95%.

While the above-listed product classification has similar performance requirements to NIOSH-approved devices, NIOSH does not have knowledge about the sustained manufacturer quality system and product quality control for these products. NIOSH also does not have knowledge about the product's handling and exposures after leaving its manufacturer's control.

This product has head bands/straps. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirators represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for <u>Crisis Capacity Strategies</u> (during known <u>shortages</u>).

Evaluation of International Respirators



Test: Modified TEB-APR-STP-0059

Date Tested: October 7, 2020

Report Prepared: October 13, 2020

Manufacturer: 3M

Item Tested: 9542V


Country of Certification: China (GB2626-2006)

Pictures have been added to the end of this report.

Filter	Flow Rate (LPM)	Initial Filter Resistance (mmH ₂ O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency (%)
1	85	13.3	4.34	4.52	95.48
2	85	12.6	2.77	2.99	97.01
3	85	10.7	1.03	1.03	98.97
4	85	14.4	0.72	0.81	99.19
5	85	14.0	1.42	1.46	98.54
6	85	12.3	1.63	1.63	98.37
7	85	14.0	1.17	1.27	98.73
8	85	13.7	0.96	1.14	98.86
9	85	13.7	1.95	2.22	97.78
10	85	13.0	1.41	1.50	98.50
Minimum Filter Efficiency: 95.48%			Maximum Filter Efficiency: 99.19%		

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of
 respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and
 therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.















佩戴说明

- 製剤口罩差鼻夹的一面。使鼻夹位于口罩上方。用手杖住口罩 000 点 点 点口罩板住下吧。
- 2. 将上方头带校过头谈, 置于头顶上方。
- 3. 将下方头带拉过头顶, 置于颈后耳朵下方。
- 林双手手指置于金属扁夹中部。一边向内按压一边暗着两头向下向时间。5次、直 屋将扁夹完全按压成鼻梁形状为止。仅用单手捏口罩鼻夹两下向前向上向口罩的密 会性。
- 5. 在进入工作区域之前,使用者必须检查口泵与胎部的密合性。
- 4 用双手罩住口罩, 避免影响口罩在脸上的位置。
- 幼株理報气、如空气从鼻梁处泄棄。后按步骤4重新调整表示,如空气从口罩边缘泄漏。 应重新调整头带、如果无法取得良好密含、不要进入污染区域。请与您的主管联系。 以知识和感觉泄漏。可进入工作区工作。
- 牢御地按压鼻夹于鼻梁上以形成良好的密合是非常重要的。



ASI

本产品只用于对某些颗粒物的呼吸防护、错误使用可导致疾病甚至死亡。位了解正确使用方法、请参见包装上误明。请示您的主管或数电 第公司。400-925-1595。

FITTING INSTRUCTIONS

- Place the respirator against your face with the nosepiece across the bridge of your nose. Hold the respirator on your face and position the respirator under the chin.
- 2. Pull the top strap over your head restling it high at the top of your head.
- 3. Pull the bottom strap over your head position it around the neck, below the ear
- 4. Place the fingertips of both hands at the top of the metal resepiece. Mold the nosepiece to the shape of the rises bridge by pushing immands while moving your fingertips down both sides of the insequence. Placking the respirator rosepiece using only one hand may result in less effective respirator fit.
- The seal of the respirator on the face should be checked by the wearar prior to entering the work area.
- a) Cover the front of the respirator with both hands, being careful not to disturb the position of the
- to inhale sharply if an leaks around the rouse bridge, readjust the nonepiece as described in step 4 flar leaks at the respirator edges, work the strips back along the side of your horst. If you CANNOT achieve a proper seal, DO NOT writer the contaminating area. See your supervisor.
- c) If no leakage is detected, then work may proceed.
- It is very important to press the nosepiece firmly to the nose bridge to form a good seal.



A WARNIN

his respendir helps protect against certain particles. Misuse can result in acknows or death. For proper use, see supervisor or liter instructions on management that are soon process.

- 9.使用过程中未遵循使用说明和警告佩戴防护口罩,并且/或未能在暴露于污染物的整个过程中 全程佩戴防护口罩,将降低防护口罩的有效性,并可导致疾病,或对健康造成永久性伤害。
- 10. 务必避免外力推击、挤压、导致口罩变形、损坏、在抽车、跑步等运动过程中佩戴斯特别注意。 11. 本产品内含有的鼻夹和呼气阀等部件如至外力作用,可能导致人身伤害。

注意:当在低温下使用本产品时要小心,过多的温气可能导致呼气阀冻结。

储在各件和有效如

未挑封产品在保证储存温度在-20°C至+30°C范围内。相对温度使于80%的条件下,保存期为从 生产之日起三年。

本产品不含天然橡胶成分。

IMPORTANT

Before use, the weater must read and understand these User instructions. Keep these instructions for reference. See information supplied by the manufacturer.

Hee For

Particles such as those from grinding, sanding, sweeping, sawing, bagging or processing minurist, sifica, coal, iron ore, heavy metal, flour, wood, pollen and certain other substances. Liquid or non-oil baced particles from sprays that do not also emit oil aerosols or harmful vapors. Metal flumes produced from we'ding, birizing, cutting and other operations involving heating of metals. Use for nuisance level organic vapors. Nuisance level refers to concentrations below the occupational exposure limit.

Do Not Use For

Do not use for gases, sandblasting, oil aerosols, or for protection against organic vapors above the occupational exposure limit for these substances; particulate concentrations that exceed either 10 times the occupational exposure limit or applicable government regulations, whichever is lower. This respirator does not supply oxygen. Do not use in atmospheres containing less than 19.5% oxygen.

Do not use for particles in concentrations that exceed the maximum allowable concentration required by the applicable regulation and standard. The maximum allowable concentration is obtained by multiplying the Assigned Protection Factor (APF) provided by GB/T 18664 "Selection, Use and Maintenance of Respiratory Protective Eguipment" and the permissible Occupational Exposure Limit (OEL) of the contaminant.

Use Limitations and Warnings

- 1. The wearer must first be trained in the proper use and fit tested before using this respirator.
- 2. Do not use this respirator when concentrations of contaminants are immediately dangerous to life or health.
- 3. Leave the contaminated area immediately if breathing becomes difficult, or dizziness or other distress occurs.
- 4. Discard and replace respirator if it becomes damaged or breathing resistance becomes excessive.
- Inspect respirator before each use to ensure it is in good working condition. Examine all the respirator parts for signs of damage including the two strags, nosefoam, nosectlp and valve.
- 6. Do not after, wash, abuse or misuse this respirator.
- 7. Do not use with bears) or other facial hair or other conditions that prevent a good seal between the face and the sealing edge of the respirator.
- 8. Can be used up to the limits specified by applicable government regulations for this product.
- 9. Faiture to follow all instructions and warnings on the use of this respirator and/or faiture to wear this respirator sturing all times of exposure can reduce respirator effectiveness and may result in illness or permanent disability.
- 10. Make sure the respirator is not crushed or squeezed by an outside force, as this may result in the deformation and damage of the respirator. Especially pay more attention when the wearer is taking, running or conducting other activities.
- 11. Outside forces acting upon the noseclip or valve on the respirator may result in injury to the wearen

CAUTION: Care should be taken when using the equipment at low temperatures as excessive moisture may cause valves to freeze.

Storage Conditions and Shelf Life

Shelf life of unopened product is 3 years from date of manufacture when stored within temperature range of -20°C to +30°C and of less than 80% relative humidity.

This respirator contains no components made from natural rubber latex.







