Manufacturer: Shanghai Yunqing Industrial Co., Ltd.

Model Tested: YQD95

Date Tested: November 19, 2020

These findings pertain to the Shanghai Yunqing Industrial Co., Ltd., model YQD95. The packaging and labeling for this product indicate that it meets GB2626-2019 (the Chinese standard for Respiratory Protective Equipment – Non-Powered Air-Purifying Particle Respirator).

Thirty respirators were submitted for evaluation. The respirators were sampled into groups of ten for evaluation. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found <a href="https://examples.com/here">here</a>.

No certificate of approval was provided with the samples received; therefore, the authenticity of the claims cannot be validated.

The maximum and minimum filter efficiency was 99.48% and 98.86%, respectively. All thirty respirators measured more than 95%.

While the above-listed product classification has similar performance requirements to NIOSH-approved devices, NIOSH does not have knowledge about the sustained manufacturer quality system and product quality control for these products. NIOSH also does not have knowledge about the product's handling and exposures after leaving its manufacturer's control.

This product has head bands/straps. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirators represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for <u>Crisis Capacity Strategies</u> (during known <u>shortages</u>).

# **Evaluation of International Respirators**



**Test:** Modified TEB-APR-STP-0059

Date Tested: November 19, 2020

Report Prepared: November 20, 2020

Manufacturer: Shanghai Yunqing Industrial Co., Ltd.

Item Tested: YQD95 (Sample Group 1 of 3)

Country of Certification: China (GB2626-2019)

Pictures have been added to the end of this report.

Filter	Flow Rate (Lpm)	Initial Filter Resistance (mmH₂O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency
1	85	9.3	0.67	0.67	99.33
2	85	9.3	0.55	0.55	99.45
3	85	7.7	0.96	0.96	99.04
4	85	12.3	0.95	0.95	99.05
5	85	8.9	0.72	0.72	99.28
6	85	7.6	0.78	0.78	99.22
7	85	9.4	0.55	0.55	99.45
8	85	9.8	0.69	0.69	99.31
9	85	7.6	0.65	0.65	99.35
10	85	10.4	0.56	0.56	99.44
Minimum Filter Efficiency: 99.04			Maximum Filter Efficiency: 99.45		

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.

## **Evaluation of International Respirators**



Test: Modified TEB-APR-STP-0059

Date Tested: November 19, 2020

Report Prepared: November 20, 2020

Manufacturer: Shanghai Yunqing Industrial Co., Ltd.

Item Tested: YQD95 (Sample Group 2 of 3)

**Country of Certification:** China (GB2626-2019)

Filter	Flow Rate (Lpm)	Initial Filter Resistance (mmH₂O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency
11	85	9.5	0.52	0.52	99.48
12	85	7.5	0.94	0.94	99.06
13	85	10.1	0.65	0.65	99.35
14	85	8.6	0.68	0.68	99.32
15	85	8.2	0.68	0.68	99.32
16	85	9.1	0.79	0.79	99.21
17	85	9.5	0.53	0.53	99.47
18	85	8.3	0.83	0.83	99.17
19	85	9.6	0.96	0.96	99.04
20	85	10.0	0.79	0.79	99.21
	Minimum Filter Eff	iciency: 99.04	Maximum Filter Efficiency: 99.48		

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of
  respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore
  cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control
  over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.

### **Evaluation of International Respirators**



Test: Modified TEB-APR-STP-0059

Date Tested: November 19, 2020

Report Prepared: November 20, 2020

Manufacturer: Shanghai Yunqing Industrial Co., Ltd.


Item Tested: YQD95 (Sample Group 3 of 3)

Country of Certification: China (GB2626-2019)

Filter	Flow Rate (Lpm)	Initial Filter Resistance (mmH₂O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency
21	85	8.4	0.70	0.70	99.30
22	85	9.5	0.62	0.62	99.38
23	85	9.3	0.85	0.85	99.15
24	85	8.6	0.98	0.98	99.02
25	85	9.6	0.74	0.74	99.26
26	85	9.4	0.77	0.77	99.23
27	85	9.8	0.63	0.63	99.37
28	85	9.4	0.68	0.68	99.32
29	85	9.3	0.88	0.88	99.12
30	85	7.7	1.14	1.14	98.86
	Minimum Filter Eff	iciency: 98.86	Maximum Filter Efficiency: 99.38		

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of
  respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore
  cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control
  over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.







#### **USE INSTRUCTIONS**

- 1. Failure to follow all instructions and limitations concerming the use of the respirator and/or failure to wear the respirator for the whole duration of exposure to contaminants can seriously reduce the performance characteristics of the respirator and lead to illness, injury or death.
- 2. Before use, the wearer must first be trained by the employer in the correct use of the respirator in accordance with applicable safety and health standards.
- 3. Discard the respirator and replace with a new one if:
  - $\cdot$  Excessive clogging of the respirator causes breathing difficullty.
  - The respirator becomes damaged.
- 4. Leave the contaminated area if dizziness, irritation or other distress occurs.

#### **USE LIMITATIONS**

- 1. Do NOT use the respirator or enter or stay in a contaminated area under the following circumstances:
  - · Atmosphere contains less than 19.5% oxygen
  - Atmosphere contains oil aerosols if using N protection class respirator
  - · For protection against gases or vapors
  - · Contaminants or their concentrations are unknown or immediately dangerous to life or health
  - · Concentrations or contaminants exceed maximum use concentrations in applicable OSHA standards or applicable government regulations or 10 times the PEL (Permissible Exposure Limit), whichever is
- $\cdot \ \, \text{For sandblasting, paint-spray operations, as best os}$
- 2. Do NOT modify or misuse the respirator.
- 3. Do NOT use the respirator with beards or other facial hair that interferes with direct contact between the face and the edge of the respirator, or any other conditions that may prevent a good face-seal.

#### **STORAGE**

- 1. The respirator should be stored at the area which is well-ventilated, avoids direct sun exposure and dry area. Keep away from fire and pollution.
- 2. Storage temperature range is -20°C to 38°C Maximum relative humidity <80%.
- 3. Shelf life: 5 years
- 4. Date of production: Certificate of conformity.



Storage temperature Store humidity





Shelf life: 5 years Refer to manufacturer information

FOR MORE INFORMATION or assistance, contact your local sales representatives of YICHITA Safety Products or YICHITA authorized distributors.





















