Manufacturer: Advoque Safeguard Model Tested: ADV001 Date Tested: December 8, 2020

These findings pertain to the Advoque Safeguard, model ADV001. The packaging and labeling indicate that it is a NIOSH-approved product, under approval number TC-84A-PH02. However, this approval was revoked on September 9, 2020, and this product is no longer considered NIOSH-approved. Additional information regarding this revocation can be found <u>here</u>.

Ten respirators were submitted for evaluation. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found <u>here</u>.

The maximum and minimum filter efficiency was 92.19% and 88.50%, respectively. All ten respirators measured less than 95% filter efficiency.

This product has head bands/straps. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirators represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for <u>Crisis Capacity Strategies (during known</u> <u>shortages)</u>.

## **Evaluation of International Respirators**

Test: Modified TEB-APR-STP-0059

Date Tested: December 8, 2020

Report Prepared: December 9, 2020

Manufacturer: Advoque Safeguard

Item Tested: ADV001

### Country of Certification: USA (42 CFR 84)

**Initial Filter Initial Percent** Maximum Filter Efficiency Flow Rate Filter Resistance Leakage Percent Leakage (LPM) (%) (mmH<sub>2</sub>O) (%) (%) 85 7.57 8.00 92.00 1 6.9 2 85 7.1 8.49 8.74 91.26 3 85 6.3 7.74 7.81 92.19 4 85 7.4 9.57 10.0 90.00 5 85 7.7 10.3 10.3 89.70 6 85 6.2 9.86 10.1 89.87 7 85 7.7 10.6 10.9 89.10 8 85 7.9 11.5 11.5 88.50 9 85 5.8 11.4 11.4 88.60 10 85 7.0 10.8 10.9 89.10 Minimum Filter Efficiency: 88.50% Maximum Filter Efficiency: 92.19%

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.



Pictures have been added to the end of this report.










P/N ADVOO1 / TC-84A-PHO2 | NIOSH N95 © ADVOQUE SAFEGUARD LLC 2020. ALL RIGHTS RESERVED. www.advoquesafeguard.com







### 

Particles such as those from sanding, grinding, sweeping, sawing, bagging, or processing minerals, coal, iron ore, flour, metal, wood, pollen, and certain other substances. Liquid or non-oil based particles from sprays that do not also emit oil aerosols or vapors. For additional information on Advogue SafeGuard use recommendations for this class of respirator please call us at 1-408-560-2990

#### DO NOT USE PRODUCT FOR

Gases and vapors, oil aerosols, asbestos, arsenic, cadmium, lead, 4,4'-methylene dianiline (MDA) or sandblasting; concentrations that exceed 10 times the PEL/OEL, specific OSHA standards, or applicable government regulations, whichever is lower. This respirator does not supply oxygen.

#### BIOLOGICAL PARTICLES

This respirator can help reduce inhalation exposures to certain airborne biological particles (e.g. mold, Bacillus anthracis, Mycobacterium tuberculosis, etc.) but cannot eliminate the risk of contracting infection, illness, or disease. OSHA and other government agencies have not established safe exposure limits for these contaminants.

#### USE INSTRUCTIONS


- 1. Failure to follow all instructions and limitations on the use of this respirator and/or failure to wear this respirator during all times of exposure can reduce respirator effectiveness and may result in sickness or death.
- Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the requirements of OSHA 29 CFR 1910.134 such as training, fit testing, medical evaluation, and applicable OSHA substance-specific standards.
- 3. The particles which can be dangerous to your health include those so small that you cannot see them.
- 4. Leave the contaminated area immediately and contact the supervisor if dizziness, irritation, or other distress occurs.
- 5. When not in use, store the respirator away from contaminated areas.
- 6. Inspect respirator before each use to ensure that it is in good operating condition. Examine all the respirator parts for signs of damage including the two headbands, staples, nose-clip, and nose-foam. The respirator should be disposed of immediately upon observation of damage or missing parts. Filtering respirators are to be inspected prior to each use to assure there are no holes in the breathing zone other than the punctures around staples and no damage has occurred. Enlarged holes resulting from ripped or torn filter material around staple punctures are considered damage. Immediately replace the respirator if damaged. Staple perforations do not affect NIOSH approval.
- Dispose of used product in accordance with applicable regulations.

#### USE LIMITATIONS

- 1. This respirator does not supply oxygen. Do not use in atmospheres containing less than 19.5% oxygen.
- 2. Do not use when concentrations of contaminants are immediately dangerous to life and health, are unknown, or when concentrations exceed 10 times the permissible exposure limit (PEL), or according to specific OSHA standards or applicable government regulations, whichever is lower.
- 3. Do not alter, abuse or misuse this respirator.
- 4. Do not use with beards, other facial hair, or other conditions that prevent a good seal between the face and the sealing surface of the respirator.
- Respirators can help protect your lungs against certain airborne contaminants, however, they will not prevent entry through other routes such as the skin, which would require additional personal protective equipment (PPE).
- 6. This respirator is designed for workplace/occupational use by adults who are properly trained in their use and limitations. This respirator is not designed to be used by children.
- Individuals with a compromised respiratory system, such as asthma or emphysema, should consult a physician and must complete a medical evaluation prior to use.

#### TIME USE LIMITATION

If respirator becomes damaged, soiled, or breathing becomes difficult, leave the contaminated area immediately and replace the respirator.



PH. The limited FFR approval will be effective until the declaration that circumstances exist justifying the authorization of the emergency use of personal respiratory protective devices during the COVID-19 outbreak is terminated under Section 564(b)(2) of the Act, or the EUA is revoked under Section 564(g) of the Act, or until such time that the approved configuration can no longer be maintained.

This respirator contains no components made from natural rubber latex.

Advoque SafeGuard | N95 Particulate Filtering Respirator / Model ADV001 | User Instructions | 6/1/2020

