

## NPPTL COVID-19 Response: International Respirator Assessment

Manufacturer: LitezAll (Wanning)

Model Tested: Protective Mask

Date Tested: April 22, 2020

These findings pertain to the LitezAll (Wanning), Protective Mask. The packaging for this product indicates that it meets GB2626-2006, the Chinese standard for Respiratory Protective Equipment – Non-Powered Air-Purifying Particle Respirator.

Ten respirators were submitted for evaluation. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found [here](#).

A certificate of approval was not provided with the samples received; therefore, the authenticity of the claims cannot be validated.

The maximum and minimum filter efficiency was 98.14% and 43.00%, respectively. Seven respirators measured less than 95%. The remaining three respirators measured more than 95%.

While the above-listed product classification has similar performance requirements to NIOSH-approved devices, NIOSH does not have knowledge about the sustained manufacturer quality system and product quality control for these products. NIOSH also does not have knowledge about the product's handling and exposures after leaving its manufacturer's control.

In addition, this product is an ear loop design. Currently, there are no NIOSH-approved products with ear loops; NIOSH-approved N95s have head bands. Furthermore, limited assessment of ear loop designs, indicate difficulty achieving a proper fit. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

**This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process.** This assessment was developed as an assessment of the filter efficiency for those respirator's represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for [Crisis Capacity Strategies \(during known shortages\)](#).

**Evaluation of International Respirators**

**Test:** Modified TEB-APR-STP-0059

**Date Tested:** April 22, 2020

**Report Prepared:** April 23, 2020

**Manufacturer:** LitezAll

**Item Tested:** Wanning Protective Mask

**Country of Certification:** China (GB2626-2006)

Pictures have been added to the end of this report.

Filter	Flow Rate (Lpm)	Initial Filter Resistance (mmH <sub>2</sub> O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency
1	85	7.8	56.10	57.00	43.00
2	85	15.2	9.89	9.89	90.11
3	85	21.9	2.99	2.99	97.01
4	85	9.4	5.69	6.01	93.99
5	85	11.1	14.30	14.60	85.40
6	85	10.4	10.35	10.35	89.65
7	85	16.8	1.86	1.86	98.14
8	85	14.2	2.99	2.99	97.01
9	85	8.1	7.31	7.53	92.47
10	85	6.0	20.24	24.74	75.26
<b>Minimum Filter Efficiency: 43.00</b>			<b>Maximum Filter Efficiency: 98.14</b>		

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.



## WARNING


1. This product does not eliminate the risk of contracting any disease or infection.
2. Failure to follow all instructions and limitations on the use of this product could reduce the effectiveness of the product and result in illness.
3. This product does not supply oxygen. Use only in adequately ventilated areas containing sufficient oxygen to support life.
4. Do not use this product when oxygen concentration is less than 19.5%
5. Do not use when concentrations of contaminants are immediately dangerous to health or life.
6. Leave work area immediately and return to fresh air if (a) breathing becomes difficult, or (b) dizziness or other distress occurs.
7. Facial hairs or beards and certain facial characteristics may reduce the effectiveness of this product.
8. Never alter or modify this product in any way.
9. Only for Single use. No maintenance necessary. Discard the used product after single use.
10. Do not leave in direct sunlight.




Please refer to the manual before use and use a mask correctly!

**844-475-7344**

[www.premierproducts.com](http://www.premierproducts.com)  
Made in China



### How to Use



1. Open the mask to make the nose clip at the top, and pull the ear straps with both your hands.
2. Hold the mask against your chin to completely cover your nose and mouth.
3. Pull the ear straps behind your ears and adjust there to make you feel comfortable.
4. Use both your hands to adjust the shape of the ip. Place your fingers in the middle of the nose clip and press it inwards while moving your fingertips along both sides of the nose clip until it is pressed to fit the bridge of your nose. (Making the nose clip get a seal with only one hand may affect the tightness of the mask).
5. Cover the mask with your hand and exhale vigorously. If you feel the air escaping from the nose clip, it is required to tighten the nose clip; if the air escapes from the edge of the mask, readjust the headband to ensure tightness.

Manufactured for  
**LitezAll**  
Comments or Questions  
Toll Free 844-475-7344  
Peru, Illinois 61354  
[www.litezall.com](http://www.litezall.com)  
©2020 Premier Products  
Made in China

**WANNING**

### 合格证 Product certificate

品名: KN95口罩  
Product Name: KN95 face mask  
执行标准: GB2626-2006  
Standard  
生产日期: 2020/3/30  
Produce Date:  
质检员: 010  
Inspector:  
有效期: 2年  
Validity: 2 Years

Made in China

NPPTL COVID-19 Response: International Respirator Assessment



NPPTL COVID-19 Response: International Respirator Assessment

