Manufacturer: Dongguan AOXING AV Equipment Co., Ltd. Model Tested: AX-KF95 Date Tested: May 26, 2021

These findings pertain to the Dongguan AOXING AV Equipment Co., Ltd., model AX-KF95. The packaging for this product indicates that it meets GB2626-2006 (the Chinese standard for Respiratory Protective Equipment – Non-Powered Air-Purifying Particle Respirator).

Thirty respirators were submitted for evaluation. The samples were tested using the modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found <u>here</u>.

No certificate of approval was provided with the samples received; therefore, the authenticity of the claims cannot be validated.

The maximum and minimum filter efficiency was 89.75% and 78.87%, respectively. All thirty respirators measured less than 95% filter efficiency.

While the above-listed product classification has similar performance requirements to NIOSH-approved devices, NIOSH does not have knowledge about the sustained manufacturer quality system and product quality control for these products. NIOSH also does not have knowledge about the product's handling and exposures after leaving its manufacturer's control.

In addition, this product is an ear loop design. Currently, there are no NIOSH-approved products with ear loops; NIOSH-approved N95s have head bands. Furthermore, limited assessment of ear loop designs indicate difficulty achieving a proper fit. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirators represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for <u>Crisis Capacity Strategies (during known</u> <u>shortages)</u>.

Evaluation of International Respirators

Test: Modified TEB-APR-STP-0059

Date Tested: May 26, 2021

Report Prepared: May 26, 2021

Manufacturer: Dongguan AOXING AV Equipment Co., Ltd.

Item Tested: AX-KF95

Country of Certification: China (GB2626-2006)

Initial Filter Initial Percent Maximum Filter Efficiency Flow Rate Filter Resistance Leakage Percent Leakage (LPM) (%) (mmH₂O) (%) (%) 85 7.6 21.13 21.13 78.87 1 2 85 8.2 17.48 17.48 82.52 3 85 7.7 15.76 15.76 84.24 4 85 6.2 15.20 15.20 84.80 5 85 7.5 15.94 15.94 84.06 6 85 6.8 17.34 17.34 82.66 7 85 14.96 14.96 85.04 6.0 85 7.5 15.42 15.42 84.58 8 9 85 7.1 11.70 11.70 88.30 10 85 8.1 14.83 14.83 85.17 Minimum Filter Efficiency: 78.87% Maximum Filter Efficiency: 88.30%

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.



Pictures have been added to the end of this report.

Evaluation of International Respirators



Test: Modified TEB-APR-STP-0059

Date Tested: May 26, 2021

Report Prepared: May 26, 2021

Manufacturer: Dongguan AOXING AV Equipment Co., Ltd.

Item Tested: AX-KF95

Country of Certification: China (GB2626-2006)

Filter	Flow Rate (LPM)	Initial Filter Resistance (mmH ₂ O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency (%)
11	85	8.2	10.46	10.46	89.54
12	85	8.0	11.31	11.31	88.69
13	85	7.6	11.94	11.94	88.06
14	85	8.0	11.39	11.39	88.61
15	85	8.6	12.00	12.00	88.00
16	85	8.7	10.78	10.78	89.22
17	85	8.4	13.30	13.30	86.70
18	85	8.6	10.25	10.25	89.75
19	85	9.4	11.80	11.80	88.20
20	85	7.4	11.93	11.93	88.07
P	Minimum Filter Efficiency: 86.70% Maximum Filter Efficiency		m Filter Efficiency:	89.75%	

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.

Evaluation of International Respirators



Test: Modified TEB-APR-STP-0059

Date Tested: May 26, 2021

Report Prepared: May 26, 2021


Manufacturer: Dongguan AOXING AV Equipment Co., Ltd.


Item Tested: AX-KF95

Country of Certification: China (GB2626-2006)

Filter	Flow Rate (LPM)	Initial Filter Resistance (mmH₂O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency (%)
21	85	7.0	18.29	18.29	81.71
22	85	7.4	11.78	11.78	88.22
23	85	7.1	13.32	13.32	86.68
24	85	7.6	10.78	10.78	89.22
25	85	9.1	11.35	11.35	88.65
26	85	8.0	13.18	13.18	86.82
27	85	8.5	12.59	12.59	87.41
28	85	8.4	11.64	11.64	88.36
29	85	7.7	12.81	12.81	87.19
30	85	6.8	13.72	13.72	86.28
Π	Minimum Filter Efficiency: 81.71%		Maximum Filter Efficiency: 89.22%		

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.






LIMITED WARRANTY IMPORTANT NOTICE TO PURCHASER

This limited warranty is made in lieu of the warranties of merchantability, fitness for particular purposes and all other warranties, expressed or implied. There are no other warranties which extend beyond the descriptions on the face hereof. The products sold meet the physical standards and specifications of AOXING.

Exclusive Remedies: The damages for the breach of this limited warranty are limited to the replacement of the same quantities of AOXING under the condition that products MUST be proved to be defectively manufactured. Except as provided above, AOXING shall not be liable or responsible for any direct and indirect loss, arising out of sale, use, or misuse, or the inability to use products by the user.

KN95 Respirator

MADE IN CHINA







P	^亡 品合格证 roduct Certificate
产品名称	KN95防护口罩 (非医用)
Name of Product	KN95 Respirator Mask(Non-medical)
型 号	AX-KF95
Model Number	AX-KF95
品 牌	AOXING
Brand Make	AOXING
规 格	160x105mm (折叠)
Specification	160x105mm(Folded)
执行标准	GB2626-2006
Standard	GB2626-2006
主要成份	熔喷布17%、无纺布53%、静电棉30%
Composition	Welt-blown17%, Nonwovens53%, Electrostatic Cotton30%
数 量	10片
Quantity	10 Pieces
生产批次	见盒底
Production Batch	See the Bottom of the Box
生产日期 Production Date	2020 -04- 18
有效期 Validity	24个月 24 Months QC PASS
检验员代码	AX-QC-02
Test Code	AX-QC-02
DONGGUAN No.2 Tiansha	莞市澳星视听器材有限公司 AOXING AV EQUIPMENT CO.,LTD 东莞市塘厦镇田沙路2号 Road,Tangxia Town,Dongguan City, angdong Province,China 中国制造 Made in China





