







- 2.1. 5 ml Cell Lysis Buffer (CLB; 50 mM Tris:50 mM EDTA, pH 8.0 + 1% Sarcosyl per tube as prepared in the "Formulas of PFGE Reagents" section)  
e. g., 5 ml x 10 tubes = 50 ml
- 2.2. 25 µl Proteinase K stock solution (20 mg/ml) is needed per tube of CLB  
e. g., 25 µl x 10 tubes = 250 µl
3. Prepare the master mix by measuring the correct volume of Cell Lysis Buffer and Proteinase K into appropriate size test tube or flask and mix well.
4. Add 5 ml of **Cell Lysis Buffer/Proteinase K** master mix to each labeled 50 ml tube.
5. Trim excess agarose from top of plugs with scalpel, razor blade or similar instrument (optional). Open reusable plug mold and transfer plugs from mold with a spatula to appropriately labeled tube. If disposable plug molds are used, remove white tape from bottom of mold and push out plug(s) into appropriately labeled tube. **Be sure plugs are under buffer and not on side of tube.**

*The excess agarose, plug mold, spatula, etc. are contaminated. Dispose of or disinfect them appropriately.*

6. Remove tape from reusable mold. Place both sections of plug mold, spatulas, and scalpel in 90% ethanol, 70% isopropanol, 1% Lysol or other suitable disinfectant. **Soak them for at least 30 minutes before washing them.** Discard disposable plug molds appropriately.
7. Place tubes in non-Styrofoam rack and incubate in a 54-55°C shaker water bath or incubator **for 15-30 min with constant and vigorous agitation** (175-200 rpm). Plugs can be lysed for longer periods of time (up to overnight).
8. Pre-heat enough sterile CLRW and TE Buffer to 54-55°C so that plugs can be washed two times with 10-15 ml water (200-250 ml for 10 tubes) and washed four times with 10-15 ml TE (400-600 ml for 10 tubes).

## Washing of Agarose Plugs After Cell Lysis

*Plug washing can be started on Day 1 and finished the morning of Day 2 after overnight refrigeration of the plugs in TE (after completion of two water and two TE washes).*

*Most laboratories will find that their plugs are sufficiently stable to perform the following washing steps at 54-55°C. However, if plugs are nicked along the edges or breaking, lower the temperature to 50°C for the following washing steps. Washes at lower temperatures should be increased to up to 30 min.*

1. Remove tubes from water bath or incubator, and carefully pour off lysis buffer into an appropriate discard container; plugs can be held in tubes with a screened cap (Bio-Rad) or spatula.
2. Add at 10-15 ml sterile Ultrapure water (CLRW) that has been pre-heated to 54-55°C to each tube and shake the tubes in a 54-55°C water bath or incubator for 10-15 minutes.
3. Pour water down drain or treat subsequent washes according to laboratory's waste disposal guidelines and repeat wash step with pre-heated water one more time.
4. Pour off water, add 10-15 ml pre-heated (54-55°C) sterile TE Buffer, and shake the tubes in 54-55°C incubator or water bath for 10-15 min.
5. Pour off TE and repeat wash step with pre-heated TE three more times.
6. Decant last wash and add 5-10 ml sterile TE. Continue with " Restriction Digestion of DNA in Agarose Plugs " section or store plugs in TE at 4°C until needed. Plugs can be transferred to smaller tubes containing TE for long-term storage.

*If restriction digestion is to be done the same day, complete Steps 1-3 of next section (Restriction Digestion) during last TE wash step for optimal use of time.*

## Restriction Digestion of DNA in Agarose Plugs

A small slice of the plug (not the entire plug) should be digested with the primary restriction enzyme (*SfiI*) because less enzyme is required and other slices of the plug can be subjected to restriction analysis with other enzymes. *NotI* is recommended as the secondary enzyme for analysis of *Vibrio cholerae* and *Vibrio parahaemolyticus* isolates. The use of a secondary enzyme is useful in situations where the PFGE patterns obtained with the primary enzyme from two or more isolates are indistinguishable.

1. Label 1.5 ml microcentrifuge tubes with culture numbers; label 3 (10-well gel) or 4 (15-well gel) tubes for *Salmonella* ser. Braenderup H9812 standards<sup>2</sup>.

The number of *Salmonella* ser. Braenderup H9812 standards is dependent on the number of samples being run on the gel. Standards should be placed on the ends and between samples as needed; no more than 4 samples should be placed in between standard lanes.

2. **Pre-Restriction Incubation Step (highly recommended):** Prepare a master mix by diluting the appropriate 10X restriction buffer (Roche/Sigma, New England Biolabs or equivalent) 1:10 with sterile CLRW according to the following table:

Reagent	µl/Plug Slice	µl/7 Plug Slices	µl/11 Plug Slices
CLRW	180 µl	1260 µl	1980 µl
10X Restriction Buffer	20 µl	140 µl	220 µl
<b>Total Volume</b>	200 µl	1400 µl	2200 µl

The appropriate restriction buffer varies between vendors and may differ between enzymes from the same vendor. **Always** use the restriction buffer recommended by the vendor for the particular restriction enzyme.

3. Add 200µl diluted restriction buffer (1X) to labeled 1.5-ml microcentrifuge tubes.
4. Carefully remove plug from TE with spatula and place in a sterile disposable Petri dish or on large glass slide.
5. Cut a 2.0 to 2.5mm wide slice from each test samples and the appropriate number of *S. ser. Braenderup H9812* standards with a single edge razor blade (or scalpel, cover slip, etc.) and transfer to tube containing diluted restriction buffer. **Be sure plug slice is under buffer.** Replace rest of plug into the original tube that contains 5 ml TE buffer and store at 4°C.

*PulseNet* recommends that the combs with larger teeth (10 mm wide teeth) be used to cast the gels because computer analysis of the gel lanes is more accurate and less tedious than analysis of gel lanes cast with combs with the smaller teeth (<10 mm). Using combs with smaller teeth is not advised. The number of slices that can be cut from the plugs will depend on the skill and experience of the operator, integrity of the plug, and whether the slices are cut vertically or horizontally (plugs made in disposable molds).

- 5.1. Incubate sample and control plug slices in a water bath for 5-10 minutes or at room temperature for 10-15 minutes.
  - 5.1.1. Incubate samples to be digested with *SfiI* at 50°C.
  - 5.1.2. Incubate samples to be digested with *NotI* and *XbaI* at 37°C.
- 5.2. After incubation, remove buffer from plug slice using a pipet fitted with 200-250µl tip all the way to bottom of tube and aspirate buffer. Be careful not to damage the plug slice with pipet tip and that plug slice is not discarded with pipet tip.
6. Prepare the restriction enzyme master mix according to the following table. Keep vials of restriction enzyme on ice or in an insulated storage box (-20°C) at all times.

<sup>2</sup> Directions for producing PFGE plugs of *Salmonella* ser. Braenderup H9812 are in PNL05.

Enzymes may be purchased in several different stock concentrations. The *SfiI* stock enzyme should be ordered in concentrated form (40 U/ $\mu$ l) rather than the less concentrated form (10 U/ $\mu$ l). Either form is acceptable for *NotI* restriction. The calculations below are based on using an enzyme at a concentration of 40 U/ $\mu$ l. If a different concentration of enzyme is used, make necessary adjustments to the volume of enzyme and water to achieve a final concentration of 40 U/ sample.

Reagent	$\mu$ l/Plug Slice	$\mu$ l/7 Plug Slices	$\mu$ l/11 Plug Slices
CLRW	178 $\mu$ l	1246 $\mu$ l	1958 $\mu$ l
10X Restriction Buffer	20 $\mu$ l	140 $\mu$ l	220 $\mu$ l
BSA (20mg/ml)	1 $\mu$ l	7 $\mu$ l	11 $\mu$ l
<i>SfiI</i> or <i>NotI</i> (40U/ $\mu$ l)	1 $\mu$ l	7 $\mu$ l	11 $\mu$ l
<b>Total Volume</b>	200 $\mu$ l	1400 $\mu$ l	2200 $\mu$ l

**Addition of Bovine Serum Albumin (BSA; highly recommended):** Several restriction enzyme vendors specifically recommend the addition of 1X BSA to enzyme restriction mixtures while others do not. PulseNet recommends adding BSA to all enzyme restriction mixtures to minimize the incidence of incomplete restriction. BSA may be purchased in different stock concentrations. The calculations above are based on using BSA at a concentration of 20 mg/ml. If a different concentration of BSA is used, make necessary adjustments to the volume of BSA and water to achieve a final concentration of 0.1 mg/ml.

7. Add 200  $\mu$ l restriction enzyme master mix to each tube. Close tube and mix by tapping gently; be sure plug slices are under enzyme mixture.
8. Incubate sample and standard (control) plug slices for 4 hours in a water bath at the appropriate temperature for the enzyme.
  - 8.1. Incubate samples digested with *SfiI* at 50°C.
  - 8.2. Incubate samples digested with *NotI* and *XbaI* at 37°C.
9. If plug slices will be loaded into the wells the same day, continue with the next section (**CASTING AN AGAROSE GEL**) approximately 1 h before restriction digest reaction is finished so the gel can solidify for at least 30 minutes before loading the restricted PFGE plugs.

## Casting an Agarose Gel

1. Confirm that water bath is equilibrated to 55-60°C.
2. Make volume of 0.5X Tris-Borate EDTA Buffer (TBE) that is needed for both the gel and electrophoresis running buffer according to one of the following tables.

### 5X TBE:

Reagent	Volume in milliliters (ml)			
	(DR-II) Running buffer + 10 well gel prep	(DR-II) Running buffer + 15 well gel prep	(CHEF Mapper/DRIII) Running buffer + 10 well gel prep	(CHEF Mapper/DRIII) Running buffer + 15 well gel prep
5X TBE	210	215	230	235
CLRW	1890	1935	2070	2115
<b>Total Volume of 0.5X TBE</b>	2100	2150	2300	2350

**10X TBE:**

Reagent	Volume in milliliters (ml)			
	(DR-II) Running buffer + 10 well gel prep	(DR-II) Running buffer + 15 well gel prep	(CHEF Mapper/DRIII) Running buffer + 10 well gel prep	(CHEF Mapper/DRIII) Running buffer + 15 well gel prep
<b>10X TBE</b>	105	107.5	115	117.5
<b>CLRW</b>	1995	2042.5	2185	2232.5
<b>Total Volume of 0.5X TBE</b>	2100	2150	2300	2350

3. Make 1% SeaKem Gold (SKG; the only acceptable agarose for PulseNet PFGE) agarose in 0.5X TBE as follows:
  - 3.1. Weigh appropriate amount of SKG into 500 ml screw-cap flask.
  - 3.2. Add appropriate amount of 0.5X TBE as prepared in the tables above; swirl gently to disperse agarose.
    - 3.2.1. Mix 1.0 g agarose with 100 ml 0.5X TBE for 14cm-wide gel form (10 wells)
    - 3.2.2. Mix 1.5 g agarose with 150 ml 0.5X TBE for 21cm-wide gel form (15 wells)
4. Loosen cap or cover loosely with clear film and microwave for 60 seconds; mix gently and repeat for 15 second intervals until agarose is completely dissolved.
5. Recap flask and return to 55-60°C water bath and equilibrate the agarose in the water bath for 15 minutes or until ready to use.

SAFETY WARNING: USE HEAT-RESISTANT GLOVES WHEN HANDLING HOT FLASKS AFTER MICROWAVING.

---

*Loading Restricted Plug Slices on the Comb (option 1)*

1. Remove restricted plug slices from 50°C or 37°C water bath. Remove enzyme/buffer mixture and add 200 µl 0.5X TBE. Incubate at room temperature for 5 min. Alternatively, digested plug slices can be kept in the refrigerator for up to three days if they are stored in 0.5X TBE.
2. Assemble gel form, place on a leveling table and adjust until perfectly leveled. Place the comb holder so the front part (side with small metal screws) and teeth face the bottom of gel frame. Make sure the comb teeth touch the gel platform, when comb is upright.
3. Remove plug slices from tubes; lay comb flat on gel form and load plug slices on the bottom edge of the teeth, including *Salmonella* serotype Braenderup H9812 standards in outside lanes (teeth) and between samples as needed (no more than 4 samples between standards).
4. Load samples on remaining teeth of the comb and note locations.
5. Remove excess buffer with tissue or kimwipe. Allow plug slices to air dry on the comb for 5-10 minutes or seal them to the comb with 1% SKG agarose (55-60°C).
6. Position comb upright in leveled gel form and confirm that the plugs slices are correctly aligned on the bottom of the comb teeth and that the lower edge of the plug slice is flush against the black platform.
7. Carefully pour the agarose (cooled to 55-60°C) into the gel form and remove any bubbles or debris.
8. Put black gel frame in electrophoresis chamber. Add 2.2 L freshly prepared 0.5X TBE (2 L for CHEF-DR II). Close cover of unit.
9. Turn on power supply, pump calibrated to a flow rate of 1 L/minute (setting of ≈70) and cooling module (14°C) approximately 30 minutes before gel is to be run.
10. Remove comb after gel solidifies, about 30 minutes.
11. Fill in wells of gel with melted and cooled (55- 60°C) 1% SKG (optional). Unscrew and remove end gates from gel form; remove excess agarose from sides and bottom of casting platform with a tissue or kimwipe. Keep gel on casting platform and carefully place gel inside black gel frame in electrophoresis chamber. Close cover of chamber.

---

## LOADING RESTRICTED PLUG SLICES INTO THE WELLS (OPTION 2)

1. Remove restricted plug slices from water bath. Remove enzyme/buffer mixture and add 200  $\mu$ l 0.5X TBE. Incubate at room temperature for 5 minutes. Alternatively, digested plug slices can be kept in the refrigerator for up to three days if they are stored in 0.5X TBE.
2. Place the gel form on a leveling table and adjust until perfectly leveled before pouring gel. Position the comb holder so that front part (side with small metal screws) and teeth face the bottom of the gel and the **bottom of the comb is 2 mm above the surface of the gel platform.**
3. Cool melted SKG agarose in 55-60°C water bath for 15-20 min; carefully pour agarose into gel form (casting stand) fitted with comb. Be sure there are no bubbles.
4. Put black gel frame in electrophoresis chamber. Add 2.2 L freshly prepared 0.5X TBE (2 L for CHEF-DR II). Close cover of unit.
5. Turn on power supply, pump calibrated to a flow rate of 1 L/minute (setting at ~70), and cooling module (14°C) approximately 30 minutes before gel is to be run.
6. Remove comb after gel solidifies, about 30 minutes.
7. Remove plug slices from tubes with tapered end of spatula and load into appropriate wells, including *Salmonella* serotype Braenderup H9812 standards in outside lanes and between samples as needed (not more than 4 samples between standards). Gently push plugs to bottom and front of wells with wide end of spatula. Manipulate position with spatula and be sure that there are no bubbles.

*Loading the plug slices can be tedious; each person has to develop his/her own technique for consistently placing the plug slices in the wells so the lanes will be straight and the bands sharp.*

8. Fill in wells of gel with melted 1% SKG Agarose (equilibrated to 55-60°C). Allow to harden for 3-5 min. Unscrew and remove end gates from gel form; remove excess agarose from sides and bottom of casting platform with a tissue or kimwipe. Keep gel on casting platform and carefully place gel inside black gel frame in electrophoresis chamber. Close cover of chamber.

*Vibrio cholerae* strains restricted with *SfiI* or *NotI*:

- Select following conditions on **CHEF Mapper with a two block program**

**Note: Press 'Enter' after each value or command is entered.**

- Press the Multi-State button
- Program with Interrupts?
  - Press 0 for "no"
- Block 1 Runtime?
  - 13 hours
- Block 1, State 1:
  - Volts = 6.0 V
  - **Angle= 60°**
  - Initial Switch Time = 2 s
  - Final Switch Time = 10 s
  - Ramping Factor = 0 (linear)
  - Continue with another state (vector)?
    - 1 = Yes
- Block 1, State 2
  - Volts = 6.0 V
  - **Angle= - 60° (Note: The angle for State 2 is Negative)**
  - Initial Switch Time = 2 s
  - Final Switch Time = 10 s
  - Ramping Factor = 0 (linear)
  - Continue with another state (vector)?
    - 0 = No
- Block 2 Runtime?
  - 6 hours
- Block 2, State 1:
  - Volts = 6.0 V
  - **Angle= 60°**
  - Initial Switch Time = 20 s
  - Final Switch Time = 25 s
  - Ramping Factor = 0 (linear)
  - Continue with another state (vector)?
    - 1 = Yes
- Block 2, State 2
  - Volts = 6.0 V
  - **Angle= - 60° (Note: The angle for State 2 is Negative)**
  - Initial Switch Time = 20 s
  - Final Switch Time = 25 s
  - Ramping Factor = 0 (linear)
  - Continue with another state (vector)?
    - 0 = No
- A Program is in memory, please enter another command
- Press the **Start Run** button

- Select following conditions on CHEF-DR III
  - Block 1
    - Initial switch time: 2 s
    - Final switch time: 10 s
    - Voltage: 6 V
    - Included Angle: 120°
    - Run time: 13 hours
  - Block 2
    - Initial switch time: 20 s
    - Final switch time: 25 s
    - Voltage: 6 V
    - Included Angle: 120°
    - Run time: 6 hours

---

*Vibrio parahaemolyticus* strains restricted with *SfiI* or *NotI*:

- Select following conditions on CHEF Mapper
  - Auto Algorithm
  - 78 kb: low MW
  - 396 kb: high MW
  - Select default values except where noted by pressing "Enter."
  - Change run time to 17-20 hours (See note below)
  - (Default values: Initial switch time = 10 s; Final switch time = 35.03 s)
- Select following conditions on CHEF DR-III
  - Initial switch time: 10 s
  - Final switch time: 35 s
  - Voltage: 6 V
  - Included Angle: 120°
  - Run time: 17-20 hours (See note below)
- Select following conditions on CHEF DR-II.
  - Initial A time: 10s
  - Final A time: 35 s
  - Start Ratio: 1.0 (if applicable)
  - Voltage: 200 V
  - Run time: 17-20 hours (See note below)

*The electrophoresis running times recommended above are based on the equipment and reagents used at the CDC. The 21 cm wide (15-well) gels require ~1 hr longer than 14 cm wide (10-well) gels. Run times may be different in your laboratory and will have to be optimized for your gels so that the lowest band in the S. ser. Braenderup H9812 standard migrates 1.0-1.5 cm from the bottom of the gel. See Appendix PNL06-1 for Vibrio cholerae block times if 19 hr total time is too short or too long.*

*Make note of the initial milliamp (mA) reading on the instrument. The initial mA should be between 110-150 mA. A reading outside of this range may indicate that the 0.5X TBE buffer was prepared improperly and the buffer should be remade.*

## Staining and Documentation of an Agarose Gel

*The following staining procedure describes the use of ethidium bromide to stain PFGE gels. Alternate DNA stains may be used. Please see the QuickTip [“20140218 Staining”](#) within the Library of PulseNet Documents forum (QuickTips/Wet Lab/PFGE) on the SharePoint site for additional information.*

1. When electrophoresis run is over, turn off equipment; remove and place gel in an approximately 14 cm x 24 cm covered container with 400 ml ethidium bromide or other approved stain (see PNL02). Larger or smaller volumes may be used for different sized containers. Stain gel, gently rocking for 20-30 min in covered container.

*Ethidium bromide (EtBr) is toxic and a mutagen. Stock solutions of 10 mg/ml EtBr in water are available from several commercial companies (see PNL01). The diluted solution can be kept in the dark at room temperature and re-used up to 15 times within 3 weeks before discarding according to your institution's guidelines for hazardous waste; **do not pour down the drain**. Aqueous solutions containing EtBr can be treated using de-staining bags from Amresco (E732-25), which effectively and safely remove EtBr from solutions. Once the EtBr is removed, the treated aqueous solutions can be discarded down the drain. Refer to the Safety Data Sheets (SDS) provided by the vendor for more details.*

*Currently, the acceptable alternative stain options are GelRed™ (Biotium, 41002), SYBR® Safe (ThermoFisher Scientific S33102) and SYBR® Gold (ThermoFisher Scientific S11494). Labs are strongly encouraged to follow manufacturer's instructions and test stains in their labs before adopting them for routine use. If one of the alternative stains is used, the de-staining steps should be omitted. However, the gels can be briefly rinsed with CLRW before imaging. Diluted GelRed™ solution can be kept in the dark at room temperature and re-used up to 10 times within 10 days before discarding down the drain. Usage parameters have not been established for other alternative stains.*

2. De-stain gel in approximately 500 ml CLRW for 60-90 min, changing water every ~20 minutes (at least 3 times). Capture image on GelDoc XR, XR+ or equivalent documentation system. If background interferes with resolution, de-stain for an additional 30-60 min (2 or 3 more washes).
3. Follow directions from the imaging equipment to save gel image as an \*.1sc (QuantityOne software) or \*.scn (ImageLab software) file; convert this file to \*.tif file for analysis with BioNumerics software program. The gel image should fill the entire window of the imaging equipment (computer) screen (without cutting off wells or lower bands). Ensure that the image is in focus and that there is little to no saturation (over-exposure) in the bands (signified by red pixilation in the software). Additional instructions are provided in PNL07 (Image Acquisition) of the PulseNet QA/QC Manual.
4. Drain buffer from electrophoresis chamber and discard. Rinse chamber with ~2 L CLRW or, if unit is not going to be used for several days, flush lines with water by letting pump run for 5-10 min **before** draining water from chamber and tubing.

**USE OF TRADE NAMES AND COMMERCIAL SOURCES IS FOR IDENTIFICATION PURPOSES ONLY AND DOES NOT IMPLY ENDORSEMENT BY CDC OR THE U.S. DEPARTMENT OF HEALTH AND HUMAN SERVICES.**

## CLIA Laboratory Procedure Manual Requirements

*Efforts have been made to assure that the procedures described in this protocol have been written in accordance with the 1988 Clinical Laboratory Improvement Amendments (CLIA) requirements for a procedure manual (42 CFR 493.1211). However, due to the format required for training, the procedures will require some modifications and additions to customize them for your particular laboratory operation.*

*Any questions regarding the CLIA requirements for a procedure manual, quality control, quality assurance, etc., should be directed to the agency or accreditation organization responsible for performing your laboratory's CLIA inspection. In*

*addition, some states and accreditation organizations may have more stringent requirements that will need to be addressed.*

## Records Management

1. Establish a record book or electronic log (for an example, see Appendix PNL06-2: PFGE Isolates and Plug Run Log) for tracking isolates, including but not limited to, the following information:
  - 1.1. Sample number (state public health lab identifier or similar)
  - 1.2. Date the isolates were received in the PFGE lab
  - 1.3. Date the plugs were made
  - 1.4. Location (box and slot) where the plugs are kept
  - 1.5. Gel and lane numbers for each sample
2. Establish a record (for an example, see Appendix PNL06-3: PFGE Plug Preparation Worksheet) of plug preparation, including but not limited to, the following information:
  - 2.1. Organism and number of samples
  - 2.2. Date plugs were made and initials of who made them
  - 2.3. Lot number and expiration dates of reagents used
3. Establish a record (for an example, see Appendix PNL06-4: PFGE Enzyme Master Mix and Gel Setup Worksheet) for each gel, including but not limited to, the following information:
  - 3.1. Gel number and date gel run
  - 3.2. Electrophoresis equipment used
  - 3.3. Running conditions/times
  - 3.4. Lot number and expiration dates of reagents used
  - 3.5. Restriction temperatures and times
  - 3.6. Order of isolates on a gel

## Formulas of Selected Reagents used in PulseNet Standardized Laboratory Protocol for PFGE

Tris: EDTA Buffer, pH 8.0: (TE, 10 mM Tris: 1 mM EDTA, pH 8.0)

10 ml of 1 M Tris, pH 8.0

2 ml of 0.5 M EDTA, pH 8.0

Dilute to 1000 ml with sterile Ultrapure water (CLRW)

CELL SUSPENSION BUFFER (100 MM TRIS:100 MM EDTA, PH 8.0)

100 ml of 1 M Tris, pH 8.0

200 ml of 0.5 M EDTA, pH 8.0

Dilute to 1000 ml with sterile Ultrapure water (CLRW)

Cell Lysis Buffer: (50 mM Tris: 50 mM EDTA, pH 8.0 + 1% Sarcosine + 0.1 mg/ml Proteinase K)

50 ml of 1 M Tris, pH 8.0

100 ml of 0.5 M EDTA, pH 8.0

100 ml 10% N-Lauroylsarcosine, Sodium salt (Sarcosyl) OR 10 g of N-Lauroylsarcosine, Sodium salt (Sarcosyl).

Dilute to 1000 ml with Sterile Ultrapure water (CLRW)

*Add 25 µl Proteinase K stock solution (20 mg/ml) per 5 ml of cell lysis buffer just before use for a final concentration in the lysis buffer of 0.1 mg/ml Proteinase K.*

*If Sarcosyl powder is added directly to the other components of this reagent, warm the solution to 50- 60°C for 30-60 minutes, or leave at room temperature for about 2 hours to completely dissolve the Sarcosyl.*

## Contacts

CDC PulseNet  
MS-C03  
1600 Clifton Rd. NE  
Atlanta, GA 30329  
Phone (404) 639-4555  
[PulseNet@cdc.gov](mailto:PulseNet@cdc.gov)

## AMENDMENTS

1. The phrase "Type I Water" has been changed to "Ultrapure Clinical Laboratory Reagent Water (CLRW)." The water composition is the same, but this reflects a change in the terminology used by the Clinical Laboratory Standards Institute (CLSI).
2. March 2013 changes:
  - Protocol was revised to combine *Vibrio cholerae* and *parahaemolyticus* instructions.
  - Corrected formula for TE buffer. TE used at CDC is 10 mM for Tris and 1 mM for EDTA.
  - Recommended disinfectant changed from 10% bleach to 1% Lysol/Amphyll or 90% ethanol.
  - Corrected 1% SKG / TE plug agarose recipe.
  - A note was added to provide guidance when working with large numbers of isolates (>10).
  - Plug washing steps can be performed at 54-55°C rather than lowering to 50°C.
  - Volume of TE needed to wash 10 plugs was corrected from 300 – 350 ml to 400 – 600 ml.
  - A statement was added to clarify that using combs with small teeth (5.5 mm) was not advised.
  - References to specific restriction buffers have been removed.
  - Moved reference to BSA out of footnotes and into the main text of the protocol. Use of pre restriction step and BSA was changed from optional to highly recommended. Calculation for including BSA in restriction enzyme master mix was added.
  - The word "Sterile" has been deleted in reference to diluting 5X or 10X TBE to 0.5X TBE. Non-sterile CLRW is acceptable.
  - A statement was included to allow the use of an alternative agarose for casting the running gel, along with recommendations strongly urging each lab to optimize the run time. Internal and external validation showed that run times could be affected by agarose type, but no trends were noted so a blanket recommendation on run times cannot be made. Additional agarose alternatives may be tested and deemed acceptable at a later date.
  - Added parameters for programming a CHEF DR-II for *Vibrio parahaemolyticus*.
  - Added a recommendation for laboratories to monitor the initial mA when electrophoresis is started.
  - A statement was included to allow the use of alternative DNA stains that are equivalent to EtBr. Labs are strongly urged to follow manufacturer's instructions as well as test stains in their own labs to gain experience using alternative agarose stains. Additional stain alternatives may be tested and deemed acceptable at a later date.

- The option to allow incubation times for restriction digestion to be increased longer than recommended was deleted.
3. April 2017 changes:
- Reformatted in accordance with new document layout (renumbered outline, removed footer and changed header.
  - Added statement recommending BHI as a growth media for *Vibrio parahaemolyticus*.
  - Updated CDC Contact information.
  - Added "Approval Signatures" section.
  - Added "Records Management" section by merging PNL08.
  - Statement allowing use of Amresco LF agarose (Amresco) was deleted. Additional testing revealed run time and normalization were negatively impacted by this agarose.
  - Added Appendices PNL06-1 through PNL06-4.

## APPENDICES

### APPENDIX PNL06-1: INDIVIDUAL BLOCK RUN TIME CALCULATION FOR PROGRAMMING *VIBRIO CHOLERAE*

**Note:** Ratio of block 1: block 2 is 2.167 +/- 0.011. Block times may be entered as either hh:mm format or decimal format.

	<b>hh:mm</b>	<b>Decimal (hr)</b>
Block 1	11:38	11.63
Block 2	5:22	5.37
<b>Total time</b>	<b>17:00</b>	<b>17</b>
Block 1	12:15	12.25
Block 2	5:39	5.65
<b>Total time</b>	<b>17:54</b>	<b>17.9</b>
Block 1	12:29	12.48
Block 2	5:46	5.76
<b>Total time</b>	<b>18:15</b>	<b>18.24</b>
Block 1	13:00	13
Block 2	6:00	6
<b>Total time</b>	<b>19:00</b>	<b>19</b>
Block 1	13:23	13.39
Block 2	6:11	6.18
<b>Total time</b>	<b>19:34</b>	<b>19.57</b>
Block 1	13:36	13.6
Block 2	6:18	6.3
<b>Total time</b>	<b>19:54</b>	<b>19.9</b>
Block 1	14:18	14.3
Block 2	6:36	6.6
<b>Total time</b>	<b>20:54</b>	<b>20.9</b>

APPENDIX PNL06-2: PFGE ISOLATES AND PLUG RUN LOG





