



POTENTIAL AREAS TO EXAMINE IF A FALSE-POSITIVE OR CONTAMINATION IS SUSPECTED



Pre-Analytical

- Were specimens collected properly?
 - › *Contact the TB clinic/hospital/healthcare provider for information*

- Were specimen containers labeled correctly?
 - › *Contact the TB nurse/clinic/hospital*

- Were specimen requisition forms completed properly and matched to specimen container labels?
 - › *Contact the laboratory accessioning department/TB clinic or compare form and container, if available*

- Were known MTBC positive patients' specimens collected on the same day in the same location?
 - › *Inquire with the TB clinic/hospital*



Analytical

- Review patient's past and current test results
- Speak with TB program or healthcare provider to determine if patient symptoms are consistent with TB disease
- Review specimen processing logs (the day of the suspected contamination and days prior to and after potential contamination) for:
 - › Clerical errors
 - › Laboratory technologist who performed the testing
 - › Order of patient specimens processed
 - › Anything unusual with the specimen, equipment, safety, processing procedure, controls, etc. that may be annotated
 - › Other known positive patient test specimens in the same processing batch
 - › Whether controls or proficiency testing specimens were processed in the same batch or biological safety cabinet (BSC)
 - › Amount of time it took for inoculated cultures to grow

Continued on next page 



- Review standard operating procedures (SOP) and workflow to rule out processes that could have led to contamination
- Ensure staff are following the SOP and/or observe technique for best practices
- Ensure identification testing passed appropriate quality control (QC)
- Review lot numbers of manufacturer reagents, media, kits, etc. and ensure all media have passed necessary sterility checks and QC
- Assess if pipettes/supplies used for specimen processing are also used with positive cultures and susceptibility testing
- Review laboratory cleaning/decontamination records
- Review BSC records to ensure BSC is operating properly
- Review autoclave records to ensure autoclave is operating properly
- Review nonconforming events (NCEs) for occurrences that may have caused a false-positive



Post-Analytical

- Review test results to ensure the correct result was reported for the correct patient
- Verify that the correct report was sent to the healthcare provider
- Determine whether laboratory results and associated reporting language were received and interpreted correctly by the healthcare provider
- Ensure that data entry or transcription errors did not occur when results were entered into the laboratory information management system (LIMS)
- Review negative specimen processing control results to ensure growth was not observed
- Ensure quality assurance procedures have been followed (e.g., secondary review)
- Cross-reference genotypes for specimens in question to any known positives processed the same day or to the H37Rv/Ra control strain