

**ATTENTION, EMERGENCY RESPONDERS:
Interim Guidance on Recommended Use of NIOSH-approved
CBRN Full Facepiece, Non-Powered Air-Purifying Respirators
(APR) for Protection against CBRN Agents**

Deleted: Author: Mike Bergman,
EG&G Technical Services, Inc.¶

This is a draft provided for comment to the InterAgency Board for Equipment Standardization and InterOperability.

This information is distributed solely for the purpose of pre dissemination peer review under applicable information quality guidelines. It has not been formally disseminated by the National Institute for Occupational Safety and Health. It does not represent and should not be construed to represent any agency determination or policy.

DEPARTMENT OF HEALTH AND HUMAN SERVICES
Centers for Disease Control and Prevention
National Institute for Occupational Safety and Health

June 2005

This document is in the public domain and may be freely copied or reprinted.

Disclaimer

Mention of any company or product does not constitute endorsement by NIOSH.

Ordering Information

DRAFT

FOREWORD

With the escalation of worldwide acts of terrorism, the personal protection equipment needs of emergency responders and follow-on support personnel continue to evolve. As a result, The National Institute for Occupational Safety and Health (NIOSH) is rapidly meeting these challenges by developing new standards for the certification of respirators for protection against chemical, biological, radiological, and nuclear (CBRN) agents.

This interim guidance document provides guidelines for the selection and use of NIOSH-approved full facepiece nonpowered, air-purifying respirators (APR) for protection against CBRN agents. Additional guidance will supplement this document in the future. Users should periodically check the NIOSH National Personal Protective Technology Laboratory website, as follows, for user guidelines for all types of CBRN respirators.

<http://www.cdc.gov/niosh/npptl/default.html>

The NIOSH CBRN APR respirator certification standards are the result of collaboration between numerous federal agencies, the United States Army Research, Development and Engineering Command (RDECOM), private certification organizations, and respirator user groups. In 1999, the InterAgency Board for Equipment Standardization and InterOperability (IAB) identified the need for the development of standards or guidelines for respiratory protection equipment as a top priority. NIOSH, the National Institute for Standards and Technology (NIST), the National Fire Protection Association (NFPA), and the Occupational Safety and Health Administration (OSHA) entered into a Memorandum of Understanding defining each agency or organization's role in developing, establishing, and enforcing standards or guidelines for respirators. The collaborative standards development and research efforts resulted in the NIOSH CBRN APR certification standard dated March 7, 2003.

The NIOSH National Personal Protective Technology Laboratory is meeting the demands for CBRN protection through the development of CBRN respirator certification standards, administering the NIOSH respirator certification program, and developing user guidelines for CBRN respirators.

As emergency responders face the new challenges and demands of the evolving emergency response environment this interim CBRN APR guidance document will be a vital tool in their respirator training programs.

CONTENTS

Foreword	3
Acknowledgements	5
Abbreviations	6
1. Overview	8
2. CBRN Agent Canister Protections	8
Chemical Agents – 110 Agents	9
Biological Agents – 13 Agents	10
Radiological and Nuclear Agents – 16 Agents	11
3. NIOSH Recommended Use of CBRN APR	11
Use Criteria	11
NIOSH Cautions and Limitations	12
CBRN Use versus Industrial Use	14
Escape Contingency	15
Canister Interchangeability	15
Occupational Exposure Limits for CBRN Chemical Agents	15
4. Service Life and Change Schedules	16
5. Canister Capacity (Cap) Selection	17
6. CBRN Respirator Use Life (CRUL)	18
7. Methods for Establishing a Change Schedule	18
CBRN Respirator Use-Life (CRUL)	19
Software	19
Manufacturer's Test Data	20
Rules of Thumb	20
Inorganic Chemicals	20
8. Developing NIOSH CBRN Respirator Guidance	21
Glossary	22
References Cited	26
Related Reading	27
Appendix: CBRN APR Canister Agent Threat Protections	28

ACKNOWLEDGEMENTS

DRAFT

ABBREVIATIONS

AEGL	acute exposure guideline level (National Research Council and Environmental Protection Agency)
AC	hydrogen cyanide
APF	assigned protection factor
APR	air-purifying respirator
CG	phosgene (choking agent)
CK	cyanogen chloride
DP	diphosgene (choking agent)
CBRN	chemical, biological, radiological, and nuclear
CDC	Centers for Disease Control and Prevention
CWA	chemical warfare agent
CRUL	CBRN respirator use life
EPA	Environmental Protection Agency
GA	Tabun (nerve agent)
GB	Sarin (nerve agent)
GD	Soman (nerve agent)
GF	Cyclohexyl Sarin (nerve agent)
H	Sulfur mustard (blister agent)
HD	distilled Sulfur mustard (blister agent)
HN-1, HN-2, HN-3	nitrogen mustards (blister agents)
IDLH	immediately dangerous to life or health
IND	improvised nuclear device
L, L-1, L-2, L-3	Lewisite (blister agents)

MUC	maximum use concentration
NIOSH	National Institute for Occupational Safety and Health
NPPTL	National Personal Protective Technology Laboratory
OSHA	Occupational Safety and Health Administration
PEL	permissible exposure limit (OSHA)
PS	chloropicrin (choking agent)
RDD	radiological dispersive device
RDECOM	U.S. Army Research, Development, and Engineering Command (formerly SBCCOM)
REL	recommended exposure limit (NIOSH)
SBCCOM	U.S. Army Soldier and Biological Chemical Command (now RDECOM)
SCBA	self-contained breathing apparatus
TIC	toxic industrial chemical
TIM	toxic industrial material
TRA	test representative agent
UI	user's instructions (respirator operations manual)
VX	VX (nerve agent)

1. Overview

The National Institute for Occupational Safety and Health (NIOSH) is testing and certifying nonpowered air-purifying respirators (APR) for use by emergency responders in atmospheres which contain chemical, biological, radiological, and nuclear (CBRN) agents. This document serves as interim guidance on the selection and use of NIOSH-approved CBRN APR until further guidance is developed.

The CBRN APR is a full facepiece, tight-fitting respirator, which uses an air-purifying canister to remove particulates, gases, and vapors from the air. This document describes NIOSH recommendations for CBRN APR use and the specific CBRN agent protections provided by the respirator. It also provides guidance for selecting the canister capacity level (also called the "Cap" level) and guidance for establishing a canister change schedule. A unique time-use limitation known as CBRN respirator use life (CRUL), which designates how long a CBRN respirator can be used in a chemical warfare agent environment, is also discussed.

A CBRN APR can be identified by its NIOSH canister label affixed to the canister. The user should also be familiar with the NIOSH matrix-style APR canister approval label and matrix-style APR respirator approval label, both of which are required to be provided with the manufacturer's user instructions. The matrix-style labels identify the unique parts and accessories which compose the NIOSH-approved respirator configuration. A detailed discussion with examples of the canister label and matrix-style approval labels is located on the NIOSH NPPTL website at:

Attention, Emergency Responders: How To Determine If Your Air-Purifying Respirator (APR) Is NIOSH Certified For CBRN Environments
<http://www.cdc.gov/niosh/npptl/topics/respirators/cbrnapproved/apr/>

2. CBRN Agent Canister Protections

CBRN agents are chemical, biological, radiological, and nuclear inhalation hazards which have the potential to be released during acts of terrorism. During the development of the certification standard for the CBRN self-contained breathing apparatus (SCBA), NIOSH collaborated with the U.S. Army Research, Development and Engineering Command (RDECOM), formerly Soldier and Biological Chemical Command (SBCCOM), to perform a hazard assessment to identify a list of potential terrorist threats. As part of the CBRN APR certification testing, NIOSH uses 11 test representative agent (TRA) chemicals (10 gases and 1 particulate oil aerosol) for canister filtration testing to represent 139 of the CBRN agents recognized by NIOSH as potential terrorist threats. As other CBRN canister protections are identified, NIOSH will inform users accordingly.

NIOSH has categorized the 139 agents into seven test representative agent (TRA) families: organic vapors, acid gases, base gases, hydrides, nitrogen oxides, particulates

(including CBRN particulates), and formaldehyde. The agents belonging to each of the seven families are listed in the appendix at the end of this document. The 139 CBRN agents are further subclassified as 110 chemical agents, 13 biological agents, and 16 radiological and nuclear agents.

2a. Chemical Agents – 110 Agents

The CBRN APR provides protection against 110 chemical gases, vapors, and solid/liquid chemical aerosols identified as potential CBRN agents by the NIOSH/U.S. Army RDECOM threat analysis.

The chemical protections include specific chemical warfare agents [including, but not limited to, GB (Sarin) and HD (sulfur mustard)], and specific toxic industrial chemicals. Among the CBRN APR chemical protections are:

- Choking agents, which include, but are not limited to, phosgene (CG), diphosgene (DP), chlorine, and chloropicrin (PS). Choking agents attack lung tissue, primarily causing pulmonary edema. They cause irritation to the bronchi, trachea, larynx, pharynx, and nose. Initial symptoms may include tears, dry throat, coughing, choking, tightness of the chest, nausea, vomiting, and headache (DOD 1995).

Additional information on choking agents can be found at the following website:
CDC Emergency Preparedness and Response, Choking/Lung/Pulmonary Agents
<http://www.bt.cdc.gov/agent/pulmonary/>

- Blood agents (cyanogens), which include, but are not limited to, hydrogen cyanide (AC) and cyanogen chloride (CK). These agents are transported by the blood to all body tissues where the agent interferes with tissue oxygenation. The effects on the brain can lead to cessation of respiration, followed by cardiovascular collapse (DOD 1995).

Additional information on blood agents can be found at the following website:
CDC Emergency Preparedness and Response, Blood Agents
<http://www.bt.cdc.gov/agent/agentlistchem-category.asp#blood>

- Nerve agents, which include, but are not limited to, GB (Sarin), GA (Tabun), GD (Soman), GF (Cyclohexyl Sarin), and VX. Nerve agents inhibit cholinesterase (ChE) enzymes. This inhibition permits acetylcholine (Ach), which transmits many nerve impulses, to collect at various sites. Physiological effects of G-series and V-series agents can include miosis (contraction of the pupil), dim vision and headache, difficulty breathing, convulsions, coma, and death (DOD 1995).

Additional information on nerve agents can be found at the following website:

CDC Emergency Preparedness and Response, Nerve Agents
<http://www.bt.cdc.gov/agent/nerve/>

- Blister agents (vesicants) which include, but are not limited to, HD (sulfur mustard), nitrogen mustard (HN-1, HN-2, HN-3) and lewisite (L, L-1, L-2, L-3). Blister agents produce reddening and blistering of the skin or any other part of the body they contact. Eye exposure results in reddening of the eyes and temporary blindness or permanent effects. Inhaled mustard damages mucous membranes and the respiratory tract (DOD 1995).

Additional information on blister agents can be found at the following website:

CDC Emergency Preparedness and Response, Blister Agents/Vesicants
<http://www.bt.cdc.gov/agent/vesicants/index.asp>

2b. Biological Agents - 13 Agents

Biological agents consist of micro-organisms such as pathogens (which include disease causing bacteria, rickettsiae, and viruses) and toxins (DOD 2005). The effects of exposure and means of airborne dissemination vary depending on the agent type. Biological agents may be disseminated as aerosols, liquid droplets (toxins only), or dry powders (Zajtchuk 1997).

The CBRN APR canister provides protection from airborne biological agents by using P100 filter media to filter agents from the air. The CBRN APR canister provides protection against the following 13 biological terrorism agents identified in the NIOSH/U.S. Army RDECOM threat analysis: anthrax, brucellosis, glanders, pneumonic plague, tularemia, query (Q) fever, smallpox, Venezuelan equine encephalitis, viral hemorrhagic fevers, T-2 mycotoxins, botulism, ricin, and Staphylococcus enterotoxin B.

Additional information on bioterrorism agents can be found at the following website:

CDC reference page for Bioterrorism Agents/Diseases
<http://www.bt.cdc.gov/agent/agentlist.asp>

2c. Radiological and Nuclear Agents - 16 Agents

Radiological refers to particulate-borne radiation dispersed by detonation of a radiological dispersive device (RDD) or "dirty bomb". The principal type of RDD combines a conventional explosive, such as dynamite, with radioactive material (NRC date). When the dynamite or other explosives are detonated, the blast disperses radioactive material into the air. The radioactive materials used in a dirty bomb would probably not create enough radiation exposure to cause immediate serious illness, except to those people who are very close to the blast site. However, the radioactive dust and smoke from the blast site could be dangerous to health if inhaled (CDC Fact sheet: *Frequently Asked Questions (FAQs) About Dirty Bombs*.
<http://www.bt.cdc.gov/radiation/dirtybombs.asp>).

Nuclear refers to particulate-borne radiation dispersed by detonation of an improvised nuclear device (IND). An IND is intended to cause a yield-producing nuclear explosion. An IND could consist of diverted nuclear weapon components, a modified nuclear weapon, or an indigenous-designed device. INDs can be categorized into two types: implosion and gun assembled. Unlike RDDs, which can be made with almost any radioactive material, INDs require fissile material—highly enriched uranium or plutonium—to produce nuclear yield (Central Intelligence Agency. Fact sheet: *Terrorism CBRN: Materials and Effects*.

http://www.cia.gov/cia/reports/terrorist_cbrn/terrorist_CBRN.htm)

The CBRN APR canister provides respiratory protection from particulate-borne radiation (liquid and solid aerosols) by capturing particulates in the canister to prevent them from being inhaled into the body and causing internal damage. The CBRN APR canister uses P100 filter media to filter particulates from the air. The protection is defined for filtering particulates, but not for radiological gases or vapors. Airborne particulates have the ability to carry radioactive alpha and beta particles released from the atomic nuclei of an unstable isotope. Additionally, some airborne particulates may be emitters of high energy gamma radiation, which are photons emitted from the atomic nuclei of a substance undergoing radioactive decay.

The CBRN APR canister *does not* shield against or filter out radiation itself (alpha and beta particles or gamma radiation), but rather *filters* the larger particulates which carry alpha and beta particles or emit gamma radiation. Radiation may still continue to be emitted through the canister by the radioactive particulates once they are collected in the canister.

The CBRN APR canister provides protection against the following 16 radiological and nuclear particulate-borne agents identified in the NIOSH/U.S. Army RDECOM threat analysis: Hydrogen 3, Carbon 14, Phosphorous 32, Cobalt 60, Nickel 63, Strontium 90, Technetium 99m, Iodine 131, Cesium 137, Promethium 147, Thallium 204, Radium 226, Thorium 232, Uranium 235 & 238, Plutonium 239 and Americium 241.

Additional information on radiological and nuclear agents can be found at the following website:

CDC reference page for Radiation Emergencies
<http://www.bt.cdc.gov/radiation/index.asp>

3. NIOSH Recommended Use of CBRN APR

3a. Use Criteria

The proper use of a respirator is a complex process. It requires knowledge about respirator selection for a specific contaminant or environment, proper respirator fit, and awareness of the protections and limitations of the respirator.

The CBRN APR is only recommended for use when **ALL** of the following criteria are met:

- The types of inhalation hazards and their concentrations have been identified.
- The CBRN APR canister is capable of filtering the identified chemical.
- Oxygen concentration is known to be at least 19.5% by volume¹.
- The contaminant concentration is less than the immediately dangerous to life or health (IDLH) limit for the particular hazard².
- The contaminant concentration is less than the maximum use concentration (MUC) of the respirator.³
- A canister change schedule has been developed if the contaminant is a gas or vapor.
- The user has received a fit test prior to being assigned a respirator.
- The user has received complete training on the operational use, protections, and limitations of their unique CBRN APR system⁴.

3b. NIOSH Cautions and Limitations

The NIOSH Cautions and Limitations for the CBRN APR are listed on the CBRN full canister label located on the APR canister housing. They are also located on the matrix-style canister label and matrix-style APR respirator approval labels required to accompany the manufacturer's user instructions. Cautions and limitations lettered as A, I, J, L, M, O, and S (located in Section 2 of these labels) apply when non-CBRN conditions are present. Cautions and limitations lettered R, T, V, W, X, Y, Z, HH, QQ, and UU (located in Section 3 of these labels) additionally apply to the non-CBRN cautions and limitations. These limitations are not all inclusive. The respirator manufacturer also may identify further cautions and limitations for their respirators. In addition, regulatory agencies also may place a limit on the use of CBRN APR in their response plans. **In all cases, all NIOSH Cautions and Limitations must be strictly followed.**

¹ NIOSH defines an oxygen-deficient atmosphere as an atmosphere which contains an oxygen partial pressure of less than 148 millimeters of mercury (19.5 percent by volume at sea level) (NIOSH 2005).

² The 2004 NIOSH Pocket Guide to Chemical Hazards lists IDLH values established by NIOSH for many chemicals.

³ The maximum use concentration (MUC) for the CBRN APR is the maximum atmospheric concentration of a hazardous substance for which an employee can be expected to be protected and is determined by the lesser of the following:

1. The assigned protection factor (APF) multiplied by the contaminant's occupational exposure limit. The APF of a CBRN APR equipped with a CBRN canister is equal to 50 when: a) the respirator user adheres to complete program requirements (such as the ones required by OSHA in 29 CFR 1910.134), b) the respirator is assembled in its NIOSH-approved configuration, and c) the user has been individually fit tested to assure a proper fit.
2. The respirator manufacturer's MUC for a hazardous substance (if any).
3. The immediately dangerous to life or health (IDLH) value of the contaminant.

⁴ CBRN APR training should include complete respirator training program requirements (such as the ones required by OSHA in 29 CFR 1910.134) including the individual manufacturer's user instructions, and additional training specific to CBRN environments and exposures.

The following NIOSH Cautions and Limitations appear in Section 2 of the full canister label and canister and respirator matrix-style approval labels:

- A Not for use in atmospheres containing less than 19.5 percent oxygen.
- I Contains electrical parts which have not been evaluated as an ignition source in flammable or explosive atmospheres by MSHA / NIOSH.
Note: Caution and limitation "I" will not be present on units which have met these evaluation requirements by MSHA/NIOSH.
- J Failure to properly use and maintain this product could result in injury or death.
- L Follow the manufacturer's User's Instructions for changing canisters.
- M All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations.
- O Refer to User's Instructions and/or maintenance manuals for information on use and maintenance of these respirators.
- S Special or critical user's instructions and/or specific limitations apply. Refer to User's Instructions before donning.

(Note: Caution and limitation "S" will only be on the NIOSH approval label if specified by the manufacturer in the user's instructions. When "S" appears on the NIOSH approval label, the corresponding cautions and limitations, that apply under "S", will be explained in a designated section of the manufacturer's user's instructions (UI).

The following NIOSH Cautions and Limitations appear in Section 3 of the full canister label and canister and respirator matrix-style approval labels and are specific to use in CBRN environments:

- R Some CBRN agents may not present immediate effects from exposure, but can result in delayed impairment, illness, or death.
- T Direct contact with CBRN agents requires proper handling of the respirator after each use and between multiple entries during the same use. Decontamination and disposal procedures must be followed. If contaminated with liquid chemical warfare agents, dispose of the respirator after decontamination.
- V Not for use in atmospheres immediately dangerous to life and health or where hazards have not been fully characterized.
- W Use replacement parts in the configuration as specified by the applicable regulations and guidance.

- X Consult manufacturer's User's Instructions for information on the use, storage, and maintenance of these respirators at various temperatures.
- Y This respirator provides respiratory protection against inhalation of radiological and nuclear dust particles. Procedures for monitoring radiation exposure and full radiation protection must be followed.
- Z If during use an unexpected hazard is encountered such as a secondary CBRN device, pockets of entrapped hazard or any unforeseen hazard, immediately leave the area for clean air.
- HH When used at defined occupational exposure limits, the rated service time cannot be exceeded. Follow established canister change schedules or observe End-of-Service-Life Indicators to ensure that canisters are replaced before breakthrough occurs.
- QQ Use in conjunction with personal protective ensembles that provide appropriate levels of protection against dermal hazard. Failure to do so may result in personal injury even when the respirator is properly fitted, used, and maintained.
- UU The respirator should not be used beyond eight (8) hours after initial exposure to chemical warfare agents to avoid possibility of agent permeation. If liquid exposure is encountered, the respirator should not be used for more than two (2) hours.

3c. CBRN Use Versus Industrial Use

At a CBRN event, respirators providing the highest level of protection, NIOSH-approved CBRN SCBA, should initially and routinely be used until hazard types and concentrations are quantified and determined to be at lower acceptable exposure limits before utilizing CBRN APR (See Section 3a Use Criteria above).

The NIOSH approval requirements for the CBRN APR were developed primarily to address the unique needs of emergency responders in CBRN environments where conditions are appropriate for air-purifying respirator use. However, manufacturers may gain dual approvals for respirator systems which utilize the same facepiece part numbers for both CBRN use and industrial (non-CBRN) use. An example of such a case could be that the same facepiece part number can be used with a CBRN canister or a canister with only industrial protections (for example, P100 protection). A properly maintained system can be used for both industrial and CBRN use if it is always assembled in the NIOSH-approved configuration specified by the assembly matrixes found in the user's instructions. All working parts and accessories must be free from damage.

NIOSH only approves canisters with either CBRN protection or industrial protections. CBRN canisters are to be used for CBRN response events only, not routine industrial use, and are to remain sealed in their original packaging until needed for such use. For both

CBRN canisters and industrial canisters, users should observe canister shelf life restrictions designated by the manufacturer and dispose of any canisters which have expired shelf life dates.

3d. Escape Contingency

The NIOSH certification test criteria for the CBRN canister includes gas testing using ten test representative agent (TRA) chemical gases at high concentrations and high flow rates. This test criterion was implemented to address high contaminant concentrations and elevated user breathing rates that emergency responders may experience while escaping from a secondary event in a CBRN environment. If unknown or high levels of a hazard are generated after entry into the work area due to a secondary hazard source (e.g., a secondary explosive device) or work activities (e.g., a spill), the user should immediately leave the area while continuing to wear the CBRN APR. After using the respirator to escape, the canister should be replaced before reusing the respirator. There is a potential for the canister's filtration capacity to become depleted during an escape due to a higher than normal breathing rate and possibly higher contaminant concentration of the escape incident.

3e. Canister Interchangeability

The design of the CBRN APR enables the interchangeable operation of canisters by standardizing the design requirements for the male thread of the canister and the interface connector on the respirator facepiece. Respirators are only in their NIOSH-approved configuration if they are assembled correctly and only with the parts listed in their individual approval matrixes. A CBRN APR respirator assembled with an interchanged canister from a different manufacturer is not in a NIOSH-approved configuration; however, the interchangeability of canisters between different manufacturers' facepieces helps to alleviate the complications of matching facepieces to canisters at an emergency response scene involving large numbers of response workers. The interchangeability provision for canisters would typically be used under emergency conditions when there is a restrictive supply of CBRN canisters for a manufacturer's facepiece. Other CBRN canisters from a different manufacturer may be used under these emergency conditions. The interchangeability provision only applies to CBRN canisters and not to other subcomponents or assemblies. **Users should not interchange canisters until they are instructed to do so by the incident commander or other designated commanding authority.**

3f. Occupational Exposure Limits for CBRN Chemical Agents

Various occupational exposure limits for CBRN chemical agents have been established by different government agencies and private organizations. Comprehensive exposure values may not be available from one particular organization. Examples of such occupational exposure limits are NIOSH recommended exposure limits (RELs), OSHA permissible exposure limits (PELs), immediately dangerous to life and health (IDLH) values established by different organizations, and acute exposure guideline levels

(AEGLs) established by the National Research Council and the Environmental Protection Agency.

OSHA and NIOSH recently collaborated on personal protective equipment (PPE) selection guidelines, including respirators, for CBRN environments. As part of this joint guideline, recommendations are made for selecting the appropriate type of respirator (self-contained breathing apparatus, powered air-purifying respirator, or air-purifying respirator) based on airborne concentrations of select nerve agents and blister agents and maximum user times in agent environments. The exposure limits presented in the joint OSHA/NIOSH guidance for nerve agents and blister agents are based on acute exposure guideline levels (AEGLs) published by the National Research Council and the Environmental Protection Agency and IDLH values either proposed or accepted by the U.S. Army. The AEGL limits characterize the risk to the general population during a one-time accident and emergency scenario with time limits not to exceed eight hours. The respirator and protective clothing recommendations made in the OSHA/NIOSH CBRN guideline are based on protecting the responder at the lowest recommended AEGL-1 level for a given exposure duration. At the AEGL-1 level an individual could experience notable discomfort, irritation, or certain asymptomatic, non-sensory effects. However, the effects are not disabling and are transient and reversible upon cessation of exposure (Federal Register 2001).

Users are referred to the joint OSHA/NIOSH guidance to aid in respirator and protective clothing selection for a CBRN event. The joint OSHA/NIOSH guidance is located on the OSHA website at:

OSHA/NIOSH Interim Guidance – August 30, 2004
Chemical - Biological - Radiological - Nuclear (CBRN)
Personal Protective Equipment Selection Matrix for Emergency Responders
<http://www.osha.gov/SLTC/emergencypreparedness/cbrnmatrix/index.html>

4. Service Life and Change Schedules

CBRN canisters have a limited time that they can remain in service before they must be changed, that is, disposed of and replaced with a new canister. A *service life* must be determined for canisters used to protect against gases or vapors. Once the service life has been determined, an appropriate *change schedule* can be implemented.

Service life is sometimes referred to as “breakthrough time.” It is the length of time required for an air-purifying element to reach a specific effluent concentration. Service life is determined by the type of substance being removed, the concentration of the substance, the ambient temperature, the specific element being tested (cartridge or canister), the flow rate resistance, and the selected breakthrough value (NIOSH 2005). The various methods of estimating service life for developing a change schedule are discussed below in Section 7 Methods for Establishing a Change Schedule.

A *change schedule*⁵ is the time interval after which a used canister is replaced with a new one. An appropriate change schedule assures that the canister will be changed before the downstream (inside the respirator) concentration exceeds a predetermined breakthrough concentration. OSHA suggests that employers apply a safety factor to the service life determination to assure that the change schedule is a conservative estimate. The change schedule also must be convenient to implement and enforce.

5. Canister Capacity (Cap) Selection

An individual knowledgeable in determining canister service life (see Section 4 above) should determine the appropriate canister capacity (Cap) level to select and when canisters need to be changed. CBRN APR canisters are approved at six different levels of canister capacity (e.g., each level is referred to as a specific Cap). Currently, at the date of this web posting (*INSERT DATE HERE*), all CBRN APR respirators have only been approved with Cap 1 canisters.

The different canister capacity levels relate to relative lengths of service time for the filtration of gases and vapors when exposure conditions are similar (similar contaminant concentration, ambient temperature, relative humidity, and user breathing rate). The higher the Cap designation, the longer the service time; conversely, the lower the Cap designation, the shorter the service time when used under similar conditions (see Table 1 below). For example, a Cap 2 canister will have about twice the service life of a Cap 1 canister when gases and vapors are concerned. If the service life is determined to be too short to be practical, a canister with a higher Cap rating should be selected and a new service life determined, or a CBRN SCBA should be used instead. **It is important to note that selecting a higher Cap designation does not change the respirator maximum use concentration (MUC).**

Service life does not vary among cap levels for the filtration of solid and liquid particulates. All CBRN canisters, regardless of Cap level designation, should be changed immediately if breathing becomes difficult due to clogging by particulates. Canisters should also be changed immediately if they become damaged. All CBRN canisters should be changed in a clean area free from contamination.

⁵ The OSHA Respiratory Protection Standard [29 CFR1910.134] requires that employers implement change schedules for canisters where end-of-service-life indicators (ESLI) do not exist or are not appropriate for the work environment. Presently, no CBRN APR canisters are approved with ESLIs.

Table 1. CBRN APR Canister Capacity Relative Gas/Vapor Service Life

CBRN Canister Capacity Designation	Relative Gas/Vapor Service Life compared to a Cap 1 Canister*
Cap 1	-----
Cap 2	2X
Cap 3	3X
Cap 4	4X
Cap 5	6X
Cap 6	8X

*Assumes similar use conditions (similar contaminant concentration, ambient temperature, relative humidity, and user breathing rate).

6. CBRN Respirator Use Life (CRUL)

CBRN respirator use life (CRUL) for the CBRN APR is a time-use limitation for the entire respirator system (including the facepiece, canister, and all accessories) when the system is used in a chemical warfare agent environment.

NIOSH caution and limitation UU states specific time-use restrictions for CBRN respirator use life for CBRN APR:

UU The respirator should not be used beyond eight (8) hours after initial exposure to chemical warfare agents to avoid the possibility of agent permeation. If liquid exposure is encountered, the respirator should not be used for more than two (2) hours.

Chemical warfare agents in the context of CRUL are presently nerve agents and blister agents. NIOSH will notify users of additional CRUL restrictions as they become available.

Confirmed contamination presence using the appropriate quantitative or qualitative methods is the key to determining the 8-hour start point of use. If exposed to a chemical warfare agent *vapor* CRUL is eight hours and 2-hours if exposed to a chemical warfare agent *liquid*. The 8-hour (vapor) and 2-hour (liquid) use life means eight continuous hours or two continuous hours in a single shift, day, or event. The time intervals are continuous and cannot be divided; for example, the 8-hour period cannot be broken into four different 2-hour periods over the course of a day.

7. Methods for Determining a Change Schedule

There are several methods available for determining a canister change schedule for gases and vapors. It is important to remember that these methods apply only to gases and vapors, not to particulates. These methods include:

- CBRN Respirator use life (CRUL)
- Software
 - using the specific chemical
 - using a NIOSH Test Representative Agent (TRA) chemical
- Manufacturers' test data (for a specific chemical or NIOSH test representative agent (TRA) chemical)
- "Rules of Thumb"

7a. CBRN Respirator Use Life (CRUL)

The CBRN Respirator use life (CRUL) limitations described in Section 6 of this document apply when a CBRN APR is used in a chemical warfare agent environment. The CRUL restrictions apply regardless of the canister Cap level. It is still important to determine canister service life in case it would indicate a change schedule more frequent than the CRUL time limitations. In all cases, the entire CBRN APR system (facepiece, canister, and accessories) should not be reused and should be decontaminated and disposed of in a manner that is consistent with the type of contamination and any government regulations regarding contaminated items.

Chemical warfare agents in the context of CRUL are presently nerve agents and blister agents. NIOSH will continue to notify users of additional CRUL restrictions as they become available.

7b. Software

Data software available on the OSHA website and also on some manufacturers' websites can determine change schedules by allowing the user to choose a chemical from a list of chemicals in the database and enter other information, including employee work rate, environmental conditions (environmental temperature and relative humidity), contaminant concentration, and canister specifications.

The respirator manufacturer may have a software program on their website which includes their specific CBRN canister. Users should contact the manufacturer for questions about using a manufacturer's software program.

The software programs on the OSHA website can only calculate service life for single organic vapors and requires information on the canister (such as the bulk density of the packed bed), which can be obtained from the manufacturer. The OSHA software also allows for chemical data (such as molecular weight and vapor pressure) to be entered if the specific chemical is not listed in the software's database.

One such programs, The Advisor Genius, was developed in 1998 and is available at: http://www.osha.gov/SLTC/etools/respiratory/advisor_genius_wood/advisor_genius.html

Breakthrough (2003), another downloadable program similar to The Advisor Genius, has been updated to include corrections for relative humidity. Breakthrough is also limited to

only calculating change schedules for single organic vapors. The program is available on the OSHA website at:

http://www.osha.gov/SLTC/etools/respiratory/advisor_genius_wood/breakthrough.html.

Only the NIOSH test representative agent (TRA) chemicals can be used as surrogates for chemicals that are not included in the software programs. The NIOSH CBRN canister protections for chemical gases and vapors are established from the NIOSH certification tests using 10 chemical gases. For example, cyclohexane can be used as the surrogate for acrylonitrile because cyclohexane is used as the TRA for CBRN organic vapor threats, which includes acrylonitrile. The NIOSH TRA chemicals for each family of agents are identified in the appendix of this document. Only the TRA which corresponds to each chemical family can be used. A surrogate cannot be used for the acid gas family. This is because there are 5 TRA that represent the acid gas family, and it is not possible to input 5 chemicals at once into the software.

7c. Manufacturers' Test Data

The manufacturer may be able to provide guidance on canister service life based on their testing of a specific CBRN chemical threat or a NIOSH Test Representative Agent (TRA) chemical gas for that threat. However, for those agents which are chemical warfare agents, the CBRN respirator use life (CRUL) limitations stated in caution and limitation 'UU' apply and service life is limited to a maximum of 8 continuous hours for a CWA gas or vapor.

7d. Rules of Thumb

The *Rules of Thumb* may provide a rough estimation of canister service life, but only for single organic vapors (AIHA 1997). However, they should **NOT** be used as the sole method of determining service life. These rules state that organic vapor cartridges will last eight hours if the organic vapor has a boiling point of greater than 70°C, the vapor's concentration is less than 200 ppm, and the worker has a breathing rate of 30 liters per minute (moderate work). These rules also indicate that service life is inversely proportional to flow rate. In other words, if the flow rate is reduced by a factor of 2, the service life is increased by a factor of 2. A summary of the "Rules of Thumb" are available on the OSHA website at:

http://www.osha.gov/SLTC/etools/respiratory/rule_of_thumb/rule_of_thumb.html.

7e. Inorganic Chemicals

The software models on the OSHA website only calculate service life for organic vapors. However, some manufacturers' service life software may include certain inorganic chemicals. For service life information on inorganic chemicals (e.g., chlorine or ammonia) users should check the manufacturers' software package or contact the manufacturer for a service-life recommendation. The manufacturer may also have test data on the contaminant in question or its NIOSH test representative agent chemical.

8. Developing NIOSH CBRN Respirator Guidance

NIOSH will continue to develop user guidance for the CBRN APR and other classes of CBRN respirators. Future CBRN APR guidance will further explain selection criteria, how to assure proper fit, cautions and limitations of use, and maintenance and cleaning guidelines. Users should continue to check the NIOSH NPPTL website for updated user guidance documents. Users are always encouraged to contact the manufacturer of their respirator for specific recommendations on CBRN APR use and establishing canister change schedules.

Glossary

Agent – A force or substance that causes change or effects on an exposed substrate.

Acute Exposure Guideline Level (AEGL) – Published by the National Research Council and Environmental Protection Agency, AEGL limits characterize the risk to the general population during a one-time accident and emergency scenario with time limits not to exceed eight hours. The AEGLs represent threshold exposure limits for the general public and are applicable to emergency exposure periods ranging from ten minutes to eight hours.

Assigned Protection Factor (APF) – The NIOSH minimum anticipated workplace level of protection provided by a properly functioning respirator or class of respirators to a percentage of properly fitted and trained personnel. An APF of 50 is assigned to the CBRN APR.

Biological Agents – Biological agents consist of micro-organisms such as pathogens (which include disease-causing bacteria and viruses) and toxins.

Blister Agents (Vesicants) – Vesicants are highly reactive chemicals that combine with proteins, DNA, and other cellular components to result in cellular changes immediately after exposure. The most commonly encountered clinical effects include dermal (skin erythema and blistering), respiratory (pharyngitis, cough, dyspnea), ocular (conjunctivitis and burns), and gastrointestinal (nausea and vomiting). Blister agents include H (sulfur mustard), HD (distilled sulfur mustard), nitrogen mustards (HN-1, HN-2, HN-3) and Lewisite (L, L-1, L-2, L-3).

CBRN Protection – A NIOSH certification term that signifies that specific respirator systems have been evaluated, reviewed, and approved/certified by NIOSH as providing a defined level of protection against chemical, biological, radiological, and nuclear agents. The approval authority is authorized by the 42 Code of Federal Regulations, Part 84 (42CFR84) and the CBRN statement of standard for a designated respirator class.

CBRN Respirator Use Life (CRUL) – Use limitation applying to the CBRN APR of a continuous 8-hour period (when exposed to a chemical warfare agent vapor) or 2-hour period (when exposed to a chemical warfare agent liquid) beginning at the time of a confirmed exposure, after which the entire CBRN APR must be decontaminated and disposed. The 8-hour (vapor) and 2-hour (liquid) use life practice means eight continuous hours (vapor) or two continuous hours (liquid) in a single shift, day, or event. The time intervals are continuous and cannot be divided; for example, the 8-hour period cannot be broken into four different 2-hour periods over the course of a day. See Section 6 of this document for details.

Chemical Warfare Agents (CWA) – Chemical warfare agents in the context of the CBRN APR cautions and limitations UU pertaining to CBRN respirator use life (CRUL) include nerve agents, blister agents, chemicals that exhibit degrading or destructive effects on respirator materials, and other chemicals for which decontamination

procedures are unable to decontaminate the respirator to a safe level for reuse. NIOSH will continue to notify users of additional CRUL restrictions as they become available.

Dirty Bomb – A conventional explosive device that has been surrounded by or contaminated with some form of radioactive material.

Exposure – Contact of the respirator wearer or APR components to a chemical, biological, radiological or nuclear agent.

Fit Test – The use of a specific measurement protocol to qualitatively or quantitatively evaluate the fit of a respirator on an individual. See also qualitative fit test (QLFT) or quantitative fit test (QNFT)

Immediately Dangerous to Life or Health (IDLH) – Conditions that pose an immediate threat to life or health or conditions that pose an immediate threat of severe exposure to contaminants, such as radioactive materials, which are likely to have adverse cumulative or delayed effects on health.

Maximum Use Concentration (MUC) – The maximum atmospheric concentration of a hazardous substance from which an employee can be expected to be protected when wearing a respirator. MUC is determined by the assigned protection factor (APF) of the respirator or class of respirators and the exposure limit of the hazardous substance. The maximum use concentration (MUC) of exposure for a respirator is generally determined by multiplying a contaminant's occupational exposure limit (OEL) by the APF assigned to a specific class or type of respirator ($MUC = OEL \times APF$). An OEL can be a NIOSH recommended exposure limit (REL), an OSHA permissible exposure limit (PEL), a short term exposure limit, ceiling limit, peak limit, or any other exposure limit for a hazardous substance.

National Institute for Occupational Safety and Health (NIOSH)/The Institute – NIOSH is the federal agency responsible for conducting research and making recommendations for the prevention of work-related injury and illness. NIOSH is part of the Department of Health and Human Services' Centers for Disease Control and Prevention.

National Personal Protective Technology Laboratory (NPPTL) – Part of NIOSH, the mission of NPPTL is "To prevent work-related injury and illness by ensuring the development, certification, deployment, and use of personal protective equipment and fully-integrated intelligent ensembles."

NIOSH-Approved – Respirator systems that have been reviewed, tested, and approved/certified by NIOSH under 42 Code of Federal Regulation, Part 84 (42CFR84). NIOSH approves respirators as complete systems assembled only with the components designated in the individual NIOSH approval matrix issued with the respirator approval.

Nerve Agents – Nerve agents consist of a group of very toxic organophosphate chemicals. Nerve agents include GB (Sarin), GA (Tabun), GD (Soman), GF (Cyclohexyl

Sarin), and VX. Nerve agents cause effects on the human body by disrupting how nerves communicate and control muscles, glands, and organs.

Oxygen Deficient Atmosphere – An atmosphere which contains an oxygen partial pressure of less than 148 millimeters of mercury (19.5 percent by volume at sea level).

Nuclear Agents – Particulate-borne radiation dispersed by detonation of an improvised nuclear device (IND) or high/low yield nuclear detonation. An IND could consist of diverted nuclear weapon components or a modified nuclear weapon. INDs require fissionable material—highly enriched uranium or plutonium—to produce nuclear yield.

Penetration – The act or process of penetrating, piercing, or entering. As it relates to CBRN approval testing, it is a term that means exposing the respirator to specific quantities of chemical warfare agent with the express intent to determine if the agent is stopped or if it penetrates through air pressure boundaries or material interfaces into the breathing zone of the respirator.

Permeation – The action of passing through the openings or interstices of a substrate at the surface level or the molecular level. As it relates to CBRN approval testing, it is a term that means exposing the respirator to a specific quantity of chemical warfare agent with the express intent to see if the agent is stopped and runs off, or if it beads up and starts to permeate through air pressure boundaries or material surfaces into the breathing zone of the respirator or respirator accessories.

Permissible Exposure Limit (PEL) – An enforceable regulatory limit set by the Occupational Safety and Health Administration (OSHA) on the amount or concentration of a substance in the air. PELs are set to protect workers against the health effects of exposure to hazardous substances and are based on an 8-hour, time weighted average exposure.

Radiological Agents – Particulate-borne radiation dispersed by detonation of a radiological dispersive device (RDD) or “dirty bomb.”

Recommended Exposure Limit (REL) – An occupational exposure limit recommended by NIOSH as being protective of worker health and safety over a working lifetime. RELs are time-weighted average concentrations for up to a 10-hour workday during a 40-hour workweek.

Respiratory Protection Program – A written program document which establishes procedures and practices to ensure that respirators are properly selected, used, fit-tested, and maintained. The program is administered by a suitable trained program administrator. The OSHA Respiratory Protection standard (29 CFR 1910.134) defines the requirements of a complete respiratory protection program.

Test Representative Agent (TRA) – Refers to any of the 11 chemicals (10 gases and 1 particulate aerosol) NIOSH uses for certification testing of the CBRN APR canister.

NIOSH uses the 11 TRA chemicals to represent 139 agents which NIOSH recognizes as CBRN threats. See the Appendix of this document for a list of the 139 CBRN threat protections of the CBRN canister and their corresponding TRA chemicals.

Toxic Industrial Chemicals (TICs)/Toxic Industrial Materials (TIMs) – a variety of chemicals used in various industrial processes which can kill, seriously injure, or incapacitate people.

User's Instructions (UI) – A NIOSH recognized manufacturer publication required to be submitted to NIOSH as part of a certification application requesting NIOSH approval. The UI are included with every new purchase of a NIOSH-approved respirator.

References Cited

Zajtchuk, R, Brigadier General, et al. [1997]. Textbook of Military Medicine: Medical Aspects of Chemical and Biological Warfare, Chapter 20, Use of Biological Weapons. Office of the Surgeon General.

CIA (Central Intelligence Agency). Fact sheet: Terrorism CBRN: Materials and Effects. http://www.cia.gov/cia/reports/terrorist_cbrn/terrorist_CBRN.htm.

CDC (Centers for Disease Control and Prevention). Fact sheet: Frequently Asked Questions (FAQs) About Dirty Bombs. <http://www.bt.cdc.gov/radiation/dirtybombs.asp>.

DOD (Department of Defense) [2005]. Potential Military Chemical/Biological Agents and Compounds. FM 3-11.9/MCRP 3-37.1B/NTRP 3-11.32/AFTTP (I) 3-2.55.

DOD (Department of Defense) [1995]. Treatment of Chemical Agent Casualties and Conventional Military Chemical Injuries. FM 8-285/Navy Medical (NAVMED) P-5041/Air Force Joint Manual (AFJMAN) 44-149/Fleet Marine Force Manual (FMFM) 11-11.

NIOSH [2005] Respirator Selection Logic 2004. Cincinnati, OH: U.S Department of Health and Human Services, Public Health Service, Centers for Disease Control, National Institute for Occupational Safety and Health, DHHS (NIOSH) Publication No. 2005-100.

NRC (Nuclear Regulatory Commission). Fact Sheet: Dirty Bombs. <http://www.nrc.gov/reading-rm/doc-collections/fact-sheets/dirty-bombs.html>.

The Occupational Environment, Its Evaluation and Control [1997], (Respiratory Protection, Sec. 5: Chapter 36). Fairfax, VA: American Industrial Hygiene Association Press.

66 Fed. Reg. 85 [2001]. Proposed AEGL Values. National Advisory Committee for Acute Exposure Guideline Levels (AEGLs) for Hazardous Substances.

Related Reading

Dower JM, Metzler RW, Palya FM, Peterson JA, Pickett-Harner M [2001]. Chemical and Biological Respiratory Protection Workshop Report. U.S Department of Health and Human Services, Public Health Service, Centers for Disease Control, National Institute for Occupational Safety and Health. Cincinnati, OH: DHHS (NIOSH) Publication No. 2000-122.

IAB (InterAgency Board). Annual Report 2003 and Standardized Equipment List 2004. Arlington: VA: <http://www.iab.gov/download/AnnualReport2003.pdf>.

Jackson B, Baker J, Ridgely M, Bartis J, Linn H [2004]. Protecting Emergency Responders, Volume 3, Safety Management in Disaster and Terrorism Response. U.S Department of Health and Human Services, Public Health Service, Centers for Disease Control, National Institute for Occupational Safety and Health, Cincinnati, OH: DHHS (NIOSH) Publication No. 2004-144.

NIOSH [1987]. Respirator Decision Logic. U.S Department of Health and Human Services, Public Health Service, Centers for Disease Control, National Institute for Occupational Safety and Health. Cincinnati, OH: DHHS (NIOSH) Publication No. 87-108.

NIOSH [2002]. Protecting Workers from Anthrax Infection. U.S Department of Health and Human Services, Public Health Service, Centers for Disease Control, National Institute for Occupational Safety and Health. Cincinnati, OH: DHHS (NIOSH) Publication No. 2002-142.

NIOSH [2002]. Protecting Workers at the World Trade Center Site. U.S Department of Health and Human Services, Public Health Service, Centers for Disease Control, National Institute for Occupational Safety and Health. Cincinnati, OH: DHHS (NIOSH) Publication No. 2002-143.

Protecting Emergency Responders, Lessons Learned from Terrorist Attacks [2002]. RAND Conference Proceedings, CF-176-OSTP.

Revoir W., Ching-Tsen Bien [1997]. *Respiratory Protection Handbook*. CRC Lewis Publishers: Boca Raton, NY.

Appendix: CBRN APR Canister Agent Threat Protections

A CBRN APR canister provides protection against a minimum of 139 identified CBRN agents, which are classified into the following 7 families: **Organic Vapors (61), Acid Gases (32), Base Gases (4), Hydrides (4), Nitrogen Oxides (5), Particulates (32)** [composed of 3 chemical, 13 biological, and 16 radiological and nuclear particulate threats], and **Formaldehyde (1)**.

The test representative agents (TRA) NIOSH uses for certification testing to represent each agent family are listed.

<u>Acid Gas Family (32 agents)</u> 5 TRA chemicals (Cyanogen Chloride, Hydrogen Cyanide, Hydrogen Sulfide, Phosgene, and Sulfur Dioxide) are used for NIOSH certification testing to represent the Acid Gas Family.	
Boron tribromide	Hydrogen chloride
Boron trichloride	Hydrogen cyanide
Boron trifluoride	Hydrogen fluoride
Bromine	Hydrogen iodide
Bromine chloride	Hydrogen sulfide
Bromine trifluoride	Phosgene
Carbonyl fluoride	Phosphorus trichloride
Chlorine	Silicon tetrafluoride
Chlorine pentafluoride	Sulfur dioxide
Chlorine trifluoride	Sulfur trioxide
Chlorosulfonic acid	Sulfuric acid
Cyanogen chloride	Sulfuryl chloride
Dichlorosilane	Titanium tetrachloride
Ethyl phosphonous dichloride	Tungsten hexafluoride
Fluorine	Bromine pentafluoride
Hydrogen bromide	Hydrogen selenide
<u>Nitrogen Oxide Family (5 agents)</u> 1 TRA chemical (nitrogen dioxide) is used for NIOSH certification testing to represent the Nitrogen Oxide Family.	
Nitric acid	Nitrogen tetraoxide
Nitric acid, fuming	Nitrogen trioxide
Nitrogen dioxide	
<u>Base Gas Family (4 agents)</u> 1 TRA chemical (Ammonia) is used for NIOSH certification testing to represent the Base Gas Family.	
Allyl amine	Dimethyl hydrazine, 1,2
Ammonia	Methyl hydrazine
<u>Hydride Family (4 agents)</u> 1 TRA chemical (Phosphine) is used for certification testing to represent the Hydride Family.	
Arsine	Phosphine
Germane	Stibine

Formaldehyde Family (1 agent)	
<i>1 TRA chemical (Formaldehyde) is used for certification testing to represent the Formaldehyde Family.</i>	
Formaldehyde	
Organic Vapor Family (61 agents)	
<i>1 TRA chemical (Cyclohexane) is used for certification testing to represent the Organic Vapor Family.</i>	
Acetone cyanohydrin	Methanesulfonyl chloride
Acrylonitrile	Methyl orthosilicate
Allyl alcohol	Methyl parathion
Allyl chlorocarbonate	Methyl phosphonic dichloride
Bromoacetone	Mustard, lewisite mixture
Bromobenzylcyanide	Nitrogen mustard HN-1
Chloroacetone	Nitrogen mustard HN-2
Chloroacetonitrile	Nitrogen mustard HN-3
Chloroacetophenone	N-propyl chloroformate
Chloroacetyl chloride	O-chlorobenzylidene malononitrile
Chloropicrin	O-ethyl-s-(2isopropylaminoethyl)methyl phosphonothiolate
Chloropivaloyl chloride	Parathion
Crotonaldehyde	Perchloromethyl mercaptan
Cyclohexyl methyphosphonate	Phenyl mercaptan
Dibenz-(b,f)-1,4-oxazepine	Phenylcarbylamine chloride
Diketene	Phenyldichloroarsine
Dimethyl sulfate	Phosgene oximedichloroformoxime
Diphenylchloroarsine	Phosphorus oxychloride
Diphenylcyanoarsine	Sarin
Diphosgene	Sec-butyl chloroformate
Distilled mustard	Soman
Ethyl chloroformate	Tabun
Ethyl chlorothioformate	Tert-octyl mercaptan
Ethyl phosphonothioicdichloride	Tetraethyl dithiopyrophosphate
Ethyl phosphorodichloridate	Tetraethyl lead
Ethylene dibromide	Tetramethyl lead
Hexachlorocyclopentadiene	Tetranitromethane
Hexaethyl tetraphosphate	Trimethoxysilane
Iso-butyl chloroformate	Trimethylacetyl chloride
Iso-propyl chloroformate	VX
Lewisite	

Particulate Family Canister Protections (32 agents)	
<i>ITRA chemical (dioctyl phthalate) is used for certification testing to represent the Particulate Family.</i>	
<u>Particulate- Chemicals (3)</u>	
Sodium azide	
Adamsite	
Sodium fluoroacetate	
<u>Particulate- Biological (13)</u>	
Anthrax	Venezuelan equine encephalitis
Glanders	Brucellosis
Tularemia	Pneumonic Plague
Smallpox	Query (Q) Fever
T-2 Mycotoxins	Viral hemorrhagic fevers
Ricin	Botulism
Staphylococcus enterotoxin B	
<u>Particulate- Radiological / Nuclear (16)</u>	
Carbon 14	Hydrogen 3
Cobalt 60	Phosphorous 32
Strontium 90	Nickel 63
Iodine 131	Technetium 99m
Promethium 147	Cesium 137
Radium 226	Thallium 204
Uranium 235 & 238	Thorium 232
Americium 241	Plutonium 239