

TECHNIQUE AND EQUIPMENT PITFALLS IN SPIROMETRY TESTING: SERIOUS THREATS TO YOUR RESPIRATORY SURVEILLANCE PROGRAM

**NIOSH Medical Surveillance Workshop
November 8, 2001**

OUTLINE

- **SPIROMETRY EVALUATION IN THE OCCUPATIONAL SETTING**
- **DATA LIMITATIONS OF MANY SPIROMETERS**
- **ERRORS THAT INFLATE TEST RESULTS**
- **ERRORS THAT REDUCE TEST RESULTS**
- **EXAMPLES**

Mary C. Townsend, Dr.P.H.

M.C. Townsend Associates

Adjunct Assistant Professor
Graduate School of Public Health
University of Pittsburgh
Pittsburgh, PA
mctownsend@aol.com

SPIROMETRY EVALUATION IN THE OCCUPATIONAL SETTING

1. Compare one test result with reference value

- Above cut-off level required by employer to perform job or wear respirator? **; or
- Normal or impaired? - periodic medical assessment

** Volumes may be under-recorded, especially for those with large lungs – test may focus on exceeding cut-off. Graphs usually saved – unfortunately may become next year’s “baseline.”

2. Examine multiple test results for rate of change over time

Technical errors cause false positives and negatives for both types of evaluation

DATA LIMITATIONS OF MANY SPIROMETERS

1. Many save only 3 “best” curves (highest sum of FEV1 + FVC); and
 2. Most report largest recorded FVC and FEV1, per ATS and OSHA.
- So.... Falsely inflated values will be reported, rather than lower accurate values.***
3. Often must delete curve so it won’t be saved, sometimes with no real-time graph;
 4. If saved, often cannot delete curve unless delete test session and start new test.

So.... Need to know that technical errors can inflate results, and Recognize which curves to delete.

ERRORS THAT INFLATE TEST RESULTS


1. Poor testing technique:

- a) Extra breath through the nose (Figure 1),
- b) Slightly submaximal expiratory effort (negative effort dependence of the FEV1) (Figure 2), and
- c) Accept/save curve with large hesitation, even when flagged by spirometer (Figure 3).

2. Flow-type spirometer malfunctions during subject testing, even though calibration successfully set or checked earlier in the day:


- a) Inaccurate zeroing of sensor, performed before each expiration (Figure 4); or
- b) Sensor characteristics change between expirations due to warming, deposition of mucus, or condensation of water vapor (Figure 5).

Figure 1: Extra Breath




| | FVC (L) | % Pred | FEV1/FVC% | % Pred |
|--------------|---------|--------|-----------|--------|
| Good Effort | 4.78 | 96 | 74 | 93 |
| Extra Breath | 5.96 | 119 | 59 | 75 |

Figure 2: Slightly Submaximal Expiratory Effort (Negative Effort Dependence of FEV1)



| | FEV1 | % Pred | FEV1/FVC | % Pred |
|-------------|------|--------|----------|--------|
| Good Effort | 2.90 | 99 | 74 | 90 |
| Submax Push | 3.13 | 106 | 79 | 96 |

Figure 3: Excessive Hesitation



| | FEV1 | % Pred | FEV1/FVC5 | % Pred |
|------------------|------|--------|-----------|--------|
| Good Effort | 3.23 | 109 | 82 | 101 |
| Large Hesitation | 3.34 | 112 | 83 | 102 |

Figure 4: Zero Flow Error: Screen Pneumotach





Figure 5: Condensation/ Mucus Deposition: Screen Pneumotach



ERRORS THAT *REDUCE* TEST RESULTS

1. Leaks in volume spirometer or breathing tubes:

- a) Reduce FVCs significantly but are *not visible in spiograms* until leak is very large (Figure 6);
- b) Checking for leaks at least daily in the calibration check is essential.

2. Errors in testing technique (usually visible in tracing):

- a) Small inspiration (Figure 7-2)
- b) Weak push (Figure 7-3)
- c) Tongue in mouthpiece (Figure 7-4)
- d) Early termination (Figure 7-6)
- e) Glottis Closure (Figure 7-7)

3. If technician corrects test error, reduced values will be replaced by higher, more accurate, results.

Figure 6: Effects of Spirometer Leaks on Test Results




Figure 7-1. GOOD EFFORT




Figure 7-2. SMALL INSPIRATION




Figure 7-3. WEAK PUSH




Figure 7-4. TONGUE IN MOUTHPIECE




Figure 7-5. GOOD EFFORT




Figure 7-6. EARLY TERMINATION






Figure 7-7. GLOTTIS CLOSURE



Example 1: FAILURE OF MEDICAL SCREENING FOR RESPIRATOR USE

- Spirometries conducted on 62 employees in a manufacturing plant.
- Employees were aged 20 - 62, median = 41.6 yr, 63% current smokers, 84 % male.
- 50% of employees tested had both FEV1 and FVC > 120 % of Predicted.
- Spirometer was frequently calibrated, plot of FVC % predicted shown below. Some probable zero flow errors, some possible moisture condensation; spirometer not holding its calibration.
- Test results are uninformative.

FVC % Pred by Date and Time of Test



EXAMPLE 2: FAILURE OF A MEDICAL SURVEILLANCE PROGRAM

Management Summary Report
Employee Spirometry Testing
12/31/99

| | | |
|--|----|--------|
| Employees receiving annual spirometry test | 31 | (100%) |
| • Annual PFT values in normal range | 5 | (15%) |
| • Annual PFT values abnormal | 1 | (3%) |
| • Annual test FVC or FEV1 decreased more than 15% in past year | 20 | (65%) |