

FDA-Approved Rapid HIV Antibody Screening Tests – Purchasing Details

March 22, 2005

	<u>Price Per Device*</u>	<u>External Controls</u>	<u># of Tests per Case</u>	<u>Catalog numbers</u>	<u>Storage Temperature</u>	<u>Operating Temperature</u>	<u>Shelf life of Test**</u>	<u>Shelf Life of Control**</u>	<u>Total Time Required to Conduct Test***</u>	<u>Window Period for Result Validity****</u>
OraQuick Rapid HIV-1 Antibody Test www.orasure.com	\$14.50	Sold Separately (\$20)	25 or 100	#1001-0052 (25 tests) #1001-0051 (100 tests) #1001-0049 (controls)	2-27°C (tests) 2-8 °C (controls)	15-27°C	8 months	8 months 21 days (after opening)	<5 minutes	20-40 min
									+ 20 min wait time	
OraQuick ADVANCE Rapid HIV-1/2 Antibody Test www.orasure.com	\$17.50	Sold Separately (\$25)	25 or 100	#1001-0079 (25 tests) #1001-0078 (100 tests) #1001-0077 (controls)	2-27°C (tests) 2-8°C (controls)	15-37°C	8 months	8 months 21 days (after opening)	<5 minutes (<10 min. for plasma)	20-40 min
									+ 20 min wait time	
Uni-Gold Recombigen HIV www.unigoldhiv.com	\$15.75	Sold Separately (\$26.25)	20	#1206503 (tests) #1206530 (controls)	2-27 °C (tests) 2-8 °C (controls)	15-27°C	12 months	12 months 1 month (after opening)	<5 minutes	10-12 min
									+ 10 minute wait time	
Reveal G-2 Rapid HIV-1 Antibody Test www.reveal-hiv.com	\$14.00	Included	20 or 60	B1057-6 (20 test kits) B1057-7 (60 test kits)	2-30°C (tests) 2-8°C (controls)	15-27°C	12 months	12 months 7 days (after reconstituting)	3-5 minutes	Result must be read immediately
									No add. wait time	
MultiSpot HIV-1/HIV-2 Rapid Test www.biorad.com	\$25.00	Included	50	#72269 (test kits)	2-8 °C Or 20-30 °C (tests and controls)	20-30°C	6-7 months (@ 2-8 °C) 3 months (@ 20-30 °C)	6-7 months	10-15 minutes	Can be read immediately or anytime up to 24 hours
									No add. wait time	

*Actual price may vary by purchasing agreements with manufacturers

** From date of manufacture, unless otherwise noted

***First time listed is estimated time required to set up test. The second time is the required wait time before reading results. Times listed exclude time needed to draw/obtain sample

**** As measured from last step of testing process

Note. Trade names are for identification purposes only and do not imply endorsement. This information was compiled from package inserts and direct calls to manufacturers.

Prepared by: Kali Stanger & Frances Margolin at HRET; Bernard Branson, Jill Clark and Margaret Lampe at CDC