



Newborn Screening Quality Assurance Program

PROFICIENCY TESTING PROGRAM FOR ANTI-HIV-1 IN DRIED BLOOD SPOTS

Quarterly Report

Quarter 2

May 2006

INTRODUCTION

The Anti-HIV-1 proficiency testing (PT) panel for Quarter 2 consisted of five individual-matrix dried-blood spot (DBS) specimens representing a variety of serostatuses. For the second quarter of 2006, we received data reports for 43 anti-HIV-1 PT panels from 39 participants. This report is the outcome of data reported for the anti-HIV-1 PT specimens that were received from participants by the designated deadline date. This quarterly report is distributed to all participants and to program colleagues by request. Each participant is asked to analyze the specimens for anti-HIV-1 using the assay schemes they routinely use and to report for each specimen the screening results along with results from any confirmatory assays performed for reactive screens.

PARTICIPANTS' RESULTS

Participants conducting screening analysis reported no misclassifications for this quarter. All but one laboratory conducting confirmatory analyses identified the reactive specimen 2461 correctly. The one laboratory reported specimen 2461 as indeterminate. In Part 1 of the report, Table 1 shows the number of laboratories using each screening method/kit; Table 2 shows the expected results and the number of laboratories reporting reactive and non-reactive results; and Table 3 summarizes EIA absorbance ranges and means

by specimen number and method of analysis (representative data are shown for three kits). Table 3 also shows the results obtained by CDC.

In Part 2 of the report, Table 4 shows the number of laboratories using each confirmatory method/kit; Table 5 shows the expected results, number of laboratories reporting reactive, non-reactive, or indeterminate results; and the number of laboratories not reporting confirmatory results. Table 6 shows the band classifications by method of PT specimens that tested positive for HIV-1 initially and after repeat EIA analysis (representative data are shown for three methods). Table 6 also shows the results obtained by CDC. ♦

The Quality Assurance Program will ship next quarter's HIV-1 DBS proficiency testing specimens on July 17, 2006, and the next major allotment of HIV-1 DBS quality control specimens on July 17, 2006. ♦

CONFERENCES

1st Annual African HIV/AIDS Clinical Update and Leadership Development Conference 15 to 16 June 2006. Panafric Hotel, Nairobi, Kenya. website: <http://www.valleyaids.org>

The XVI International AIDS Conference will be held August 13-18, 2006, in Toronto, Ontario, Canada.

It is the world's largest, most comprehensive HIV/AIDS conference.

The theme will be "Time to Deliver."

Website: <http://www.aids2006.org>

SpotLight

Training opportunity

Assuring the Quality of HIV Prevention Counseling: Practical Approaches for Supervisors
City: Philadelphia, PA, August 1-3, 2006 or St. Louis, MO September 12-14, 2006 Sponsor: Centers for Disease Control & Prevention (CDC), Division of HIV/AIDS, Prevention Capacity Building Branch.

Contact: For additional information or to register access the Website at www.cdc.gov/hiv/cba/ or contact Theresa Folsom by phone: 404-639-0982. ♦

Quarterly publication for colleagues and participants of the Performance Evaluation Program for Anti-HIV-1 in Dried Blood Spots.

Editorial and Distribution and Inquiries:
Centers for Disease Control and Prevention
(CDC)

Newborn Screening Quality
Assurance Program
4770 Buford Highway NE, MS/F43
Atlanta, GA 30341-3724

Phone: 770-488-7897
Fax: 770-488-4255
E-mail: NMeredith@cdc.gov
Editor: Nancy K. Meredith
Production: Connie Singleton



DEPARTMENT OF HEALTH AND HUMAN SERVICES
CENTERS FOR DISEASE CONTROL AND PREVENTION
SAFER • HEALTHIER • PEOPLE



Quarter 2, 2006

PART 1. SCREENING

Table 1.
Screening Methods

Kit Source	Total Participants
Genetic Systems rLAV EIA (Bio-Rad)	9
bioMerieux Vironostika HIV-1 Microelisa System	11
bioMerieux Vironostika Uni-Form II <i>plus</i> O	7
Fujirebio Serodia-HIV	2
Abbott Murex HIV 1+2 Gacelisa	1
Murex HIV 1.2.0	3
Tecnosuma (Cuba) UMELISA HIV 1+2	0
Other	9
In House	1
Total	43

Table 2.
Number of Labs Reporting Screening Results

Specimen Number	Expected Results	Reactive	Non-reactive	No Interpretation
2641	Reactive	43	0	0
2642	Non-reactive	0	43	0
2643	Non-reactive	0	43	0
2644	Non-reactive	0	43	0
2645	Non-reactive	0	43	0

Table 3.
Kit Source

		Specimen Number					Cutoff Value
		2641	2642	2643	2644	2645	
Genetic Systems CDC Results	Range	1.459-1.876	0.044-0.163	0.035-0.110	0.036-0.129	0.058-0.219	0.268-0.306
	Mean	1.698	0.104	0.079	0.096	0.130	0.280
		1.812	0.167	0.155	0.225	0.268	0.297
bioMerieux Vironostika CDC Results	Range	2.287-2.739	0.270-0.378	0.226-0.298	0.254-0.350	0.272-0.395	0.464-0.519
	Mean	2.287	0.270	0.226	0.254	0.272	0.464
		2.244	0.235	0.195	0.200	0.223	0.429
bioMerieux Vironostika Uni-Form II <i>plus</i> O	Range	1.003-1.945	0.057-0.158	0.055-0.158	0.051-0.197	0.059-0.178	0.148-0.207
	Mean	1.346	0.093	0.096	0.108	0.111	0.183

Use of trade names and commercial sources is for identification only and does not imply endorsement by the U.S. Department of Health and Human Services.

Quarter 2, 2006

PART 2. CONFIRMATORY

Table 4.
Confirmatory Methods

Kit Source	Total Participants
Genetic Systems HIV-1 WB (Bio-Rad)	15
Bio-Rad New LAV Blot I	3
OraSure HIV-1 WB Kit	2
Genelab Diagnostics HIV 2.2 WB	1
Cambridge Biotech HIV-1 WB Kit (Calypte)	1
In House	1
Other	1
Total	24

Table 5.
Number of Labs Reporting Western Blot Results

Specimen Number	Expected Results	Reactive	Non-reactive	Indeterminate	Not Tested
2641	Reactive	22	0	1	0
2642	Non-reactive	0	7	0	17
2643	Non-reactive	0	7	0	17
2644	Non-reactive	0	7	0	17
2645	Non-reactive	0	7	0	17

Table 6.

Specimen 2641

Methods	gp160		gp120		p66		p55		p51		gp41		p31		p24		p18	
	+	-	+	-	+	-	+	-	+	-	+	-	+	-	+	-	+	-
CDC (GS HIV-1 Kit, Bio-Rad)	✓		✓		✓		✓		✓		✓		✓		✓		✓	
Genetic Systems HIV-1 Kit (Bio-Rad)	15	0	15	0	9	6	13	2	10	5	12	3	5	10	15	0	13	2
Bio-Rad New LAV Blot I	2	1	1	2	1	2	3	0	1	2	1	2	1	2	2	1	3	0
OraSure HIV-1 WB Kit	2	0	2	0	1	1	0	2	1	1	2	0	0	2	2	0	2	0

This **NEWBORN SCREENING QUALITY ASSURANCE PROGRAM** report is an internal publication distributed to program participants and selected program colleagues. The laboratory quality assurance program is a project cosponsored by the Centers for Disease Control and Prevention (CDC) and the Association of Public Health Laboratories.

CENTERS FOR DISEASE CONTROL AND PREVENTION (CDC)
ATLANTA, GA 30341

Director

Julie Louise Gerberding, M.D., M.P.H.

Director

National Center for Environmental Health

Howard Frumkin, M.D., Dr.P.H., M.P.H.

Director

Division of Laboratory Sciences

Eric J. Sampson, Ph.D.

Chief

Newborn Screening Branch

W. Harry Hannon, Ph.D.



Contributors:

Barbara W. Adam
Carol Bell
Paul Dantonio
Marie C. Earley, Ph.D.
F. Hugh Gardner
Sherri Hall
L. Omar Henderson, Ph.D.
Lisa Kalman, Ph.D.
Lixia Li, Ph.D.
Timothy Lim, Ph.D.
Elizabeth McCown
Joanne Mei, Ph.D.
Nancy Meredith
Hetal Patel
Jarad Schiffer
Anand Swamy, Ph.D.
Robert Vogt, Ph.D.
Yingtao Zhou

Production:

Sarah Brown
Felicia Manning
Connie Singleton

ASSOCIATION OF PUBLIC HEALTH LABORATORIES
WASHINGTON, DC 20036-3320

President

Katherine Kelley, Dr.P.H.

Chairman, Newborn Screening and Genetics in Public Health Committee

William Becker, D.O., M.P.H.

Chairman, Newborn Screening Quality Assurance Subcommittee

John Sherwin, Ph.D.



INQUIRIES TO:

Carol Bell, Editor • Centers for Disease Control and Prevention (CDC)
Newborn Screening Quality Assurance Program • Mailstop F-43
4770 Buford Highway, N.E. • Atlanta, GA 30341-3724
Phone (770) 488-4582 • FAX (770) 488-4255 • E-mail: CBell@cdc.gov