


CDC Workplace
Safety and Health


NIOSH


NPPTL Research to Practice
through Partnerships


Click Here


Permeation Calculator Version 2.4

















Permeation Calculator

File Help

Import Data View Graph Modify Settings Save Print

Choice of Variable

* All active fields are required data fields unless noted as optional

Manually Select Data Columns

Analyzer Response Format

- ☒ Option 1: Use Concentration (in $\mu\text{g/L}$)
- ☐ Option 2: Use Concentration (in ppm)
- ☐ Option 3: Use Other Analyzer Output Reading

Time Format

- ☒ Time in Minutes
- ☐ YYYY/MM/DD HH:MM:SS
- ☐ MM/DD/YYYY HH:MM:SS ##

Choose System Type

- ☐ Open Loop System (OL)
 - ☒ Constant Flow Rate of Fresh Collection Medium (F in ASTM F 739):
 L / min
Analytical Method Detection Limit: $\mu\text{g/mL}$ (optional)
 - ☐ Variable Flow Rate:
Minimum detectable permeation rate: $\mu\text{g}/(\text{cm}^2 \cdot \text{min})$
- ☒ Closed Loop System (CL)
 - Total Volume of the Collection Medium (V_t in ASTM F 739):
 L
 - ☒ Continuous Sampling
 - ☐ Discrete Sampling
 - ☒ Sample Volume NOT replaced, enter Volume Removed (V_s in ASTM F 739)
 L
 - ☐ Sample Volume IS replaced, enter Volume Replaced (V_s in ASTM F 739)
 L
 - Minimum detectable mass permeated: $\mu\text{g}/\text{cm}^2$

Cancel Next

Permeation Calculator

File Help

Import Data View Graph Modify Settings Save Print

Choice of Variable

* All active fields are required data fields unless noted as optional

Manually Select Data Columns

Choose System Type

☐ Open Loop System (OL)

☒ Constant Flow Rate of Fresh Collection Medium (F in ASTM F 739):

Analyzer Response Format

☒ Option 1: Use Concentration (in

☐ Option 2: Use Concentration (in

☐ Option 3: Use Other Analyzer Output Reading

Select Data Columns

Analyzer Output: Column

Time: Column

Flow Rate (optional) Column

OK Clear Cancel

Time Format

☒ Time in Minutes

☐ YYYY/MM/DD HH:MM:SS

☐ MM/DD/YYYY HH:MM:SS ##

☐ Sample Volume IS replaced, enter Volume Replaced (Vs in ASTM F 739)

L

Minimum detectable mass permeated: $\mu\text{g}/\text{cm}^2$

Cancel Next

$\mu\text{g}/\text{mL}$ (optional)

$\mu\text{g}/(\text{cm}^2 \cdot \text{min})$

ASTM F 739):

Volume Removed (Vs in ASTM F 739)

Permeation Calculator

File Help

Import Data View Graph Modify Settings Save Print

Choice of Variable

* All active fields are required data fields unless noted as optional

Manual Data Columns Selected

Analyzer Response Format

- ☒ Option 1: Use Concentration (in $\mu\text{g/L}$)
- ☐ Option 2: Use Concentration (in ppm)
- ☐ Option 3: Use Other Analyzer Output Reading

Time Format

- ☒ Time in Minutes
- ☐ YYYY/MM/DD HH:MM:SS
- ☐ MM/DD/YYYY HH:MM:SS ##

Choose System Type

☐ Open Loop System (OL)

- ☒ Constant Flow Rate of Fresh Collection Medium (F in ASTM F 739):

 L / min

 Analytical Method Detection Limit: $\mu\text{g/mL}$ (optional)
- ☐ Variable Flow Rate:

 Minimum detectable permeation rate: $\mu\text{g}/(\text{cm}^2\cdot\text{min})$

☒ Closed Loop System (CL)

Total Volume of the Collection Medium (V_t in ASTM F 739):

L

☒ Continuous Sampling

☐ Discrete Sampling

- ☒ Sample Volume NOT replaced, enter Volume Removed (V_s in ASTM F 739)

 L
- ☐ Sample Volume IS replaced, enter Volume Replaced (V_s in ASTM F 739)

 L

Minimum detectable mass permeated: $\mu\text{g}/\text{cm}^2$

Cancel Next

Permeation Calculator

File Help

Import Data View Graph Modify Settings Save Print

Data Input

* All active fields are required data fields

Cumulative Permeation vs Time
Cumulative Permeation ($\mu\text{g}/\text{cm}^2$)

Switch Exposure Size (for A in ASTM F 739)

☒ Diameter: 1.00 ☒ inch
☐ Area: 5.07 ☐ cm

Specimen Weight

1.00 ☒ grams
☐ grams/m²

Cumulative Permeation for: 60 min

Cumulative Permeation Mass target:
150 $\mu\text{g}/\text{cm}^2$


Enter times T_{i-1}: 0.00 and T_i: 120.83
(in Minutes) for Average Permeation Rate

Time (Minutes)
View Data Graph

☐ Steady State Permeation was not established
(Program will report the maximum Permeation Rate)

Back Cancel Next

0



Permeation Calculator

File Help

Import Data View Graph Modify Settings Save Print

Additional Data Input

* All fields are optional data fields (values entered here will not affect the results)

Report Title: Acetone Against Neoprene

Date: 4/3/2007 MM/DD/YYYY

Material Type (Manufacture/Product): Neoprene

Average Material Thickness (mm): 0.685

Chemical Name: 100% Acetone

Physical State:
☒ Liquid
☐ Gas

Nominal Test Temperature: 23.5 degree Celsius

Temperature Range: 22.2 to 25.1 degree Celsius

Comments: This is to compare decontamination methods.

Project Number: PR-1234

Operator: Jane Doe

Experiment Setting

Instrument Type (e.g., MIRAN IR, GC, etc.): Miran IR

Instrument Settings: Wavelength 8.5, Pathlength 2l

Collection Medium: air

Instrument ID Number: CDC 1236

Pump ID Number: Wr-156p

Data Sampling Interval (second): 3

Back Cancel Finish

Permeation Calculator

File Help

Import Data View Graph Modify Settings Save Print

Report Title: Acetone Against Neoprene

Project Number: PR-1234

Experiment type: Closed Loop, Continuous Sampling.

--Results--

Breakthrough Time--

- ASTM F739 Standardized Breakthrough Time (permeation mass at $0.25 \mu\text{g}/\text{cm}^2$): 35.8 min
- Breakthrough Detection Time (BDT): 36.06489 min
- ISO 6529 Normalized Breakthrough Time (permeation mass at $2.5 \mu\text{g}/\text{cm}^2$): 38.9 min


Steady-State Permeation Rate (SSPR)--

- SSPR: $2.74 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$
- Correlation Factor (R^2) in the steady-state region from 93.51667 to 108.8667 min: 0.999718

Cumulative Permeation--

- Elapsed Time for Cumulative Permeation Mass of $150 \mu\text{g}/\text{cm}^2$: 98.3 minutes
- Cumulative Permeation for 60 minutes: $39.9 \mu\text{g}/\text{cm}^2$
- Average Permeation Rate from 0.00 to 120.83 min: $1.73 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$

Operator: Jane Doe



Permeation Calculator

File Help

Import Data View Graph **Modify Settings** Save Print

Report Title: Acetone Against Neoprene

Project Number: PR-1234

Experiment type: Closed Loop, Continuous Sampling.

--Results--

Breakthrough Time--

- ASTM F739 Standardized Breakthrough Time (permeation mass at $0.25 \mu\text{g}/\text{cm}^2$): 35.8 min
- Breakthrough Detection Time (BDT): 36.06489 min
- ISO 6529 Normalized Breakthrough Time (permeation mass at $2.5 \mu\text{g}/\text{cm}^2$): 38.9 min

Steady-State Permeation Rate (SSPR)--

- SSPR: $2.74 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$
- Correlation Factor (R^2) in the steady-state region from 93.51667 to 108.8667 min: 0.999718

Cumulative Permeation--

- Elapsed Time for Cumulative Permeation Mass of $150 \mu\text{g}/\text{cm}^2$: 98.3 minutes
- Cumulative Permeation for 60 minutes: $39.9 \mu\text{g}/\text{cm}^2$
- Average Permeation Rate from 0.00 to 120.83 min: $1.73 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$

Operator: Jane Doe

Permeation Calculator

File Help

Import Data View Graph Modify Settings Save Print

Report Title: Acetone Against Neoprene

Choice of Variable

* All active fields are required data fields unless noted as optional

Manual Data Columns Selected

Analyzer Response Format

- ☒ Option 1: Use Concentration (in $\mu\text{g/L}$)
- ☐ Option 2: Use Concentration (in ppm)
- ☐ Option 3: Use Other Analyzer Output Reading


Time Format

- ☒ Time in Minutes
- ☐ YYYY/MM/DD HH:MM:SS
- ☐ MM/DD/YYYY HH:MM:SS ##

Choose System Type

- ☐ Open Loop System (OL)
 - ☒ Constant Flow Rate of Fresh Collection Medium (F in ASTM F 739):
 L / min
Analytical Method Detection Limit: $\mu\text{g/mL}$ (optional)
 - ☐ Variable Flow Rate:
Minimum detectable permeation rate: $\mu\text{g}/(\text{cm}^2 \cdot \text{min})$
- ☒ Closed Loop System (CL)
 - Total Volume of the Collection Medium (V_t in ASTM F 739):
 L
 - ☒ Continuous Sampling
 - ☐ Discrete Sampling
 - ☒ Sample Volume NOT replaced, enter Volume Removed (V_s in ASTM F 739)
 L
 - ☐ Sample Volume IS replaced, enter Volume Replaced (V_s in ASTM F 739)
 L
 - Minimum detectable mass permeated: $\mu\text{g}/\text{cm}^2$

Cancel Next



Permeation Calculator

File Help

Import Data View Graph **Modify Settings** Save Print

Report Title: Acetone Against Neoprene

Project Number: PR-1234

Experiment type: Closed Loop, Continuous Sampling.

--Results--

Breakthrough Time--

- ASTM F739 Standardized Breakthrough Time (permeation mass at $0.25 \mu\text{g}/\text{cm}^2$): 35.8 min
- Breakthrough Detection Time (BDT): 36.06489 min
- ISO 6529 Normalized Breakthrough Time (permeation mass at $2.5 \mu\text{g}/\text{cm}^2$): 38.9 min


Steady-State Permeation Rate (SSPR)--


- SSPR: $2.74 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$
- Correlation Factor (R^2) in the steady-state region from 93.51667 to 108.8667 min: 0.999718


Cumulative Permeation--


- Elapsed Time for Cumulative Permeation Mass of $150 \mu\text{g}/\text{cm}^2$: 98.3 minutes
- Cumulative Permeation for 60 minutes: $39.9 \mu\text{g}/\text{cm}^2$
- Average Permeation Rate from 0.00 to 120.83 min: $1.73 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$


Operator: Jane Doe











Permeation Calculator

File Help

Import Data View Graph Modify Settings **Save** Print

Report Title: Acetone Against Neoprene

Project Number: PR-1234

Experiment type: Closed Loop, Continuous Sampling.

--Results--

Breakthrough Time--

- ASTM F739 Standardized Breakthrough Time (permeation mass at $0.25 \mu\text{g}/\text{cm}^2$): 35.8 min
- Breakthrough Detection Time (BDT): 36.06489 min
- ISO 6529 Normalized Breakthrough Time (permeation mass at $2.5 \mu\text{g}/\text{cm}^2$): 38.9 min

Steady-State Permeation Rate (SSPR)--

- SSPR: $2.74 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$
- Correlation Factor (R^2) in the steady-state region from 93.51667 to 108.8667 min: 0.999718

Cumulative Permeation--

- Elapsed Time for Cumulative Permeation Mass of $150 \mu\text{g}/\text{cm}^2$: 98.3 minutes
- Cumulative Permeation for 60 minutes: $39.9 \mu\text{g}/\text{cm}^2$
- Average Permeation Rate from 0.00 to 120.83 min: $1.73 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$

Operator: Jane Doe

Permeation Calculator

File Help

Import Data View Graph Modify Settings Save Print

Report Title: Acetone Against Neoprene

Project Number: PR-1234

Experiment type: Closed Loop, Continuous Sampling.

--Results--

Breakthrough Time--
ASTM F739 Standard
Breakthrough Detection
ISO 6529 Normalized

Steady-State Permeation Rate
SSPR: 2.74 $\mu\text{g}/(\text{cm}^2 \cdot \text{min})$
Correlation Factor (R^2) in the steady-state region from 93.51667 to 108.8667 min: 0.999718

Cumulative Permeation--
Elapsed Time for Cumulative Permeation Mass of 150 $\mu\text{g}/\text{cm}^2$: 98.3 minutes
Cumulative Permeation for 60 minutes: 39.9 $\mu\text{g}/\text{cm}^2$
Average Permeation Rate from 0.00 to 120.83 min: 1.73 $\mu\text{g}/(\text{cm}^2 \cdot \text{min})$

Operator: Jane Doe


Save File As


Save

☒ Save As an Excel File
☐ Save As Text File

OK Cancel

0.25 $\mu\text{g}/\text{cm}^2$: 35.8 min
g/cm²): 38.9 min





Microsoft Excel - Results for Data File 1

File Edit View Insert Format Tools Data Window Help

Type a question for help

MS Sans Serif 10 B I U \$ % , .00 .00

Reply with Changes... End Review...

H2O

	A	B	C	D	E	F	G	H	I	J	K	L
1	Permeation Calculator	Report Title:	Acetone Against Neoprene	Result:	--							
2												
3	Operator:	Jane Doe										
4	Date:	4/3/2007										
5	Data Filename:	Data File 1.xls										
6	Project Number:	PR-1234										
7		ASTM F739 Standardized Breakthrough Time (permeation mazz at 0.25 $\mu\text{g}/\text{cm}^2$):		35.8 min								
8		Breakthrough Detection Time (BDT):		36.06489 min								
9		ISO 6529 Normalized Breakthrough Time (permeation mazz at 2.5 $\mu\text{g}/\text{cm}^2$):		38.9 min								
10		Steady-State Permeation Rate (SSPR):		2.74 $\mu\text{g}/(\text{cm}^2\cdot\text{min})$								
11		Correlation Factor (R ²) in the Steady-State Region:		0.999718								
12		Start Time in the Steady-State Region:		93.51667 min								
13		End Time in the Steady-State Region:		108.8667 min								
14		Elapsed Time for Cumulative Permeation Mazz of 150 $\mu\text{g}/\text{cm}^2$:		98.3 minutes								
15		Cumulative Permeation for 60 minutes:		39.9 $\mu\text{g}/\text{cm}^2$								
16		Average Permeation Rate from 0.00 to 120.83 minutes:		1.73 $\mu\text{g}/(\text{cm}^2\cdot\text{min})$								
17												
18												
19	Experiment Information											
20		Test Duration:		2.01 hours								
21	Material	Manufacturer & Product:		Neoprene								
22		Average Thickness:		0.635 mm								
23		Exposure Area:		5.07 cm^2								
24		Weight Per Unit Area of Specimen:		1970 g/m^2								
25	Chemical(s)	Physical State:		Liquid								
26		Components with Concentration (%):		100% Acetone								
27	Temperature	Nominal Test:		23.5 degrees Celsius								
28		Range:		22.2 to 25.1 degrees Celsius								
29	Analytical Technique	Instrument Type:		Miran IR								
30		Instrument ID Number:		ODC 1236								
31		Instrument Settings:		Wavelength 8.5, Pathlength 20.5								
32		Sampling Pump ID:		W-156p								
33	Collection System	Medium:		air								
34		Total Volume of the Collection Medium (Vt):		5.64 L								
35	System Type:	Closed Loop, Continuous Sampling:										
36		Data Sampling Interval (seconds):		3								
37	Additional Comment:	This is to compare decontamination methods.										
38												

Results

Draw AutoShapes

Ready

NUM

For open loop testing under a constant flow rate, there is an option to enter a value for the “Analytical Method Detection Limit” to calculate the “Minimum Detectable Permeation Rate” (see the report on the next slide)

The screenshot shows the 'Permeation Calculator' software window. The 'Report Title' is 'Neoprene - Acetone'. A 'Choice of Variable' dialog box is open, allowing the user to select system type and data format. The dialog includes sections for 'Manually Select Data Columns', 'Analyzer Response Format', 'Time Format', and 'Choose System Type'. The 'Open Loop System (OL)' is selected, with 'Constant Flow Rate of Fresh Collection Medium' chosen. The flow rate is 3.94 L/min, and the 'Analytical Method Detection Limit' is set to 1000 µg/mL. The 'Closed Loop System (CL)' is also an option, with 'Continuous Sampling' and 'Sample Volume NOT replaced' selected. The 'Next' button is highlighted with a red arrow.

Permeation Calculator

File Help

Import Data View Graph Modify Settings Save Print

Report Title: Neoprene - Acetone

Choice of Variable

* All active fields are required data fields unless noted as optional

Manually Select Data Columns

Analyzer Response Format

- ☒ Option 1: Use Concentration (in µg/L)
- ☐ Option 2: Use Concentration (in ppm)
- ☐ Option 3: Use Other Analyzer Output Reading

Time Format

- ☒ Time in Minutes
- ☐ YYYY/MM/DD HH:MM:SS
- ☐ MM/DD/YYYY HH:MM:SS ##

Choose System Type

- ☒ Open Loop System (OL)
 - ☒ Constant Flow Rate of Fresh Collection Medium (F in ASTM F 739):
3.94 L / min
Analytical Method Detection Limit: 1000 µg/mL (optional)
 - ☐ Variable Flow Rate.
Minimum detectable permeation rate: 0.10 µg/(cm²*min)
- ☐ Closed Loop System (CL)
 - Total Volume of the Collection Medium (Vt in ASTM F 739):
5.64 L
 - ☒ Continuous Sampling
 - ☐ Discrete Sampling
 - ☒ Sample Volume NOT replaced, enter Volume Removed (Vs in ASTM F 739)
0.05 L
 - ☐ Sample Volume IS replaced, enter Volume Replaced (Vs in ASTM F 739)
0.05 L

Minimum detectable mass permeated: 0.25 µg/cm²

Cancel Next

Permeation Calculator

File Help

Import Data View Graph Modify Settings Save Print

Project Number: PR-1234

Experiment type: Open Loop, Constant Flow Rate.

--Results--

Breakthrough Time--

- ASTM F739 Standardized Breakthrough Time (permeation rate at $0.1 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$): 3.07 min
- Breakthrough Detection Time (BDT): 3.287359 min
- ISO 6529 Normalized Breakthrough Time (permeation rate at $1.0 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$): 3.70 min

Steady-State Permeation Rate (SSPR)--

- SSPR: $117 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$
- Determined around 13.15000 (average of: 13.20000; 13.15000; 13.10000) min
- Maximum Permeation Rate: $117 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$
- Determined at 13.20000 min

Cumulative Permeation--

- Elapsed Time for Cumulative Permeation Mass of $150 \mu\text{g}/\text{cm}^2$: 7.76 minutes
- Cumulative Permeation for 13.20 minutes: $726 \mu\text{g}/\text{cm}^2$
- Average Permeation Rate from 0.00 to 13.20 min: $55.4 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$

Minimum detectable permeation rate based on analytical method detection limit: $0.78 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$