# It's NOT What You Think!



An N95 respirator and a surgical mask provide the same protection.



#### 

I can decorate my N95 respirator to look stylish!





### 

There's no way to verify if an N95 respirator is actually NIOSH-approved.

### FACT:

A NIOSH-approved N95 respirator forms a seal against the user's face, preventing particle penetration around the edges. The filter has passed NIOSH tests to determine that it protects against at least 95% of airborne particles.

Because surgical masks do not seal against the face and the filters have not been tested, the same level of protection against airborne particles cannot be guaranteed.

## FACT:



Any alteration of the N95 will void the NIOSH approval because respirator can no longer be guaranteed to provide the necessary level of protection.

#### FACT:





A list of approved N95 respirators is updated every month. Find it at Knowits.niosh.gov under the **"Approved Particulate Filtering** Facepiece Respirators" tab.

