



NCID Organism to SNOMED Mapping Application: A Practical Tool to Achieve Lab Vocabulary's PHIN Compliance

Wenkai Li, James Tolson, Heather Sun, Jyothi Kondamuri, Malini Kesavan, Sheila Abner, Dale Nordenberg

Office of Informatics (OI)

National Center for Infectious Diseases (NCID)

Centers for Disease Control and Prevention (CDC)

SAFER•HEALTHIER•PEOPLE[™]

Presentation Outline

- Introduction
 - NCID STARLIMS implementation
 - SNOMED CT as a national recommended standard
- NCID organism mapping issues
- NCID organism to SNOMED mapping tool
- Summary

Benefits of Implementing Standardized Vocabularies

- Support information sharing/interoperability
 - Compliant with regulations/requirements
- Capture institutional knowledge
- Help to reengineer business process
- Make data analysis easier
 - Document management
 - Knowledge management


STARLIMS NCID Implementation

- STARLIMS is a commercial off the shelf (COTS)
 laboratory information management system
- NCID STARLIMS Implementation
 - Started in 2004
 - An enterprise laboratory information system
 - Adopt PHIN compliant Standards
 - HL7
 - SNOMED CT
 - LOINC


SNOMED CT

(Systemized Nomenclature of Medicine Clinical Terms)

- A comprehensive clinical terminology that used to code, retrieve, and analyze clinical data
- Comprised of concepts, terms and relationships to precisely represent clinical information across the scope of health care
- Concepts are the basic units of SNOMED CT
- SNOMED CT is organized into hierarchies of concepts
- The January 2005 release contains 364,461 concepts, 399,755 preferred terms and 188,859 synonyms
- There are 25,888 organism concepts and 71,246 organism terms in January 2005 release



SNOMED CT Database



SNOMED CT Concepts Table

ConceptID	Concept Status	FullySpecificedName	CTV3ID	SNOMEDID	IsPrimitive
23496000	0	Fungus (organism)	X79pZ	L-40000	1
233604007	0	Pneumonia (disorder)	X100E	D2-0007F	0
274103002	6	Pneumonia NOS (disorder)	XM0rv	D2-0009B	1

SNOMED CT Descriptions Table

DescriptionID	Description Status	ConceptID	Term	Initial capital Status	Description Type	Language Code
753022014	0	23496000	Fungus (organism)	0	3	en
39449010	0	23496000	Fungus	0	1	en
1223107017	0	23496000	Fungi	0	2	en
39451014	0	23496000	Mycete	0	2	en

SNOMED CT Relationships Table

RelationshipID	ConceptID1	Relationship Type	ConceptID2	Characteristic Type	Refinability	Relation- ship Group	
2689692020	23496000	116680003	36272005	0	0	0	
	Fungus	IS-A relationship	Infectious agent				
2689693026	23496000	116680003	264395009	0	0	0	
	Fungus	IS-A relationship	Microorganism				

NCID Organism Mapping Challenges and Solutions

- NCID proprietary organism terms
 - came from legacy systems like LITS Plus
 - Not well formatted
 - Not systematically classified
- SNOMED CT not a prefect standardized vocabulary
- Mapping results reviewed by practicing lab microbiologists/NCID specialty expert group
- Informatics specialists, instead of making decisions, function as process facilitators in the final mapping decision making process

Examples of a Proprietary Term Set Mapping Issues

- Absidia Corymbifera
 - Inferred as a bacterium
 - A fungal concept in SNOMED CT
- Amycolatopsis orentalis
 - In SNOMED CT: amycolatopsis orientalis ss orientalis
- Gordonia rubropertincta
 - In SNOMED CT: Gordona rubropertinctus

Examples of a Proprietary Term Set Mapping Issues

- Disease vs. causative agent: Tuberculosis
 - It is a disease concept in SNOMED CT
 - The causative agent should be mycobacterium tuberculosis
- Duplicated proprietary terms:
 - Nocardia dassonville (synonym in SNOMED CT)
 - Nocardiopsis Dassovillei (preferred term in SNOMED CT)
 - SNOMED CT concept: nocardiopsis dassonvillei


Concept Classification Examples from a NCID Lab:

- Clonorchis/Opisthorchis/Heterophyes/ Metagonimus
 - SNOMED CT: Family Opisthorchiidae??
- Egg
- Charcot-Leyden crystals
- OTHER_specify
- Remark_No Parasites Found
- Remark_Unacceptable Specimen


Supporting Tool is needed for the Organism Mapping Process

- To handle the larger number of terms
- To manage NCID legacy organism terms
- To maximize automation of the mapping process
- To document issues, notes and decisions during the mapping process
 - Discard proprietary terms
 - Submit new terms or concepts to SNOMED CT
- Support researches on the mapping results

Mapping Process Diagram



NCID Organism Mapping Database ER Diagram



Functions of the Mapping Tool

- Bulk load the subset of SNOMED CT organism description terms
- Bulk load and register an organism value set from a NCID laboratory
- Manually review and revise a NCID organism term and document the change as needed
- Automatically map a NCID organism value set against SNOMED CT organism subset


Functions of the Mapping Tool

- Manually map NCID organism terms that can not be automatically matched against the SNOMED CT database
- Check the mapping results
- Track the organism term submission to SNOMED CT
- Assign a interim NCID organism ID to a organism term submitted to SNOMED CT
- Generate a SNOMED mapped organism value set for a particular laboratory module in the NCID StarLIMS implementation


Technical Specifications of the Mapping Tool

- Access through NCID network
 - Limited security accesses by authorized users'
 CDC network user IDs and passwords
- Backend database maintained in SQL server 2000
- Front end interfaces created by MS Access project (.ADP)

Main Page




Register a NCID Lab




View Raw Data

2 N	1icrosoft Access - [dbo_NCID_Raw_Organism_List]								_ 6	l ×
: =8	<u>File E</u> dit <u>I</u> nsert	<u>R</u> ecords <u>W</u> indow <u>H</u> elp					Ту	pe a question f	or help	5	× 4
	Raw_Organism_ID	Local_Organism_Name	Initial_Cleaned_Name	Fin	Lo	Sta	Lab_	Organism_0	Classifica	ation_IC	J_
	123	Achromobacter spp	Achromobacter Species	3			1				
	124	Achromobacter denitrificans					1				
	125	Achromobacter insolitus					1				
	126	Achromobacter (Alcaligenes) piechau	ıdii				1				
	127	Achromobacter ruhlandii					1				
	128	Achromobacter spanius					1				
	129	Achromoabacter xylosoxidans					1				
	130	Achromobacter group B					1				
	131	Achromobacter group E					1				
	132	Acidovorax spp	Acidovorax Species				1				
	133	Acidovorax delafieldii					1				
	134	Acidovorax faecilis					1				
	135	Acidovorax temperans					1				
	136	Acinetobacter spp	Acinetobacter Species				1				
	137	Acinetobacter baumannii					1				
	138	Acinetobacter baylyi					1				
	139	Acinetobacter bouvetii					1				
	140	Acinetobacter calcoaceticus					1				
	141	Acinetobacter gerneri					1				
	142	Acinetobacter grimontii					1				
	143	Acinetobacter haemolyticus					1				
	144	Acinetobacter johnsonii					1				
	145	145 Acinetobacter junii					1				
	146	Acinetobacter lwoffi					1				
	147	Acinetobacter parvus					1				
	148	Acinetobacter radioresistens					1				
	149	Acinetobacter schindleri					1				- 1
Rec	ord: I4 4	1 ▶ ▶I ▶* 🐼 ▶! of 722	<u>'</u>								
	Start 📗 🚮 🏉 💲	I O	8 MainForm : Form	8 db	o_N	CID	_Raw	_Organ	Q	44	EN

Check Duplicated Mappings



Clue Browser View



Keys to Standard Implementations

- Full understanding of an application's functional requirements
- Familiar with characteristics of relevant source standards
- Strong technical expertise
- Willing to play in a multidisciplinary team sport

Summary

- The NCID Organism to SNOMED mapping application is a practical tool to make lab vocabulary compliant with PHIN specifications
 - Provided a centralized management of the organism to SNOMED mapping process
 - Increased the organism to SNOMED mapping efficiency and accuracy
 - Can be adapted with a little efforts to any other domains for SNOMED CT concept mappings
 - Piloting stage
- Further applicant improvements and function enhancements are underway

Acknowledgments

CDC Information Technology Services (CITS)

Electronic Data Systems (EDS)

Lockheed Martin

Northrop Grumman

Emergint

Business Computer Application (BCA)

NCID PHIN/Interoperability Working Group CDC Public Health Informatics Fellowship