

Suggested citation for this article: Diederien BMW, Veenendaal D, Jansen R, Herpers BH, Ligtoet EEJ, IJzerman EPF. Rapid antigen test for pandemic (H1N1) 2009 virus [letter]. *Emerg Infect Dis.* 2010 May; [Epub ahead of print]

Rapid Antigen Test for Pandemic (H1N1) 2009 Virus

To the Editor: Drexler et al. recently compared the sensitivity of the BinaxNOW Influenza A & B Rapid Test (BinaxNOW; Inverness Medical, Cologne, Germany) with that of a real-time reverse transcription–PCR (RT-PCR) assay specific for influenza A pandemic (H1N1) 2009 virus (*1*). Of 1,838 clinical specimens tested, 221 were confirmed as positive for pandemic (H1N1) 2009 by RT-PCR. When 144 of these 221 specimens were evaluated by using the BinaxNOW, results were positive for only 16 (11%).

At onset of the pandemic, we evaluated the first 135 nasopharyngeal aspirates submitted to the Regional Laboratory of Public Health Haarlem, the Netherlands. We compared the performance of the BinaxNOW for diagnosing influenza A (H1N1) virus by using molecular detection of influenza virus as the reference standard. Samples were analyzed with a general influenza A assay targeting the matrix gene (the RespiFinder assay) (PathoFinder B.V., Maastricht, the Netherlands [*2*]) and a pandemic (H1N1) 2009–specific RT-PCR assay targeting the neuraminidase gene (*3*). We tested 135 patient samples (76 from male patients); mean age of patients was 32 years (range 0–81 years). Samples from 38 (28%) patients had positive results in both RT-PCRs, and samples from 97 (72%) patients had negative results in the matrix gene RT-PCR and neuraminidase RT-PCR assays. Sensitivity and specificity were estimated to be 47% (18/38, 95% confidence interval [CI] 32%–62%) and 95% (92/97, 95% CI 88%–98%), respectively, for the BinaxNOW antigen test. Patients' ages did not significantly differ between rapid test–positive and –negative results.

Our results largely agree with those of Vasoo et al. (*4*) and the Centers for Disease Control and Prevention (*5*). Those studies determined that the sensitivity of the BinaxNOW compared with nucleic acid amplification tests is ≈40%. The lower sensitivity observed by Drexler et al. (*1*) might be because of differences in study type (retrospective evaluation

compared with a prospective cohort in our study), sample size, technical factors (with regard to specimen collection, specimen transport, and specimen storage), differences in the test kit, and differences between individual patients (multiple categories of age and stages of illness, differences in virus shedding).

Many clinicians are not aware of the performance of specific test devices and rely on test results to make clinical decisions. Because negative results cannot rule out influenza, this test is of little use in a clinical setting with appreciation of the limitations of the test. However, because the BinaxNOW has reasonable specificity, it might prove useful in clinical or epidemiologic situations in which test sensitivity is not critical, e.g., in facility outbreaks in which multiple specimens are collected to rapidly identify the causative organism.

Bram M.W. Diederén, Dick Veenendaal, Ruud Jansen, Bjorn H. Herpers,
Eric E.J. Ligtoet, and Ed P.F. IJerman

Author affiliation: Regional Laboratory of Public Health Haarlem, Haarlem, the Netherlands

References

1. Drexler JF, Helmer A, Kirberg H, Reber U, Panning M, Müller M, et al. Poor clinical sensitivity of rapid antigen test for influenza A pandemic (H1N1) 2009 virus. *Emerg Infect Dis*. 2009;15:1662–4. [PubMed DOI: 10.3201/eid1502.081028](#)
2. Reijans M, Dingemans G, Klaassen CH, Meis JF, Keijdenier J, Mulders B, et al. RespiFinder: a new multiparameter test to differentially identify fifteen respiratory viruses. *J Clin Microbiol*. 2008;46:1232–40. [PubMed DOI: 10.1128/JCM.02294-07](#)
3. RIVM Laboratory Protocol Library. Influenza A PCR light cycler–probe test A-Matrix-H1-H1v-H3-H5-N1-N1v-N2. Updated June 2009 [cited 2009 Aug 14]. http://www.rivm.nl/cib/binaries/Influenza_diagnostic_qPCR_RIVM_tcm92-61120.pdf
4. Vasoo S, Stevens J, Singh K. Rapid antigen tests for diagnosis of pandemic (swine) influenza A/H1N1. *Clin Infect Dis*. 2009;49:1090–3. [PubMed DOI: 10.1086/644743](#)
5. Centers for Disease Control and Prevention. Evaluation of rapid influenza diagnostic tests for detection of novel influenza A (H1N1) virus—United States. August 7, 2009. *MMWR Morb Mortal Wkly Rep*. 2009;58:826–9. [PubMed](#)

Address for correspondence: B.M.W. Diederén, Regional Laboratory of Public Health Haarlem, Boerhaavelaan 26, 2035 RC Haarlem, the Netherlands; email: bramdiederén@gmail.com

Suggested citation for this article: Drexler JF, Drosten C, Eis-Hübinger AM. Rapid antigen test for pandemic (H1N1) 2009 virus [letter]. *Emerg Infect Dis.* 2010 May; [Epub ahead of print]

In Response: We read with interest the report by Diederer et al. (1) showing a 47% sensitivity of the BinaxNOW (Inverness Medical, Cologne, Germany) antigen-based rapid influenza diagnostic test (RIDT) for the clinical detection of pandemic (H1N1) 2009 virus. We agree that RIDTs may be of little benefit in situations where a timely diagnosis by reverse transcription–PCR (RT-PCR) or optimized direct fluorescent antibody (DFA) tests can be achieved.

Our recent study yielded even lower sensitivity for RIDT: 11.1% (2). RIDT sensitivity is greatly influenced by differences in the level of virus shedding between children and adults, making studies difficult to compare (3). In general, age profiles and virus concentrations should be provided and considered when comparing cohorts examined by any virus detection method. Moreover, quality and origin of specimens can influence the sensitivity of RT-PCR– and antigen-based tests. One important example is the use of flocked swabs for collecting respiratory samples. Under optimal conditions, for instance, a DFA test was recently shown to yield high diagnostic sensitivity comparable with that of RT-PCR for pandemic (H1N1) 2009 virus (4). Another critical factor, especially for RIDT, may be the compatibility of test monoclonal antibodies with the novel virus. Lower sensitivities of such tests for pandemic (H1N1) 2009 virus in comparison with seasonal influenza viruses have been reported (3,5). Adaptation of RIDT antibody selection to pandemic (H1N1) 2009 virus may thus be necessary. Finally, we would like to emphasize the medical risks associated with use of RIDTs by untrained operators, e.g., lesions from inadequate sampling and false interpretation of test results. Such use may be specifically promoted by ready availability of such tests on the Internet or at pharmacies.

Jan Felix Drexler, Christian Drosten, and Anna Maria Eis-Hübinger

Author affiliation: University of Bonn Medical Centre, Bonn, Germany

References

1. Diederer BMW, Veenendaal D, Jansen R, Herpers BH, Ligtvoet EEJ, IJzerman EPF. Rapid antigen test for pandemic (H1N1) 2009 virus [letter]. *Emerg Infect Dis*. 2010 May; [Epub ahead of print].
2. Drexler JF, Helmer A, Kirberg H, Reber U, Panning M, Muller M, et al. Poor clinical sensitivity of rapid antigen test for influenza A pandemic (H1N1) 2009 virus. *Emerg Infect Dis*. 2009;15:1662–4.
3. Centers for Disease Control and Prevention. Evaluation of rapid influenza diagnostic tests for detection of novel influenza A (H1N1) virus—United States, 2009. *MMWR Morb Mortal Wkly Rep*. 2009;58:826–9.
4. Pollock NR, Duong S, Cheng A, Han LL, Smole S, Kirby JE. Ruling out novel H1N1 influenza virus infection with direct fluorescent antigen testing. *Clin Infect Dis*. 2009;49:e66–8.
5. Kok J, Blyth CC, Foo H, Patterson J, Taylor J, McPhie K, et al. Comparison of a rapid antigen test with nucleic acid testing during cocirculation of pandemic influenza A/H1N1 2009 and seasonal influenza A/H3N2. *J Clin Microbiol*. 2010;48:290–1.

Address for correspondence: Christian Drosten, Institute of Virology, University of Bonn Medical Centre, 53127 Bonn, Germany; email: drosten@virology-bonn.de