

**TABLE 2. Tests for hepatitis C virus (HCV) infection**

Test/Type	Application	Comments
<b>Hepatitis C virus antibody (anti-HCV)</b>		
<ul style="list-style-type: none"> <li>EIA (enzyme immunoassay)</li> <li>Supplemental assay (i.e., recombinant immunoblot assay [RIBA™])</li> </ul>	<ul style="list-style-type: none"> <li>Indicates past or present infection, but does not differentiate between acute, chronic, or resolved infection</li> <li>All positive EIA results should be verified with a supplemental assay</li> </ul>	<ul style="list-style-type: none"> <li>Sensitivity <math>\geq 97\%</math></li> <li>EIA alone has low-positive predictive value in low-prevalence populations</li> </ul>
<b>HCV RNA (hepatitis C virus ribonucleic acid)</b>		
<b>Qualitative tests*†</b>		
<ul style="list-style-type: none"> <li>Reverse transcriptase polymerase chain reaction (RT-PCR) amplification of HCV RNA by in-house or commercial assays (e.g., Amplicor HCV™)</li> </ul>	<ul style="list-style-type: none"> <li>Detect presence of circulating HCV RNA</li> <li>Monitor patients on antiviral therapy</li> </ul>	<ul style="list-style-type: none"> <li>Detect virus as early as 1–2 weeks after exposure</li> <li>Detection of HCV RNA during course of infection might be intermittent; a single negative RT-PCR is not conclusive</li> <li>False-positive and false-negative results might occur</li> </ul>
<b>Quantitative tests*†</b>		
<ul style="list-style-type: none"> <li>RT-PCR amplification of HCV RNA by in-house or commercial assays (e.g., Amplicor HCV Monitor™)</li> <li>Branched chain DNA<sup>§</sup> (bDNA) assays (e.g., Quantiplex™ HCV RNA Assay)</li> </ul>	<ul style="list-style-type: none"> <li>Determine concentration of HCV RNA</li> <li>Might be useful for assessing the likelihood of response to antiviral therapy</li> </ul>	<ul style="list-style-type: none"> <li>Less sensitive than qualitative RT-PCR</li> <li>Should not be used to exclude the diagnosis of HCV infection or to determine treatment endpoint</li> </ul>
<b>Genotype*†</b>		
<ul style="list-style-type: none"> <li>Several methodologies available (e.g., hybridization, sequencing)</li> </ul>	<ul style="list-style-type: none"> <li>Group isolates of HCV based on genetic differences, into 6 genotypes and &gt;90 subtypes</li> <li>With new therapies, length of treatment might vary based on genotype</li> </ul>	<ul style="list-style-type: none"> <li>Genotype 1 (subtypes 1a and 1b) most common in United States and associated with lower response to antiviral therapy</li> </ul>
<b>Serotype*</b>		
<ul style="list-style-type: none"> <li>EIA based on immunoreactivity to synthetic peptides (e.g., Murex HCV Serotyping 1–6 Assay)</li> </ul>	<ul style="list-style-type: none"> <li>No clinical utility</li> </ul>	<ul style="list-style-type: none"> <li>Cannot distinguish between subtypes</li> <li>Dual infections often observed</li> </ul>

\* Currently not U.S. Food and Drug Administration approved; lack standardization.

† Samples require special handling (e.g., serum must be separated within 2–4 hours of collection and stored frozen [-20 C or -70 C]; frozen samples should be shipped on dry ice).

§ Deoxyribonucleic acid.