

Newborn Screening Quality Assurance Program

Quality Control Specimen Certification

Set 1—January 2015

Acylcarnitines Method: MS/MS Non-Derivatized - MS/MS non-kit

ENRICHMENT LEVELS (endogenous levels not included)

| <i>Analyte (μmol/L whole blood)</i> | <i>Lot</i> | <i>Base</i> | <i>Lot</i> | <i>Low</i> | <i>Lot</i> | <i>Intermediate</i> | <i>Lot</i> | <i>High</i> |
|--|------------|-------------|------------|------------|------------|---------------------|------------|-------------|
| Free carnitine (C0) | 1465 | 0 | 1466 | 10.0 | 1467 | 20.0 | 1468 | 30.0 |
| Acetylcarnitine (C2) | 1465 | 0 | 1466 | 10.0 | 1467 | 20.0 | 1468 | 30.0 |
| Propionylcarnitine (C3) | 1465 | 0 | 1466 | 3.0 | 1467 | 7.5 | 1468 | 12.0 |
| Butyrylcarnitine (C4) | 1465 | 0 | 1466 | 1.0 | 1467 | 2.5 | 1468 | 5.0 |
| Malonylcarnitine + Hydroxybutyrylcarnitine (C3DC + C4OH) | 1465 | 0 | 1466 | 1.0 | 1467 | 2.5 | 1468 | 5.5 |
| Isovalerylcarnitine (C5) | 1465 | 0 | 1466 | 0.5 | 1467 | 1.5 | 1468 | 3.0 |
| Glutaryl carnitine (C5DC) | 1465 | 0 | 1466 | 0.5 | 1467 | 1.0 | 1468 | 2.5 |
| Hydroxyisovalerylcarnitine (C5OH) | 1465 | 0 | 1466 | 0.5 | 1467 | 1.5 | 1468 | 2.5 |
| Hexanoylcarnitine (C6) | 1465 | 0 | 1466 | 0.5 | 1467 | 1.0 | 1468 | 2.5 |
| Octanoylcarnitine (C8) | 1465 | 0 | 1466 | 0.5 | 1467 | 1.0 | 1468 | 2.5 |
| Decanoylcarnitine (C10) | 1465 | 0 | 1466 | 0.5 | 1467 | 1.0 | 1468 | 2.5 |
| Dodecanoylcarnitine (C12) | 1465 | 0 | 1466 | 0.5 | 1467 | 1.0 | 1468 | 2.5 |
| Myristoylcarnitine (C14) | 1465 | 0 | 1466 | 0.5 | 1467 | 1.5 | 1468 | 3.0 |
| Palmitoylcarnitine (C16) | 1465 | 0 | 1466 | 3.0 | 1467 | 8.0 | 1468 | 12.0 |
| Hydroxypalmitoylcarnitine (C16OH) | 1465 | 0 | 1466 | 0.25 | 1467 | 0.5 | 1468 | 1.0 |
| Stearoylcarnitine (C18) | 1465 | 0 | 1466 | 1.0 | 1467 | 2.0 | 1468 | 5.0 |
| Hydroxystearoylcarnitine (C18OH) | 1465 | 0 | 1466 | 0.5 | 1467 | 1.0 | 1468 | 1.5 |

Newborn Screening Quality Assurance Program

Quality Control Specimen Certification

Set 1— January 2015

Acylcarnitines Method: MS/MS Non-Derivatized - MS/MS non-kit

ANALYTICAL INFORMATION Lot Numbers, Mean Values (\bar{x} , $\mu\text{mol/L}$ whole blood), and 95% Confidence Limits (CL)

| Analyte | Lot | Mean/ 95% CL | Lot | Mean/ 95% CL | Lot | Mean/ 95% CL | Lot | Mean/ 95% CL |
|----------------|------|------------------------------------|------|------------------------------------|------|------------------------------------|------|------------------------------------|
| C0 | 1465 | $\bar{x} = 13.6$ CL = 11.7–15.4 | 1466 | $\bar{x} = 23.0$ CL = 20.8–25.3 | 1467 | $\bar{x} = 32.9$ CL = 29.5–36.4 | 1468 | $\bar{x} = 42.5$ CL = 37.7–47.4 |
| C2 | 1465 | $\bar{x} = 7.8$ CL = 6.7–8.9 | 1466 | $\bar{x} = 16.8$ CL = 15.0–18.5 | 1467 | $\bar{x} = 25.7$ CL = 23.3–28.1 | 1468 | $\bar{x} = 34.0$ CL = 29.9–38.2 |
| C3 | 1465 | $\bar{x} = 0.8$ CL = 0.6–1.0 | 1466 | $\bar{x} = 3.4$ CL = 2.5–4.3 | 1467 | $\bar{x} = 7.5$ CL = 5.6–9.5 | 1468 | $\bar{x} = 11.5$ CL = 9.0–14.0 |
| C4 | 1465 | $\bar{x} = 0.1$ CL = 0.1–0.2 | 1466 | $\bar{x} = 0.9$ CL = 0.6–1.1 | 1467 | $\bar{x} = 2.0$ CL = 1.4–2.6 | 1468 | $\bar{x} = 4.0$ CL = 3.0–4.9 |
| C3DC + C4OH | 1465 | $\bar{x} = 0.1$ CL = 0.0–0.1 | 1466 | $\bar{x} = 0.3$ CL = 0.2–0.4 | 1467 | $\bar{x} = 0.6$ CL = 0.4–0.8 | 1468 | $\bar{x} = 1.4$ CL = 1.0–1.9 |
| C5 | 1465 | $\bar{x} = 0.1$ CL = 0.0–0.1 | 1466 | $\bar{x} = 0.5$ CL = 0.4–0.6 | 1467 | $\bar{x} = 1.3$ CL = 1.0–1.6 | 1468 | $\bar{x} = 2.5$ CL = 1.9–3.2 |
| C5DC | 1465 | $\bar{x} = 0.1$ CL = 0.0–0.1 | 1466 | $\bar{x} = 0.6$ CL = 0.3–0.8 | 1467 | $\bar{x} = 1.0$ CL = 0.7–1.4 | 1468 | $\bar{x} = 2.5$ CL = 1.7–3.2 |
| C5OH | 1465 | $\bar{x} = 0.6$ CL = 0.4–0.8 | 1466 | $\bar{x} = 1.0$ CL = 0.8–1.3 | 1467 | $\bar{x} = 1.9$ CL = 1.4–2.4 | 1468 | $\bar{x} = 2.7$ CL = 2.1–3.4 |
| C6 | 1465 | $\bar{x} = 0.0$ CL = 0.0–0.0 | 1466 | $\bar{x} = 0.4$ CL = 0.3–0.5 | 1467 | $\bar{x} = 0.8$ CL = 0.6–1.0 | 1468 | $\bar{x} = 2.0$ CL = 1.4–2.6 |
| C8 | 1465 | $\bar{x} = 0.0$ CL = 0.0–0.0 | 1466 | $\bar{x} = 0.4$ CL = 0.4–0.5 | 1467 | $\bar{x} = 0.9$ CL = 0.7–1.1 | 1468 | $\bar{x} = 2.2$ CL = 1.8–2.7 |
| C10 | 1465 | $\bar{x} = 0.0$ CL = 0.0–0.0 | 1466 | $\bar{x} = 0.5$ CL = 0.4–0.7 | 1467 | $\bar{x} = 1.0$ CL = 0.8–1.2 | 1468 | $\bar{x} = 2.7$ CL = 2.0–3.5 |
| C12 | 1465 | $\bar{x} = 0.0$ CL = 0.0–0.0 | 1466 | $\bar{x} = 0.4$ CL = 0.3–0.6 | 1467 | $\bar{x} = 0.8$ CL = 0.6–1.0 | 1468 | $\bar{x} = 2.1$ CL = 1.6–2.5 |
| C14 | 1465 | $\bar{x} = 0.0$ CL = 0.0–0.1 | 1466 | $\bar{x} = 0.5$ CL = 0.4–0.5 | 1467 | $\bar{x} = 1.4$ CL = 1.1–1.6 | 1468 | $\bar{x} = 2.7$ CL = 2.2–3.3 |
| C16 | 1465 | $\bar{x} = 0.7$ CL = 0.6–0.8 | 1466 | $\bar{x} = 3.3$ CL = 2.7–3.8 | 1467 | $\bar{x} = 7.6$ CL = 6.5–8.7 | 1468 | $\bar{x} = 10.7$ CL = 9.2–12.2 |
| C16OH | 1465 | $\bar{x} = 0.0$ CL = 0.0–0.0 | 1466 | $\bar{x} = 0.1$ CL = 0.1–0.2 | 1467 | $\bar{x} = 0.3$ CL = 0.2–0.3 | 1468 | $\bar{x} = 0.6$ CL = 0.4–0.7 |
| C18 | 1465 | $\bar{x} = 0.6$ CL = 0.5–0.7 | 1466 | $\bar{x} = 1.5$ CL = 1.3–1.6 | 1467 | $\bar{x} = 2.4$ CL = 2.0–2.8 | 1468 | $\bar{x} = 5.0$ CL = 4.1–5.9 |
| C18OH | 1465 | $\bar{x} = 0.0$ CL = 0.0–0.0 | 1466 | $\bar{x} = 0.3$ CL = 0.2–0.3 | 1467 | $\bar{x} = 0.5$ CL = 0.4–0.6 | 1468 | $\bar{x} = 0.8$ CL = 0.6–0.9 |

Note: The values provided in the above tables are for reference use only. The mean value and confidence limits (CL) are determined by CDC for each Quality Control (QC) lot. Each participating laboratory must establish its own mean values and CL for its test method with these QC materials. Temporary estimates of mean values and CL can be determined after 10 successive, independent measurements. Slazyk WE, Hannon WH. *Quality assurance in the newborn screening laboratory*. In: Therrell BL Jr, editor. *Laboratory methods for neonatal screening*. Washington (DC): American Public Health Association, 1993:23-46.

Newborn Screening Quality Assurance Program
 Acylcarnitines Quality Control Specimen Certification
 Set 1— January 2015 Transition Set
Previous Lot Transition Materials (parallel testing)

Acylcarnitines Method: MS/MS Non-Derivatized - MS/MS non-kit

ENRICHMENT LEVELS (endogenous levels not included)

| <i>Analyte (μmol/L whole blood)</i> | <i>Lot</i> | <i>Base</i> | <i>Lot</i> | <i>Low</i> | <i>Lot</i> | <i>Intermediate</i> | <i>Lot</i> | <i>High</i> |
|--|------------|-------------|------------|------------|------------|---------------------|------------|-------------|
| Free carnitine (C0) | 1461 | 0 | 1462 | 10.0 | 1463 | 20.0 | 1464 | 30.0 |
| Acetylcarnitine (C2) | 1461 | 0 | 1462 | 10.0 | 1463 | 20.0 | 1464 | 30.0 |
| Propionylcarnitine (C3) | 1461 | 0 | 1462 | 3.0 | 1463 | 7.5 | 1464 | 12.0 |
| Butyrylcarnitine (C4) | 1461 | 0 | 1462 | 1.0 | 1463 | 2.5 | 1464 | 5.0 |
| Malonylcarnitine + Hydroxybutyrylcarnitine (C3DC + C4OH) | 1461 | 0 | 1462 | 1.0 | 1463 | 2.5 | 1464 | 5.5 |
| Isovalerylcarnitine (C5) | 1461 | 0 | 1462 | 0.5 | 1463 | 1.5 | 1464 | 3.0 |
| Glutaryl carnitine (C5DC) | 1461 | 0 | 1462 | 0.5 | 1463 | 1.0 | 1464 | 2.5 |
| Hydroxyisovalerylcarnitine (C5OH) | 1461 | 0 | 1462 | 0.5 | 1463 | 1.5 | 1464 | 2.5 |
| Hexanoylcarnitine (C6) | 1461 | 0 | 1462 | 0.5 | 1463 | 1.0 | 1464 | 2.5 |
| Octanoylcarnitine (C8) | 1461 | 0 | 1462 | 0.5 | 1463 | 1.0 | 1464 | 2.5 |
| Decanoylcarnitine (C10) | 1461 | 0 | 1462 | 0.5 | 1463 | 1.0 | 1464 | 2.5 |
| Dodecanoylcarnitine (C12) | 1461 | 0 | 1462 | 0.5 | 1463 | 1.0 | 1464 | 2.5 |
| Myristoylcarnitine (C14) | 1461 | 0 | 1462 | 0.5 | 1463 | 1.5 | 1464 | 3.0 |
| Palmitoylcarnitine (C16) | 1461 | 0 | 1462 | 3.0 | 1463 | 8.0 | 1464 | 12.0 |
| Hydroxypalmitoylcarnitine (C16OH) | 1461 | 0 | 1462 | 0.1 | 1463 | 0.5 | 1464 | 1.0 |
| Stearoylcarnitine (C18) | 1461 | 0 | 1462 | 1.0 | 1463 | 2.0 | 1464 | 5.0 |
| Hydroxystearoylcarnitine (C18OH) | 1461 | 0 | 1462 | 0.5 | 1463 | 1.0 | 1464 | 1.5 |

Newborn Screening Quality Assurance Program
 Acylcarnitines Quality Control Specimen Certification
 Set 1— January 2015 Transition Set
Previous Lot Transition Materials (parallel testing)

Acylcarnitines Method: MS/MS Non-Derivatized - MS/MS non-kit

ANALYTICAL INFORMATION *Lot Numbers, Mean Values (\bar{x} , $\mu\text{mol/L}$ whole blood), and 95% Confidence Limits (CL)*

| <i>Analyte</i> | <i>Lot</i> | <i>Mean/ 95% CL</i> | <i>Lot</i> | <i>Mean/ 95% CL</i> | <i>Lot</i> | <i>Mean/ 95% CL</i> | <i>Lot</i> | <i>Mean/ 95% CL</i> |
|----------------|------------|------------------------------------|------------|------------------------------------|------------|------------------------------------|------------|------------------------------------|
| C0 | 1461 | $\bar{x} = 13.0$ CL = 11.2–14.7 | 1462 | $\bar{x} = 21.1$ CL = 18.6–23.6 | 1463 | $\bar{x} = 29.3$ CL = 25.9–32.6 | 1464 | $\bar{x} = 38.1$ CL = 33.1–43.1 |
| C2 | 1461 | $\bar{x} = 12.4$ CL = 11.0–13.7 | 1462 | $\bar{x} = 22.0$ CL = 19.4–24.6 | 1463 | $\bar{x} = 32.2$ CL = 28.6–35.9 | 1464 | $\bar{x} = 41.8$ CL = 36.4–47.2 |
| C3 | 1461 | $\bar{x} = 1.2$ CL = 0.9–1.5 | 1462 | $\bar{x} = 3.9$ CL = 3.1–4.8 | 1463 | $\bar{x} = 7.9$ CL = 5.5–10.2 | 1464 | $\bar{x} = 11.6$ CL = 8.9–14.3 |
| C4 | 1461 | $\bar{x} = 0.1$ CL = 0.1–0.2 | 1462 | $\bar{x} = 1.0$ CL = 0.7–1.2 | 1463 | $\bar{x} = 2.2$ CL = 1.8–2.7 | 1464 | $\bar{x} = 4.4$ CL = 3.2–5.5 |
| C3DC + C4OH | 1461 | $\bar{x} = 0.1$ CL = 0.0–0.1 | 1462 | $\bar{x} = 0.3$ CL = 0.2–0.4 | 1463 | $\bar{x} = 0.5$ CL = 0.4–0.7 | 1464 | $\bar{x} = 1.1$ CL = 0.8–1.4 |
| C5 | 1461 | $\bar{x} = 0.1$ CL = 0.1–0.1 | 1462 | $\bar{x} = 0.6$ CL = 0.4–0.7 | 1463 | $\bar{x} = 1.4$ CL = 1.1–1.7 | 1464 | $\bar{x} = 2.7$ CL = 2.1–3.3 |
| C5DC | 1461 | $\bar{x} = 0.1$ CL = 0.0–0.2 | 1462 | $\bar{x} = 0.7$ CL = 0.4–1.0 | 1463 | $\bar{x} = 1.2$ CL = 0.7–1.8 | 1464 | $\bar{x} = 2.8$ CL = 2.0–3.6 |
| C5OH | 1461 | $\bar{x} = 0.8$ CL = 0.5–1.0 | 1462 | $\bar{x} = 1.2$ CL = 0.8–1.6 | 1463 | $\bar{x} = 2.1$ CL = 1.4–2.8 | 1464 | $\bar{x} = 3.1$ CL = 2.4–3.9 |
| C6 | 1461 | $\bar{x} = 0.0$ CL = 0.0–0.1 | 1462 | $\bar{x} = 0.5$ CL = 0.3–0.6 | 1463 | $\bar{x} = 0.9$ CL = 0.7–1.2 | 1464 | $\bar{x} = 2.2$ CL = 1.8–2.7 |
| C8 | 1461 | $\bar{x} = 0.0$ CL = 0.0–0.1 | 1462 | $\bar{x} = 0.6$ CL = 0.4–0.8 | 1463 | $\bar{x} = 1.1$ CL = 0.9–1.4 | 1464 | $\bar{x} = 2.6$ CL = 2.0–3.1 |
| C10 | 1461 | $\bar{x} = 0.1$ CL = 0.0–0.1 | 1462 | $\bar{x} = 0.7$ CL = 0.4–0.9 | 1463 | $\bar{x} = 1.1$ CL = 0.8–1.3 | 1464 | $\bar{x} = 2.7$ CL = 2.0–3.4 |
| C12 | 1461 | $\bar{x} = 0.0$ CL = 0.0–0.0 | 1462 | $\bar{x} = 0.4$ CL = 0.4–0.5 | 1463 | $\bar{x} = 0.9$ CL = 0.7–1.1 | 1464 | $\bar{x} = 2.2$ CL = 1.7–2.8 |
| C14 | 1461 | $\bar{x} = 0.0$ CL = 0.0–0.1 | 1462 | $\bar{x} = 0.6$ CL = 0.4–0.7 | 1463 | $\bar{x} = 1.4$ CL = 1.1–1.7 | 1464 | $\bar{x} = 2.8$ CL = 2.2–3.5 |
| C16 | 1461 | $\bar{x} = 1.0$ CL = 0.8–1.1 | 1462 | $\bar{x} = 3.3$ CL = 2.6–3.9 | 1463 | $\bar{x} = 7.2$ CL = 5.6–8.8 | 1464 | $\bar{x} = 10.0$ CL = 7.9–12.1 |
| C16OH | 1461 | $\bar{x} = 0.0$ CL = 0.0–0.0 | 1462 | $\bar{x} = 0.1$ CL = 0.0–0.1 | 1463 | $\bar{x} = 0.3$ CL = 0.2–0.4 | 1464 | $\bar{x} = 0.5$ CL = 0.3–0.7 |
| C18 | 1461 | $\bar{x} = 0.7$ CL = 0.5–0.9 | 1462 | $\bar{x} = 1.5$ CL = 1.1–2.0 | 1463 | $\bar{x} = 2.3$ CL = 1.8–2.8 | 1464 | $\bar{x} = 4.9$ CL = 3.7–6.0 |
| C18OH | 1461 | $\bar{x} = 0.0$ CL = 0.0–0.0 | 1462 | $\bar{x} = 0.3$ CL = 0.1–0.4 | 1463 | $\bar{x} = 0.5$ CL = 0.2–0.8 | 1464 | $\bar{x} = 0.8$ CL = 0.3–1.2 |

Note: The values provided in the above tables are for reference use only. The mean value and confidence limits (CL) are determined by CDC for each Quality Control (QC) lot. Each participating laboratory must establish its own mean values and CL for its test method with these QC materials. Temporary estimates of mean values and CL can be determined after 10 successive, independent measurements. *Slazyk WE, Hannon WH. Quality assurance in the newborn screening laboratory. In: Therrell BL Jr, editor. Laboratory methods for neonatal screening. Washington (DC): American Public Health Association, 1993:23-46*