



# LABORATORY INFORMATION TO MONITOR WHEN RESULTS HAVE THE POTENTIAL TO BE FALSE-POSITIVE

Review of patient's past results and current test orders at time of testing

---

Review of patient's past results at time of reporting

---

For newly diagnosed patients, check if there were other positive specimens in the same processing batch

---

Periodic monitoring of all test results

---

Documentation to determine possible false-positive causes (e.g., specimen processing logs and lot numbers)

---

Monitoring of results triggered by external inquiries

---

Use the below section to fill in additional laboratory information to monitor

---

---

---

---

---

---

